Instruction/
Operation Manual

Restek Methanizer for Agilent GCs

(Cat.# 22650)

RESTEK

Pure Chromatography

www.restek.com

Table of Contents

1.0 Overview	3
2.0 Specifications	3
3.0 Parts Identification Guide	4
4.0 Plumbing the Hydrogen Source Supply Line.....	5
5.0 Installation Guide	6–9
5.1 Methanizer Controller Box and Heater Chamber Wires	6
5.2 Inside the GC Oven	6–8
5.3 Column Installation	8
5.4 Methanizer Controller Box Connections.....	9
6.0 Operating Procedure	10
7.0 Front Panel Buttons	10
8.0 LCD Screen Modes	11
8.1 Main Screens	11
8.2 Fault Screens	12
9.0 Replacing the Catalyst Tube	13
10.0 Cautions	13
11.0 Maintenance and Service	13
12.0 Back Label Explanation.....	14
13.0 EC Declaration of Conformity	43


Restek Methanizer for Agilent GCs

Installation/Operation Guide*

1.0 Overview

The Restek methanizer enables a flame ionization detector (FID) to detect low levels of carbon monoxide (CO) and carbon dioxide (CO₂). The gas mixture is first separated using a GC column, and then it passes over a hot catalyst that (in the presence of hydrogen) converts any CO and CO₂ into methane. Methane (CH₄) is readily detected by the FID. A typical setup is shown in Figure 1.

Figure 1: Restek Methanizer Installed in a GC-FID

2.0 Specifications

Input Power Rating: 100-240 VAC, 50/60 Hz, 2.0 A

Output Power Rating: 24 V $\overline{\text{DC}}$, 3.5 A

Operating Temperature Range for Control Box: 32–120 °F (0–48 °C)

Operating Temperature Range for Catalyst and Heater Chamber Assembly: 350–450 °C (factory preset at 380 °C)

Controller Temperature Accuracy: ± 1 °C

Warranty: 1 year

Certification/Compliance: CE, Canadian ICES-001, WEEE, RoHS 2, China RoHS 2, UKCA

Indoor Use Only

Altitude up to 2000 m

Maximum relative humidity: 80% for temperatures up to 31 °C (decreases linearly to 50% relative humidity at 40 °C)

Main supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage

Transient overvoltage up to the levels of overvoltage category II

Pollution degree 2

NOTE: If equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

* Cannot be used with 6850 GCs

3.0 Parts Identification Guide

- Methanizer control box
- Heater chamber (cat.# 27210)
- Catalyst tube (cat.# 27209)
- Hydrogen supply line (cat.# 27208)
- Adaptor for Agilent capillary FID (cat.# 27211)
- Adaptor for Agilent packed FID (cat.# 27212)
- 1/8" Swagelok nut (cat.# 23151, not shown)
- Ferrules (cat.# 20215, 20219, 20220, 27067, not shown)
- Right angle wrench (cat.# 22649, not shown)
- Transition tubing (not shown)

Parts Manufactured by Restek Corporation

Methanizer control box

Heater chamber
(cat.# 27210)

Hydrogen supply line
(cat.# 27208)

Catalyst tube
(cat.# 27209)

Adaptor for Agilent capillary FID
(cat.# 27211)


Adaptor for Agilent packed FID
(cat.# 27212)

Recommended Related Products (Available at www.restek.com)

- Flowmeter (cat.# 22656) and FID flow measuring adaptor (cat.# 21000)
- Leak detector (cat.# 22655)
- Hi-Duty tubing cutter (cat.# 22356)
- Replacement triple gas filters (cat.# 22020)
- Filter baseplate (cat.# 22344)
- Installation kit (cat.# 27213)
 - 1/8" Metering valve (cat.# 23204)
 - 1/8" Stainless-steel gas supply line (10 feet) (cat.# 29026)
 - 1/8" Tee (cat.# 23172)
 - 1/8"-1/16" Tube end reducer (cat.# 23176)
 - 1/8" Stainless-steel ball valve (cat.# 23200)
 - MINICYL regulator with 1/8" stainless-steel fittings (cat.# 27215)

4.0 Plumbing the Hydrogen Source Supply Line

Figure 2: Components and Setup

Warning: Use caution when working with hydrogen. Hydrogen mixed with air can be both flammable and an explosion hazard. Turn off the hydrogen supply at its source before working on the instrument.

Note: Use new, clean tubing for the installation of the gas line. Hydrogen can remove deposits on the inside of old tubing, and this may affect the operation of the methanizer.

Note: Use stainless-steel tubing for the hydrogen gas line. Copper can become brittle when exposed to hydrogen gas, which could create a safety hazard.

Note: Use hydrogen that is of carrier gas purity (99.999% or greater) and use gas filters to remove hydrocarbons, water, and oxygen.

Figure 2 shows the ideal use of parts to plumb the hydrogen source line to the hydrogen supply line of the methanizer. If the operator chooses to run a dedicated gas supply line from a gas bottle or generator for the methanizer, then the tee is not needed. If the operator chooses to splice into the existing gas line supplying hydrogen to the FID, the following procedure should be used.

- 4.1 Verify hydrogen is switched off at its source.
- 4.2 Locate the hydrogen source supply line at the back of the GC and install a 1/8" tee after the filter on the hydrogen gas line.
- 4.3 Install the ball valve after the tee using a short length of 1/8" stainless-steel tubing.
- 4.4 Install the MINICYL regulator after the ball valve using a short length of 1/8" stainless-steel tubing.
- 4.5 Install the metering gas valve using a short length of 1/8" stainless-steel tubing.
- 4.6 Install the 1/8" to 1/16" tube end reducer into the output end of the metering gas valve.

NOTE: A hydrogen supply line (cat.# 27208) is included with the methanizer. It will be connected to the tube end reducer in the installation section of this manual.

5.0 Installation Guide

5.1 Methanizer Controller Box and Heater Chamber Wires

5.1.1 Switch off all of the supply gases at their sources.

5.1.2 Cool down the GC oven and the detector.

5.1.3 Place the methanizer controller box on the top right side of the GC oven. This will make the power switch and connectors more accessible.

NOTE: Do not place the unit where it will be difficult to access the power cord or the on/off switch. Also, when positioning the unit, ensure that the ventilation slots are not blocked.

5.1.4 Wearing gloves, remove one of the cutouts at the top of the GC oven and move insulation aside (Figure 3).

5.1.5 Route the two connectors on the heating chamber out through both the cutout hole and the insulation toward the back of the methanizer controller box. Do NOT connect them to the control box at this time. Pull any excess wire through the hole (Figure 4).

Figure 3: Removing GC Oven Cutout and Insulation


Figure 4: Routing Connectors

5.2 Inside the GC Oven


5.2.1 Prepare the methanizer by inserting the catalyst tube into the heater chamber as shown in Figure 5. The arrow on the heater chamber indicates the top of the unit.

5.2.2 Make sure the support ring on the catalyst tube is at the bottom of the heater chamber.

Figure 5: Heater Chamber and Catalyst Tube


Figure 6: Adaptor and FID Weldment

5.2.3 Install the adaptor for FID as follows.

- For capillary FID weldment (Figure 6):
 - Place the adaptor for Agilent capillary FID and a 0.8 mm ferrule (cat.# 27067) over the transition tubing. Wipe the end of the transition tubing clean with a tissue dampened with isopropanol to remove fingerprints and ferrule pieces before inserting it into the FID.
 - Thread the adaptor for Agilent capillary FID and the transition tubing into the FID detector base.
 - Insert the transition tubing all the way into the detector jet until it stops, then withdraw it 1-2 mm. Tighten the adaptor until finger tight. Use a wrench to tighten the adaptor an additional $\frac{1}{4}$ turn.
- For packed column FID weldment:
 - Place the back and front ferrules over the adaptor for Agilent packed FID (cat.# 27212). Insert the adaptor all the way into the detector.
 - Tighten the adaptor nut until finger tight. Use a wrench to tighten the adaptor nut an additional $\frac{1}{2}$ to $\frac{3}{4}$ turn.

Figure 7: Component Parts inside GC Oven

- 5.2.4 Connect the top fitting of the catalyst tube to the adaptor used above. If using the adaptor for capillary FID, a small part of the transition tubing will extend into the catalyst tube. Tighten the nut until finger tight. (See Figure 7 for reference.)
- 5.2.5 Using two wrenches against each other, tighten the nut an additional $\frac{1}{2}$ to $\frac{3}{4}$ turn to be sure the catalyst tube is tightly connected to the adaptor nut.
- 5.2.6 Feed the methanizer hydrogen supply line (cat.# 27208) out of the GC through the same cutout opening through which the wires of the heater chamber were routed.
- 5.2.7 Connect the hydrogen supply line (cat.# 27208) to the $\frac{1}{8}$ " to $\frac{1}{16}$ " tube end reducer connected to the metering valve.
- 5.2.8 Connect the other end of the hydrogen supply line (cat.# 27208) to the $\frac{1}{8}$ " fitting on the bottom of the catalyst tube. Tighten the nut until finger tight, then, using two wrenches against each other to avoid twisting the hydrogen supply line tubing (cat.# 27208), tighten the nut an additional $\frac{1}{2}$ to $\frac{3}{4}$ turn.

5.3 Column Installation

- 5.3.1 Install the analytical column between the inlet side of the GC oven (or sampling valve) and the methanizer column inlet. (An installed column is shown in Figure 10.)
- 5.3.2 Use the column fitting and the reducing ferrules that are provided (Table I).
- Packed column installation: Place the column nut and the ferrule over one end of the column. Insert the column all the way into the methanizer column inlet until it stops and then withdraw it 1-2 mm. Tighten the column nut until finger tight and then, using two wrenches against each other to avoid twisting the hydrogen supply line (cat.# 27208) tube, tighten the nut an additional $\frac{1}{4}$ to $\frac{1}{2}$ turn.
 - Capillary and micropacked column installation: Place the column nut and ferrule over one end of the column. The optimal insertion distance of the column into the methanizer inlet should be 1" from the bottom of the column nut to the top of the column. With the column in the proper position, finger tighten the column nut. Then, using two wrenches against each other to avoid twisting the hydrogen supply line (cat.# 27208), tighten the nut an additional $\frac{1}{4}$ to $\frac{1}{2}$ turn.

Table I: Fitting/Ferrule/Column Compatibility*

Fitting	Ferrule ID	Fits Column	Quantity per Pack	Catalog Number
$\frac{1}{8}$ "	0.8 mm	0.53 mm ID	5	20215
$\frac{1}{8}$ "	$\frac{1}{16}$ "	1 mm ID	10	20220
$\frac{1}{8}$ "	$\frac{1}{8}$ "	$\frac{1}{8}$ " OD	10	20219

*If you are using other column dimensions, contact Technical Service for fitting/ferrule recommendations.

5.4 Methanizer Controller Box Connections (Figures 8 and 9)

5.4.1 Be sure the power switch is turned off on the power entry module of methanizer controller box. (O = off, I = on)

5.4.2 Connect the H1 output to H1. Align the key on the connector's insert and turn it clockwise ¼ turn.

5.4.3 Connect the TC1 input from the thermocouple to TC1. Push the connector in until it stops.

5.4.4 Plug the VAC cord into the power entry module of methanizer controller box.

5.4.5 Leave the power OFF until you are ready to perform a test.

5.4.6 Connect the power cord to a properly grounded outlet.

5.4.7 DO NOT block the ventilation slots on either side of the methanizer.

Figure 8

H1= Output to heater
TC1= Input from thermocouple

Figure 9



Figure 10: Restek Methanizer Installed in a GC-FID

6.0 Operating Procedure

- 6.1 Start the system by initiating the flow of hydrogen. Set the MINICYL regulator to 50 psi and adjust the flow to 25 mL/min with the metering valve. Verify the flow using a Restek flowmeter (cat.#22656) and an FID flow measuring adaptor (cat.# 21000).

Warning: Leak check all the new connections using a Restek leak detector (cat.# 22655). If a leak is detected, tighten the connection an additional 1/8 turn.

- 6.2 Turn the carrier gas on and set it to an appropriate value.
- 6.3 Let the system purge using the methanizer hydrogen supply and carrier gas for at least 10 minutes.
- 6.4 Turn on the methanizer by switching the power button located above the power cord on the back of the unit to the ON position.
- 6.5 The methanizer is factory preset to 380 °C for optimal catalytic efficiency. However, this set point temperature can be adjusted using the up and down arrow buttons to the right of the digital display (Figure 11).
- 6.6 Push the RUN button on the methanizer unit to start heating the catalyst tube. The actual temperature changes are displayed on the digital display (Figures 12 and 13).
- 6.7 Allow system to stabilize for 1 hour.
- 6.8 During the stabilization period, adjust the inlet and FID temperature settings and light the FID (Table II).
- 6.9 Condition the catalyst tube by setting the GC oven temperature to 200 °C. Do not exceed the maximum temperature of your GC column.

Table II: Recommended FID Gases for Standard Installation

Gas	Flow Rate (mL/min)
Carrier	Gas and column dependent
Hydrogen - FID	15 (25 mL/min is already coming from the methanizer supply gas, giving a total flow of around 40 mL/min [10:1 air to hydrogen])
Air - FID	400

7.0 Front Panel Buttons

- 7.1 The buttons are capacitive touch buttons (similar to the touch screen on your smart phone). They do not have a mechanical detent when they are touched. When a button is touched, the lights on the front panel will flash, indicating the touch was recognized.
- 7.2 The MODE button allows you to change screens.
- 7.3 The RUN button turns the heater on and off.
- 7.4 The up and down arrows change the set point of the methanizer.

8.0 LCD Screen Modes

8.1 Main Screens

Figure 11: Heater Off


Figure 12: Heating



Figure 13: Set point Achieved

8.2 Fault Screens

- 8.2.1 If a fault is indicated, the unit will turn off the heater and will not again operate until the fault has been cleared by fixing the problem. The unit will not automatically restart.
- 8.2.2 When the problem has been fixed, the fault screen will change to NO FAULTS (Figure 14).
- 8.2.3 Touch the MODE button to return to the main screen and then touch the RUN button to turn the heater on and continue normal operation.
- 8.2.4 The Over Temp Condition fault indicator is shown in Figure 15. When this displays, check the thermocouple connector for possible thermocouple malfunction.
- 8.2.5 The Heater Open fault indicator is shown in Figure 16. When this displays, check the heater connection.

Figure 14: No Faults


Figure 15: Over Temp Condition Fault Indicator



Figure 16: Heater Open Fault Indicator

9.0 Replacing the Catalyst Tube

When conversion performance of the catalyst material becomes significantly degraded, replace the catalyst tube.

WARNING: The catalyst chamber and all associated fittings can be extremely HOT! Allow the actual temperature on the controller box and the GC oven to cool down to 50 °C or less before attempting to touch any part of this assembly.

- 9.1 Set the GC oven to 25 °C and cool down the methanizer heater chamber by pushing the RUN button on the controller box to turn the heater off. Observe the actual temperature displayed on the controller box display (Figure 11) and do not touch until the GC oven and catalyst heater chamber have cooled down to 50 °C or less.
- 9.2 Locate the hydrogen supply line and turn off the flow of hydrogen.
- 9.3 Using two wrenches against each other, remove the analytical column from the catalyst tube.
- 9.4 Using two wrenches against each other, loosen the bottom nut of the catalyst tube that is connected to the methanizer's hydrogen supply line. Remove the hydrogen supply line.
- 9.5 Using two wrenches, loosen the top nut of the catalyst tube connecting the adaptor for capillary (or packed) FID. Remove the heater chamber and slide the catalyst tube out of the bottom of the heater chamber.
- 9.6 Replace the catalyst tube with a new one, keeping the support ring of the catalyst tube towards the bottom of the heater chamber.
- 9.7 Reassemble the methanizer components following the directions in section 5.0 (Installation Guide).


10.0 Cautions

- 10.1 Even small amounts of sulfur compounds can poison the nickel catalyst, and a poisoned catalyst must be replaced as it cannot be regenerated. If sulfur compounds are present in a sample, the catalyst should be bypassed using a switching valve or by backflushing the column.
- 10.2 Hydrogen gas in the presence of the hot nickel catalyst is a reducing agent. Operators should bypass the catalyst when analyzing samples containing compounds that are easily reduced, such as carboxylic acids or unsaturated hydrocarbons.

11.0 Maintenance and Service

- 11.1 All maintenance and servicing of the controller box will be performed at Restek.
- 11.2 For your safety, DO NOT open the controller box. Opening the controller box will void the warranty. If the unit is opened, first disconnect it from the electrical supply in order to prevent a shock hazard.
- 11.3 For questions, problems, or repair services, contact Restek Technical Service at 1-800-356-1688 or 1-814-353-1300, ext. 4 within the U.S. Outside the U.S. please contact your Restek representative.
- 11.4 DO NOT replace main supply cords with inadequately rated cords. Only use cords supplied by Restek.
- 11.5 To prevent injury and/or burning yourself when replacing the catalyst, confirm that the catalyst chamber and GC oven have cooled to 50 °C or less prior to touching any parts. Note that the actual temperature of the catalyst is displayed on the front of the controller box.

12.0 Back Label Explanation

1. Product name.
2. Manufacturer company name, address, and contact information.
3. Country of origin
4. Product catalog number.
5. Product serial number.
6. Electrical parameters.
7. This Class A digital apparatus complies with Canadian ICES-001.
8. This complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
 - 1) This device may not cause harmful interference.
 - 2) This device must accept any interference received, including interference that may cause undesired operation.
9. This unit conforms to EU/EMC Directive 2004/108/EC; standards, to which conformity is declared include:
 - EN61326-1:2006 w/A3 Class A.
 - Unit is WEEE, RoHS 2, China RoHS 2, UKCA compliant
 - UKCA replaces CE mark for Britain.

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MODEL / TYPE: 22650

DESCRIPTION: Methanizer used in gas chromatography for trace analysis of CO and/or CO²

SERIAL NUMBER: TBD

REPORTS: F2LQ7444B-02S, F2LQ7444B-03S, F2LQ7444B-04S, F2LQ7444B-05S & F2LQ7444B-01E

DIRECTIVES: Low Voltage Directive 2006/95/EC, EMC Directive 2004/108/EC

STANDARDS CONSIDERED, FULL OR IN PART:

- **EN 61010-1:2010** – Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements
- **EN 61010-2-010:2014** – Safety requirement for electrical equipment for measurement and control and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials.
- **EN 61326-1:2013** – Electrical equipment for measurement, control, and laboratory use – EMC requirements – Part 1: General requirements
- **EN 55011:2009+A1:2010** – Industrial, scientific, and medical equipment – Radio-frequency disturbance characteristics – Limits and methods of measurement
- **EN 61000-3-2:2006, inc. A1:2009 & A2:2009** – Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current $\leq$ 16A per phase)
- **EN 61000-3-3:2013** – Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current $\leq$ 16A per phase and not subject to conditional connection
- **UKCA** – Supply of Machinery (Safety) Regulations 2008, the Electromagnetic Compatibility Regulations 2016 and the Radio Equipment Regulations 2017

THIRD PARTY TESTING:

26501 Ridge Road
Damascus, MD 20872
Phone: 301-253-4500
Fax: 301-253-5179
Toll Free: 877-405-1580

16740 Peters Road
Middlefield, OH 44062
Phone: 440-632-5541
Fax: 440-632-5542

www.f2labs.com

Restek Corporation declares under its sole responsibility that the Methanizer system, model 22650 is in conformity with the Low Voltage Directive 2006/95/EC and the EMC Directive 2004/108/EC.

Authorized by:

Name: Gerald Johnston

Date: Jun 15, 2021

Title: Product Development Engineer

For the most up-to-date information, see our website www.restek.com

For questions, problems, or repair services:

Within the U.S.:

Call Restek Technical Service at 1-800-356-1688 or 1-814-353-1300, ext. 4.

Outside the U.S.:

Contact your local Restek representative.

Questions about this or any other Restek product?
Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2021 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com

#203-04-003 Rev. date: 06/21

RESTEK
Pure Chromatography