


Split Vent Trap Kit

cat.# 20698, 20699

Installation Procedure

Restek Recommends: To ensure optimal bed packing and performance, and to minimize the likelihood of material “channeling,” gas filters should be installed vertically.

1. Connect the 9-inch piece of copper tubing to the GC split vent outlet.
2. Connect the split vent trap, in either direction, to the 9-inch piece of copper tubing.
3. Connect the 3-inch piece of copper tubing to the other end of the split vent trap.
4. Position the split vent trap vertically. The trap can be mounted on the side of the GC with the enclosed Velcro strip.
5. Adjust the column split vent flow rates. (The split vent trap slightly increases the backpressure on the split vent outlet causing the flow rate through the column and split vent to decrease approximately 5mL/min.)



Questions about this or any other Restek product?
Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com

#204-08-015 Rev. date: 12/24



RESTEK
Pure Chromatography