


# Installing Ashcroft Gauges on Canisters Using Stainless Steel Connectors


cat.# 22121, 24268, 24285



1. Verify you have received the following parts: a 1/4" x 1/4" NPT female connector with a nut and ferrules, a 1/4" nut, a 1/4" port connector.



2. Place the nut on the port connector against the base of the flared end. Leave the nut and ferrules on the female connector.



3. Install the port connector and nut onto the valve and tighten 1/4 turn past finger-tight.



4. Wrap three layers of PTFE tape around the threads on the gauge. Thread the female connector onto the gauge and tighten.



5. Finally, slide the gauge with the female connector attached onto the port connector and tighten 1/4 turn past finger-tight. You are now ready to open the valve and take vacuum and pressure readings.

**Questions about this or any other Restek product?**  
**Contact us or your local Restek representative ([www.restek.com/contact-us](http://www.restek.com/contact-us)).**

Restek patents and trademarks are the property of Restek Corporation. (See [www.restek.com/Patents-Trademarks](http://www.restek.com/Patents-Trademarks) for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

[www.restek.com](http://www.restek.com)

**#500-09-001** Rev. date: 08/24



**RESTEK**  
Pure Chromatography