



# Tubing Assembly for ASE 200/300 Systems

cat.# 26251, 26248

Restek's tubing assembly for ASE 200/300 systems eliminates the need for an adaptor fitting and allows  $\frac{1}{32}$ " tubing to be inserted directly into either static or pressure relief valves. This simplifies installation and removes a potential source of leaks.

## Installation Instructions

1. Remove the old tubing assembly and adaptor from the static or pressure relief valve. The adaptor is not needed when using Restek's tubing assembly for ASE 200/300 systems (Figure 1).
2. Install the  $\frac{1}{16}$ " tubing side of Restek's tubing assembly for ASE 200/300 systems directly into the static or pressure relief valve, whichever applies (Figure 2).
3. Install the  $\frac{1}{32}$ " tubing side of Restek's tubing assembly for ASE 200/300 systems into the upper or lower autoseal tip, whichever applies.



**Questions about this or any other Restek product? Contact us or your local Restek representative ([www.restek.com/contact-us](http://www.restek.com/contact-us)).**

Restek patents and trademarks are the property of Restek Corporation. (See [www.restek.com/Patents-Trademarks](http://www.restek.com/Patents-Trademarks) for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2020 Restek Corporation. All rights reserved. Printed in the U.S.A.

[www.restek.com](http://www.restek.com)

#803-04-002 Rev. date: 07/20



**RESTEK**  
Pure Chromatography