



# A Novel Solution for Vitamin K<sub>1</sub> and K<sub>2</sub> Analysis in Human Plasma by LC-MS/MS

By Shun-Hsin Liang and Frances Carroll


## Abstract

Vitamin K<sub>1</sub> and K<sub>2</sub> analysis is typically complex and time-consuming because these lipophilic vitamins occur at low levels and are subject to matrix interference, particularly from phospholipids. The new method established here provides fast, accurate, and precise vitamin K analysis in plasma using a simple phospholipid removal procedure and a sensitive, 4-minute LC-MS/MS analysis on a Raptor Biphenyl column, making it suitable for both research and high-throughput testing.

## Introduction

Vitamin K is a group of fat-soluble vitamins divided into Vitamin K<sub>1</sub> (one compound, phyloquinone) and K<sub>2</sub> (a group of compounds, menaquinones). Among the K<sub>2</sub> family, MK4 and MK7 are the most nutritionally recognized. While K<sub>1</sub> plays an important role in controlling blood clotting, studies have shown that MK4 and MK7 have distinct biological function in the regulation of bone metabolism and vascular calcification. As interest in their biological action in extra-hepatic tissues is increasing, an accurate and simple measurement of vitamin K status remains a critical issue for both clinical research and diagnostics. Vitamin K<sub>1</sub> and K<sub>2</sub> analysis is extremely challenging because they are the most lipophilic and least abundant of the fat-soluble vitamins. Common strategies are complex and time-consuming: typical methods use liquid-liquid extraction paired with silica-based SPE to separate polar and nonpolar lipids, followed by post-column derivatization and fluorescence detection. In this study, a simple and fast phospholipid removal sample preparation procedure was developed and used in combination with an LC-MS/MS analysis using a Raptor Biphenyl column. The method established here provides a novel solution for vitamin K measurement for both research and high-throughput clinical diagnostic labs.

**Figure 1:** Structures of Vitamin K<sub>1</sub> and K<sub>2</sub>



## Experimental

### Calibration Standards and Quality Control Samples

SeraFlx BIOMATRIX (BioIVT) was fortified with the three vitamin K analytes to prepare calibration standards and QC samples. Seven-point calibration curves established a 0.10-10 ng/mL linear range. Three QC levels were prepared at 0.25, 0.75, and 8.00 ng/mL for all three vitamin K analytes. The fortified standards and QC samples were subjected to sample preparation procedure as described further on.

### Human Plasma Analysis

To ensure method applicability across a range of patient populations that vary in endogenous vitamin K levels, multiple human plasma samples (BioIVT) were analyzed. These samples included pooled plasma (from 20 people), plasma from an individual taking vitamin K<sub>2</sub> medication, and plasma from an individual of a Japanese descent because the typical Japanese diet is high in MK7-containing foods.

### Sample Preparation

A 500 µL aliquot of sample was mixed with 5 µL of internal standard solution (K<sub>1</sub>-d7, MK4-d7, and MK7-d7 at 100 ng/mL in methanol) and 1.5 mL of acetonitrile followed by vortexing for 20 seconds at 3000 rpm. After centrifugation at 4300 rpm for 10 minutes, the supernatant was loaded onto a Biotage ISOLUTE PLD+ 96-well plate (50 mg) and vacuum was applied to collect the eluate. The eluate was then evaporated to dryness at 50 °C under a gentle stream of nitrogen. The dried extract was reconstituted with 100 µL of 15:85 water:methanol and 5 µL of sample was injected for analysis.

### Instrument Conditions

LC-MS/MS analysis of vitamin K<sub>1</sub> and K<sub>2</sub> in human plasma was performed on a Waters Xevo TQ-S mass spectrometer paired with an ACQUITY UPLC system. Instrument conditions were as follows, and analyte transitions are provided in Table I.

Analytical column: Raptor Biphenyl 2.7 µm,  
50 mm x 2.1 mm (cat.# 9309A52)  
Mobile phase A: Water, 0.1% formic acid,  
5 mM ammonium formate  
Mobile phase B: Methanol, 0.1% formic acid  
Gradient: 

Time (min)	%B
0.00	90
1.00	100
3.00	100
3.01	90
4.00	90

  
Flow rate: 0.4 mL/min  
Injection volume: 5 µL  
Column temp.: 40 °C  
Ion mode: Positive ESI

**Table I: Analyte Transitions**

Compound	Precursor Ion	Product Ion Quantifier	Product Ion Qualifier
MK4	445.5	187.2	81.1
K <sub>1</sub>	451.5	187.2	128.2
MK7	649.7	187.2	81.1
MK4-d7	452.5	194.2	—
K <sub>1</sub> -d7	458.5	194.2	—
MK7-d7	656.7	194.2	—


## Results and Discussion

### Chromatographic Performance


A fast, 4-minute chromatographic vitamin K<sub>1</sub> and vitamin K<sub>2</sub> analysis was achieved with injection of reconstituted solution obtained from a simple phospholipid removal procedure. No chromatographic interferences were observed from the analysis of blank plasma samples (Figure 2), indicating the technique provided sufficient removal of matrix components.

**Figure 2: Analysis of Blank and Fortified SeraFlx BIOMATRIX**

**Blank**




**0.5 ng/mL Fortified**



### Linearity

Using 1/x weighted linear regression, all three analytes showed acceptable linearity with  $r^2$  values of 0.995 or greater, and deviations of <15% from the nominal concentrations. The established limit of quantitation was 0.10 ng/mL. Representative standard curves are shown in Figure 3.

**Figure 3: Standard Curves**



### Accuracy and Precision

Precision and accuracy testing for vitamin K<sub>1</sub> and K<sub>2</sub> analysis according to the method established here was performed on three different days. The method was demonstrated to be accurate with recovery values falling within 10% of the nominal concentration for all QC levels. The %RSD was 0.207-7.83% and 3.39-5.75% for intraday and interday assessments, respectively, indicating acceptable method precision (Table II).


**Table II:** Accuracy and Precision of QC Samples for Vitamin K<sub>1</sub> and K<sub>2</sub> Analysis

Analyte	QC-1 (0.250 ng/mL)			QC-2 (0.750 ng/mL)			QC-3 (8.00 ng/mL)		
	Avg. Conc. (ng/mL)	Avg. Accuracy (%)	Precision (%RSD)	Avg. Conc. (ng/mL)	Avg. Accuracy (%)	Precision (%RSD)	Avg. Conc. (ng/mL)	Avg. Accuracy (%)	Precision (%RSD)
K <sub>1</sub>	0.246	98.2	5.75	0.763	102	3.39	8.23	103	4.36
MK4	0.249	100	4.55	0.729	97.2	3.83	8.11	101	5.17
MK7	0.245	98.0	3.88	0.715	95.3	4.98	8.23	103	4.76

### Human Plasma Analysis

Distinct endogenous vitamin K profiles were revealed by analyzing vitamin K<sub>1</sub> and vitamin K<sub>2</sub> in multiple plasma samples (Figure 4). The analysis of pooled human plasma is a good indication of the average vitamin K level of a collected population. The result showed clear detection of K<sub>1</sub> and MK4 in pooled plasma with no detectable MK7, reflecting that the collected population has very low consumption of the MK7-containing food. Conversely, a plasma sample from an individual of Japanese descent showed higher MK7 and no detectable MK4. The plasma sample from an individual taking vitamin K<sub>2</sub> medication showed strong K<sub>1</sub>, MK4, and MK7 signals. Overall, these results demonstrated that our established method can measure the variation of vitamin K<sub>1</sub> and K<sub>2</sub> concentrations seen in different populations and is thus applicable for a broad range of clinical samples.

**Figure 4:** Analysis of Endogenous Vitamin K<sub>1</sub> and K<sub>2</sub> in Human Plasma



## Conclusion

A method was developed for fast, accurate, and precise analysis of vitamin K<sub>1</sub> and K<sub>2</sub> (MK4 and MK7) in plasma. With quick and easy phospholipid removal sample preparation, 4 minutes of chromatographic analysis using a Raptor Biphenyl column, and sensitive MS/MS detection, the established method provides a reliable and high-throughput assay for clinical research and assessment of vitamin K<sub>1</sub> and K<sub>2</sub>.



## Raptor Biphenyl LC Columns (USP L11)

The innovative Biphenyl is Restek's most popular LC stationary phase because it is particularly adept at separating compounds that are hard to resolve or that elute early on C18 and other phenyl chemistries. As a result, the rugged Raptor Biphenyl column is extremely useful for fast separations in bioanalytical testing applications like drug and metabolite analyses, especially those that require a mass spectrometer (MS). Increasing retention of early-eluting compounds can limit ionization suppression, and the heightened selectivity helps eliminate the need for complex mobile phases that are not well suited for MS detection.

Length	2.1 mm cat.#	3.0 mm cat.#	4.6 mm cat.#
<b>1.8 µm Columns</b>			
30 mm	9309232	—	—
50 mm	9309252	930925E	—
100 mm	9309212	930921E	—
150 mm	9309262	—	—
<b>2.7 µm Columns</b>			
30 mm	9309A32	9309A3E	9309A35
50 mm	9309A52	9309A5E	9309A55
100 mm	9309A12	9309A1E	9309A15
150 mm	9309A62	9309A6E	9309A65
<b>5 µm Columns</b>			
30 mm	—	930953E	—
50 mm	9309552	930955E	9309555
100 mm	9309512	930951E	9309515
150 mm	9309562	930956E	9309565
250 mm	—	—	9309575

### Column Characteristics:

Stationary Phase Category: Phenyl (L11)  
 Ligand Type: Biphenyl  
 Particle: 1.8 µm, 2.7 µm, or 5 µm superficially porous silica (SPP or "core-shell")  
 Pore Size: 90 Å  
 Carbon Load: 7% (1.8 µm), 7% (2.7 µm), 5% (5 µm)  
 End-Cap: yes  
 Surface Area: 125 m<sup>2</sup>/g (1.8 µm), 130 m<sup>2</sup>/g (2.7 µm), or 100 m<sup>2</sup>/g (5 µm)

#### Recommended Usage:

pH Range: 2.0 to 8.0  
 Maximum Temperature: 80 °C  
 Maximum Pressure: 1,034 bar/15,000 psi\* (1.8 µm), 600 bar/8,700 psi (2.7 µm); 400 bar/5,800 psi (5 µm)


\* For maximum lifetime, recommended maximum pressure for 1.8 µm particles is 830 bar/12,000 psi.

#### Properties:

- Increased retention for dipolar, unsaturated, or conjugated solutes.
- Enhanced selectivity when used with methanolic mobile phase.
- Ideal for increasing sensitivity and selectivity in LC-MS analyses.

#### Switch to a Biphenyl when:

- Limited selectivity is observed on a C18.
- You need to increase retention of hydrophilic aromatics.





**Questions? Contact us or your local Restek representative ([www.restek.com/contact-us](http://www.restek.com/contact-us)).**

Restek patents and trademarks are the property of Restek Corporation. (See [www.restek.com/Patents-Trademarks](http://www.restek.com/Patents-Trademarks) for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit [www.restek.com/subscribe](http://www.restek.com/subscribe). To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2018 Restek Corporation. All rights reserved. Printed in the U.S.A.

**[www.restek.com](http://www.restek.com)**



**Lit. Cat.# CFAN2934-UNV**