

Optimizing Method Standard HJ 834-2017 for GC-MS Semivolatiles Analysis

Increase Speed and Certainty with Highly Inert RMX-5Sil MS Columns

By Erica Pack, Grace Idowu, Chris English, Ramkumar Dhandapani, Colton Myers

Key Highlights

- Next-generation TriMax GC column deactivation creates a robust and exceptionally neutral sample flow path.
- Maximum symmetry reduces peak tailing for a wide range of challenging semivolatiles—including acids, bases, and neutrals—improving identification and quantification.
- Optimized GC-MS conditions speed up sample throughput while still meeting method criteria.



Abstract

RMX-5Sil MS columns feature a groundbreaking deactivation technology that creates a highly inert surface that is broadly effective across a wide range of semivolatile compound classes. In this study, we evaluated the columns performance against the data quality criteria in Method Standard HJ 834-2017. Method criteria for system suitability, linearity, and recovery were easily met with optimized conditions that reduced sample analysis time by 20 minutes.

Introduction

Semivolatile compounds are monitored worldwide to protect human health and the environment. Analysis can be done with a variety of detectors, with GC-MS techniques being among the most popular. The People's Republic of China Method Standard HJ 834-2017 is used for GC-MS analysis of semivolatiles in soil and sediment samples. To ensure laboratories generate accurate results, Method Standard HJ 834-2017 specifies criteria for key parameters, including system suitability, linearity, and recovery. Since semivolatiles comprise a wide range of compound chemistries, including reactive compounds, it is essential that the sample flow path be highly inert. Activity in the flow path may contribute to poor peak shape and drifting retention times, which can make it difficult to meet data quality requirements.

In this study, we demonstrate how exceptionally inert RMX-5Sil MS columns can help laboratories meet data quality objectives for Method Standard HJ 834-2017 semivolatiles analysis. RMX-5Sil MS columns undergo a unique TriMax surface deactivation that eliminates active sites (e.g., silanols) and is broadly effective across compound classes. The RMX-5Sil MS column contains a traditional 5sil polymer, so it is a direct 5sil replacement, but the highly neutral surface improves peak shape for a wide variety of compounds, making it easier to meet data quality requirements. In this study, in addition to demonstrating the effectiveness of the column for GC-MS semivolatiles analysis, we also developed optimized instrument conditions that speed up run times and improve sample throughput.

Experimental

Standard and Sample Preparation

Multicomponent semivolatiles standards were prepared at 1, 5, 10, 20, and 50 ppm in dichloromethane for calibration. Internal standards and surrogate standards were present at 40 ppm in each calibration standard.

Instrument Conditions

GC-MS semivolatiles analysis was performed following the conditions in Method Standard HJ 834-2017. In addition, semivolatiles analysis was run under optimized conditions to allow faster analysis times. Table I summarizes the instrument conditions for both approaches.

Table I: Instrument Conditions for GC-MS Semivolatiles Analysis

Parameter	Method Standard HJ 834-2017 Reference Conditions	Optimized Conditions
Column	RMX-5Sil MS, 30 m x 0.25 mm ID x 0.25 µm (cat.# 17323)	RMX-5Sil MS, 30 m x 0.25 mm ID x 0.25 µm (cat.# 17323)
Injection	1 µL, splitless, 280 °C	1 µL, splitless, 280 °C
Carrier gas	Helium at 1 mL/min (constant flow)	Helium at 1 mL/min (constant flow)
Oven	35 °C (hold 2 min) to 150 °C at 15 °C/min (hold 5 min) to 290 °C at 3 °C/min (hold 2 min)	Faster oven conditions speed up analysis: 35 °C (hold 1 min) to 290 °C at 8 °C/min to 350 °C at 12 °C/min (hold 5 min)
Detector	El; 70 eV; 230 °C source; 150 °C quad; scan 35-450 m/z; solvent delay 5 min; 280 °C transfer line	El; 70 eV; 230 °C source; 150 °C quad; scan 35-450 m/z; solvent delay 2 min; 280 °C transfer line

Results and Discussion


System Suitability

The instrument was tuned with DFTPP in accordance with Method Standard HJ 834-2017, and then system suitability was assessed by evaluating the degradation of DDT to DDE and DDD. The requirement to not exceed 15% was easily met with DDT degradation to both DDE and DDD being <1%.

Chromatographic Performance

As shown in Figure 1, the RMX-5Sil MS column produced sharp, symmetrical peaks for 64 monitored semivolatiles, six surrogates, and six internal standards (76 compounds total) at 5 ppm under both the original Method Standard HJ 834-2017 conditions and the modified GC-MS conditions. While some coelutions occurred, these compounds could easily be differentiated by their ion m/z ratios. The standard method conditions produced good results, but the >60-minute run time limits the number of samples that can be analyzed. The improved method conditions resulted in a faster 40-minute run time, which is beneficial to high-throughput labs that have many samples to analyze under tight deadlines.


Figure 1: TIC Chromatogram of Method Standard HJ 834-2017 Conditions (>60 minutes)



Peaks	tR (min)	Peaks	tR (min)	Peaks	tR (min)
1. Bis(N-methoxy-N-methylamino)methane	3.989	26. p-Chloroaniline	9.630	51. Azobenzene	17.714
2. Phenol, 2-fluoro-	5.769	27. 1,3-Butadiene, 1,1,2,3,4,4-hexachloro-	9.767	52. Phenol, 2,4,6-tribromo-	17.972
3. Phenol-d6-	7.057	28. Phenol, 2-chloro-5-methyl-	10.486	53. Benzene, 1-bromo-4-phenoxy-	19.798
4. Phenol	7.069	29. Naphthalene, 1-methyl-	10.766	54. Benzene, hexachloro-	19.924
5. Bis(2-chloroethyl) ether	7.194	30. Hexachlorocyclopentadiene	11.140	55. Phenol, pentachloro-	21.347
6. Phenol, 2-chloro-	7.269	31. Phenol, 2,4,6-trichloro-	11.412	56. Phenanthrene-D10	22.411
7. Benzene, 1,3-dichloro-	7.496	32. Phenol, 2,4,5-trichloro-	11.486	57. 9H-Fluorene, 9-methylene-	22.561
8. 1,4-Dichlorobenzene-D4	7.579	33. 1,1'-Biphenyl, 2-fluoro-	11.640	58. Anthracene	22.911
9. Benzene, 1,4-dichloro-	7.602	34. Naphthalene, 1-chloro-	11.915	59. Carbazole	24.334
10. Benzene, 1,2-dichloro-	7.825	35. Dimethyl (2-nitroanilino)maleate	12.221	60. Dibutyl phthalate	28.107
11. Phenol, 2-methyl-	7.956	36. Dimethyl phthalate	12.926	61. Fluoranthene	31.814
12. Bis(2-chloro-1-methylethyl) ether	7.995	37. Benzene, 2-methyl-1,3-dinitro-	13.078	62. Pyrene	32.727
13. 1-Propanamine, N-nitroso-N-propyl-	8.183	38. Biphenylene	13.231	63. p-Terphenyl-d14	34.881
14. p-Cresol	8.183	39. m-Nitroaniline	13.680	64. Benzyl butyl phthalate	39.418
15. Ethane, hexachloro-	8.336	40. Acenaphthene-d10	13.813	65. Benz[a]anthracene	42.103
16. Nitrobenzene-D5	8.393	41. Acenaphthene	13.947	66. Chrysene-D12	42.103
17. Benzene, nitro-	8.421	42. Phenol, 2,4-dinitro-	14.149	67. Chrysene	42.298
18. Isophorone	8.795	43. Phenol, 4-nitro-	14.626	68. Bis(2-ethylhexyl) phthalate	44.571
19. Phenol, 2-nitro-	8.915	44. Dibenzofuran	14.763	69. Phthalic acid, hept-4-yl octyl ester	49.346
20. Phenol, 2,3-dimethyl-	9.002	45. Benzene, 1-methyl-2,4-dinitro-	14.820	70. Benzo[b]fluoranthene	49.728
21. Methane, bis(2-chloroethoxy)-	9.146	46. Diethyl phthalate	16.434	71. Benzo[k]fluoranthene	49.915
22. Phenol, 2,4-dichloro-	9.290	47. Fluorene	16.576	72. Benzo[a]pyrene	51.738
23. Benzene, 1,2,4-trichloro-	9.427	48. Benzene, 1-chloro-4-phenoxy-	16.802	73. Perylene - D12	52.132
24. Naphthalene-D8	9.508	49. Benzene, 1-chloro-3-phenoxy-	16.802	74. Indeno[1,2,3-cd]pyrene	58.454
25. Naphthalene	9.540	50. Phenol, 2-methyl, 4,6-dinitro-	17.043	75. Diben[a,h]anthracene	58.786
				76. Benzo[ghi]perylene	59.786

Column	RMX-5Sil MS, 30 m, 0.25 mm ID, 0.25 µm (cat.# 17323)		
Standard/Sample	1000 ppm HJ 834-2017 VOC and SVOCs mixture 155 (LGC) 1000 ppm HJ 834-2017 substitutes mixture 156 (LGC) 1000 ppm HJ 834-2017 internal standard mixture 174 (internal standard) Dichloromethane 5 µg/mL		
Diluent:			
Conc.:			
Injection			
Inj. Vol.:	1 µL splitless (hold 1 min)		
Liner:	Topaz 4.0 mm ID single taper inlet liner w/wool (cat.# 23303)		
Inj. Temp.:	280 °C		
Oven			
Oven Temp.:	35 °C (hold 2 min) to 150 °C at 15 °C/min (hold 5 min) to 290 °C at 3 °C/min (hold 2 min)		
Carrier Gas	He, constant flow		
Flow Rate:	1 mL/min @ 35 °C		
Detector	MS		
Mode:	Scan		
Scan Program:			
	Group 1	Start Time (min) 2.5	Scan Range (amu) 35-450
			Scan Rate (scans/sec)
Transfer Line Temp.:	280 °C		
Source Temp.:	230 °C		
Quad Temp.:	150 °C		
Electron Energy:	70 eV		
Tune Type:	PFTBA		
Ionization Mode:	EI		
Instrument	Agilent 7890B GC & 5977B MSD		
Sample Preparation	Reference standards were diluted to 5 ppm in dichloromethane.		

Figure 2: TIC Chromatogram of Optimized Method Conditions (40 minutes)



Peaks	t _r (min)	Peaks	t _r (min)	Peaks	t _r (min)
1. Bis(N-methoxy-N-methylamino)methane	3.941	27. 1,3-Butadiene, 1,1,2,3,4,4-hexachloro-	12.859	53. Benzene, 1-bromo-4-phenoxy-	20.596
2. Phenol, 2-fluoro-	6.189	28. Phenol, 2-chloro-5-methyl-	14.122	54. Benzene, hexachloro-	20.697
3. Phenol-d6-	8.226	29. Naphthalene, 1-methyl-	14.409	55. Phenol, pentachloro-	21.316
4. Phenol	8.252	30. Hexachlorocyclopentadiene	14.927	56. Phenanthrene-D10	21.822
5. Bis(2-chloroethyl) ether	8.417	31. Phenol, 2,4,6-trichloro-	15.284	57. 9H-Fluorene, 9-methylene-	21.88
6. Phenol, 2-chloro-	8.502	32. Phenol, 2,4,5-trichloro-	15.362	58. Anthracene	22.025
7. Benzene, 1,3-dichloro-	8.861	33. 1,1'-Biphenyl, 2-fluoro-	15.564	59. Carbazole	22.57
8. 1,4-Dichlorobenzene-D4	9.01	34. Naphthalene, 1-chloro-	15.806	60. Dibutyl phthalate	23.889
9. Benzene, 1,4-dichloro-	9.048	35. Dimethyl (2-nitroanilino)maleate	16.158	61. Fluoranthene	25.396
10. Benzene, 1,2-dichloro-	9.411	36. Dimethyl phthalate	16.837	62. Pyrene	26.035
11. Phenol, 2-methyl-	9.726	37. Benzene, 2-methyl-1,3-dinitro-	16.945	63. p-Terphenyl-d14	26.717
12. Bis(2-chloro-1-methylethyl) ether	9.784	38. Biphenylene	16.995	64. Benzyl butyl phthalate	28.309
13. 1-Propanamine, N-nitroso-N-propyl-	10.101	39. m-Nitroaniline	17.381	65. Benz[a]anthracene	29.638
14. p-Cresol	10.133	40. Acenaphthene-d10	17.45	66. Chrysene-D12	29.684
15. Ethane, hexachloro-	10.273	41. Acenaphthene	17.533	67. Chrysene	29.744
16. Nitrobenzene-D5	10.421	42. Phenol, 2,4-dinitro-	17.692	68. Bis(2-ethylhexyl) phthalate	30.159
17. Benzene, nitro-	10.468	43. Phenol, 4-nitro-	18.011	69. Phthalic acid, hept-4-yl octyl ester	31.991
18. Isophorone	11.137	44. Dibenzofuran	18.05	70. Benzo[b]fluoranthene	32.608
19. Phenol, 2-nitro-	11.322	45. Benzene, 1-methyl-2,4-dinitro-	18.104	71. Benzo[k]fluoranthene	32.683
20. Phenol, 2,3-dimethyl-	11.561	46. Diethyl phthalate	18.972	72. Benzo[a]pyrene	33.417
21. Methane, bis(2-chloroethoxy)-	11.822	47. Fluorene	19.061	73. Perylene - D12	33.583
22. Phenol, 2,4-dichloro-	12.006	48. Benzene, 1-chloro-4-phenoxy-	19.168	74. Indeno[1,2,3-cd]pyrene	35.743
23. Benzene, 1,2,4-trichloro-	12.23	49. Benzene, 1-chloro-3-phenoxy-	19.175	75. Diben[a,h]anthracene	35.798
24. Naphthalene-D8	12.36	50. Phenol, 2-methyl, 4,6-dinitro-	19.301	76. Benzo[ghi]perylene	36.197
25. Naphthalene	12.411	51. Azobenzene	19.626		
26. p-Chloroaniline	12.622	52. Phenol, 2,4,6-tribromo-	19.78		

Column RMX-5SIL MS, 30 m, 0.25 mm ID, 0.25 µm (cat.# 17323)

Standard/Sample 1000 ppm HJ 834-2017 VOC and SVOCs mixture 155 (LGC)

1000 ppm HJ 834-2017 substitutes mixture 156 (LGC)

1000 ppm HJ 834-2017 internal standard mixture 174 (internal standard)

Diluent: Dichloromethane

Conc.: 5 µg/mL

Injection

Inj. Vol.: 1 µL splitless (hold 1 min)

Liner: Topaz 4.0 mm ID single taper inlet liner w/wool (cat.# 23303)

Inj. Temp.: 280 °C

Oven

Oven Temp.: 35 °C (hold 1 min) to 290 °C at 8 °C/min to 350 °C at 12 °C/min (hold 5 min)

Carrier Gas

Flow Rate: He, constant flow

1 mL/min @ 35 °C

Detector

Mode: MS

Scan Program: Scan

Group 1

Start Time (min) 2.5

Scan Range (amu) 35-450

Scan Rate (scans/sec)

Transfer Line Temp.: 280 °C

Source Temp.: 230 °C

Quad Temp.: 150 °C

Electron Energy: 70 eV

Tune Type: PFTBA

Ionization Mode: EI

Instrument Agilent 7890B GC & 5977B MSD

Sample Preparation Reference standards were diluted to 5 ppm in dichloromethane.

Want to customize a method?

Enter your target analytes into Restek's free Pro EZGC chromatogram modeling software and instantly generate optimized method conditions for your specific compound list.

Try it now!



Quantitative Performance

Good peak shape improves signal-to-noise ratios, which allows reliable integration at low concentrations and minimizes the need for time-consuming user intervention. Accurate integration across the calibration range increases confidence that method performance criteria will be met, even for difficult reactive compounds. Method performance was assessed by evaluating linearity (unweighted, $R > 0.990$); relative response factor (RSD $< 30\%$); and recovery of a midpoint standard (70-130% recovery at 10 ppm).

Limits of detection and quantitation were determined using generally accepted methods where $LOD = 3.3 \times$ standard deviation of the lowest calibration point/calibration slope, and $LOQ = 10 \times$ standard deviation of the lowest calibration point/calibration slope. Values were compared to the LOD and LOQ values described in the method. The method requirements show LOD and LOQ as sample concentrations. Since matrix was not used in this study, sample concentrations were converted to extract concentrations based on a theoretical sample mass of 20 g, and final extract volume of 1 mL (Table II). Without matrix, the differences between column performance can be separated from the efficacy of the sample preparation.

For quantitative analysis, calibration was performed using the optimized instrument conditions. Summarized results are presented in Table II and the values for individual semivolatiles are given in Table III. All compounds (excluding internal standards and surrogates) passed linearity criteria of RRF RSD $< 30\%$ and $R > 0.990$. RSD values ranged 2.1-22.5% and R ranged 0.999-1.00. While not required for HJ 834, R^2 was also evaluated because it is often employed in other methods. In this study R^2 was > 0.995 for each compound, ranging 0.997-1.0.

Table II: Results Summary for Method Standard HJ 834-2017 GC-MS Semivolatiles Analysis (Optimized Conditions)

	RRF %RSD ($< 30\%$)	R^2 (> 0.990)	R (> 0.990)	LOD (Method HJ 834-2017)	LOQ (Method HJ 834-2017)	%Recovery
Acceptable	64	64	64	64	64	64
Unacceptable	0	0	0	0	0	0
Maximum	23%	1.000	1.000	0.89	2.98	109%
Minimum	2%	0.997	0.999	0.02	0.05	75%

* R^2 is not specified in Method Standard HJ 834-2017. It is included because it is a commonly used linearity metric in other methods for GC-MS semivolatiles analysis.

Table III: Individual Results by Compound for Method Standard HJ 834-2017 GC-MS Semivolatiles Analysis (Optimized Conditions)

Compound	RRF %RSD (<30%)	R ²	R	Method HJ 834-2017 LOD (ppm)	Method HJ 834-2017 LOQ (ppm)	LOD (ppm)	LOQ (ppm)	%Recovery of 10 ppm
Bis(N-methoxy-N-methylamino)methane	2%	1.000	1.000	1.60	6.40	0.36	1.22	105%
Phenol	4%	0.998	0.999	2.00	8.00	0.33	1.11	95%
Bis(2-chloroethyl) ether	3%	0.999	1.000	2.00	8.00	0.17	0.58	99%
Phenol, 2-chloro-	4%	0.999	0.999	2.00	8.00	0.16	0.55	96%
Benzene, 1,3-dichloro-	5%	0.999	0.999	1.80	7.20	0.20	0.68	100%
Benzene, 1,4-dichloro-	5%	1.000	1.000	1.20	4.80	0.19	0.62	97%
Benzene, 1,2-dichloro-	7%	1.000	1.000	1.60	6.40	0.13	0.45	96%
Phenol, 2-methyl-	5%	0.998	0.999	1.60	6.40	0.25	0.82	101%
Bis(2-chloro-1-methylethyl) ether	6%	0.999	0.999	1.60	6.40	0.34	1.13	104%
1-Propanamine, N-nitroso-N-propyl-	3%	1.000	1.000	2.00	8.00	0.43	1.43	109%
p-Cresol	6%	1.000	1.000	2.00	8.00	0.11	0.38	100%
Ethane, hexachloro-	10%	1.000	1.000	2.00	8.00	0.53	1.77	97%
Benzene, nitro-	7%	0.999	1.000	1.40	5.60	0.23	0.76	102%
Isophorone	6%	1.000	1.000	2.00	8.00	0.09	0.30	99%
Phenol, 2-nitro-	10%	1.000	1.000	2.00	8.00	0.15	0.50	96%
Phenol, 2,3-dimethyl-	6%	1.000	1.000	1.80	7.20	0.25	0.84	99%
Methane, bis(2-chloroethoxy)-	5%	1.000	1.000	1.40	5.60	0.09	0.31	100%
Phenol, 2,4-dichloro-	13%	1.000	1.000	4.00	16.00	0.09	0.30	98%
Benzene, 1,2,4-trichloro-	11%	1.000	1.000	1.80	7.20	0.16	0.54	98%
Naphthalene	9%	0.998	0.999	1.60	6.40	0.16	0.54	99%
p-Chloroaniline	8%	0.999	0.999	1.40	5.60	0.17	0.57	94%
1,3-Butadiene, 1,1,2,3,4,4-hexachloro-	5%	1.000	1.000	1.40	5.60	0.29	0.98	97%
Phenol, 2-chloro-5-methyl-	6%	1.000	1.000	1.80	7.20	0.23	0.77	93%
Naphthalene, 1-methyl-	10%	0.999	0.999	1.80	7.20	0.07	0.22	98%
Hexachlorocyclopentadiene	22%	0.998	0.999	1.20	4.80	0.27	0.92	89%
Phenol, 2,4,6-trichloro-	21%	0.999	1.000	1.20	4.80	0.13	0.44	92%
Phenol, 2,4,5-trichloro-	21%	0.999	0.999	1.60	6.40	0.27	0.91	92%
Naphthalene, 1-chloro-	15%	0.999	1.000	2.00	8.00	0.10	0.32	95%
Dimethyl (2-nitroanilino)maleate	21%	0.999	1.000	2.00	8.00	0.62	2.05	89%
Dimethyl phthalate	15%	1.000	1.000	2.00	8.00	0.12	0.39	93%
Benzene, 2-methyl-1,3-dinitro-	20%	0.998	0.999	2.00	8.00	0.56	1.87	90%
Biphenylene	18%	1.000	1.000	2.00	8.00	0.10	0.32	99%
m-Nitroaniline	22%	0.999	1.000	1.60	6.40	0.89	2.98	90%
Acenaphthene	10%	0.998	0.999	1.80	7.20	0.02	0.08	99%
Phenol, 2,4-dinitro-	16%	0.997	0.999	1.40	5.60	0.69	2.32	75%
Phenol, 4-nitro-	13%	1.000	1.000	1.60	7.20	0.23	0.75	85%
Dibenzofuran	14%	1.000	1.000	2.00	8.00	0.09	0.29	99%
Benzene, 1-methyl-2,4-dinitro-	23%	1.000	1.000	2.00	8.00	0.13	0.42	91%
Diethyl Phthalate	13%	1.000	1.000	2.00	8.00	0.05	0.18	96%
Fluorene	15%	0.999	1.000	1.80	7.20	0.04	0.14	98%
Benzene, 1-chloro-4-phenoxy-	8%	1.000	1.000	1.80	7.20	0.12	0.39	103%
Benzene, 1-chloro-3-phenoxy-	13%	1.000	1.000	4.00	16.00	0.33	1.11	95%
Phenol, 2-methyl, 4,6-dinitro-	22%	0.998	0.999	1.60	6.40	0.39	1.30	103%
Azobenzene	6%	0.999	0.999	6.00	24.00	0.15	0.50	104%

Table III: Individual Results by Compound for Method Standard HJ 834-2017 GC-MS Semivolatiles Analysis (Optimized Conditions)

Compound	RRF %RSD (<30%)	R ² *	R	Method HJ 834-2017 LOD (ppm)	Method HJ 834-2017 LOQ (ppm)	LOD (ppm)	LOQ (ppm)	%Recovery of 10 ppm
Benzene, 1-bromo-4-phenoxy-	13%	1.000	1.000	2.00	8.00	0.20	0.67	95%
Benzene, hexachloro-	8%	1.000	1.000	2.00	8.00	0.02	0.07	99%
Phenol, pentachloro-	11%	0.998	0.999	2.00	8.00	0.39	1.30	89%
9H-Fluorene, 9-methylene-	8%	0.998	0.999	2.00	8.00	0.03	0.11	101%
Anthracene	13%	1.000	1.000	4.00	16.00	0.03	0.09	104%
Carbazole	11%	0.998	0.999	2.00	8.00	0.07	0.24	96%
Dibutyl phthalate	12%	0.999	0.999	2.00	8.00	0.03	0.10	99%
Fluoranthene	10%	1.000	1.000	4.00	16.00	0.07	0.22	100%
Pyrene	9%	1.000	1.000	2.00	8.00	0.03	0.11	101%
Benzyl butyl phthalate	12%	0.998	0.999	2.00	8.00	0.20	0.68	94%
Benz[a]anthracene	5%	0.999	0.999	2.00	8.00	0.08	0.26	95%
Triphenylene	7%	0.999	1.000	2.00	8.00	0.13	0.42	106%
Bis(2-ethylhexyl) phthalate	15%	1.000	1.000	4.00	16.00	0.09	0.29	102%
Phthalic acid, hept-4-yl octyl ester	18%	1.000	1.000	2.00	8.00	0.02	0.05	103%
Benzo[b]fluoranthene	17%	0.999	1.000	2.00	8.00	0.09	0.31	96%
Benzo[k]fluoranthene	16%	1.000	1.000	4.00	16.00	0.09	0.32	102%
Benzo[a]pyrene	15%	0.999	0.999	2.00	8.00	0.08	0.26	104%
Indeno[1,2,3-cd]pyrene	11%	0.999	1.000	2.00	8.00	0.05	0.18	96%
Diben[a,h]anthracene	5%	0.999	0.999	2.00	8.00	0.10	0.35	107%
Benzo[ghi]perylene	7%	1.000	1.000	4.00	16.00	0.04	0.15	106%

*R² is not specified in Method Standard HJ 834-2017, but it is included here because it is a statistical metric that is commonly used to demonstrate linearity in other methods for GC-MS semivolatiles analysis.

Conclusion

This study demonstrated that the highly inert RMX-5Sil MS column produced method-compliant results for a wide range of challenging semivolatiles. The inertness of the column surface produced sharp, symmetrical peaks that simplified integration and resulted in excellent results for system suitability; linearity (RRF %RSD and R²); LOD; LOQ; and recovery. In addition to meeting Method Standard HJ 834-2017 requirements and producing highly linear calibration curves, the optimized GC-MS conditions that were developed here reduced sample analysis time from 60 to 40 minutes, allowing high-throughput laboratories to improve sample throughput.

Need support with applications or product selection?

Contact Technical Service!



Want help with quoting, purchasing, and more?

Contact Us!



Featured Products

RMX-5Sil MS GC Capillary Column

Catalog No.	Product Name	Units
17323	RMX-5Sil MS GC Capillary Column, 30 m, 0.25 mm ID, 0.25 µm	ea.



Topaz Single Taper Inlet Liner

Catalog No.	Product Name	Units
23303	Topaz, Single Taper Inlet Liner, 4.0 mm x 6.5 x 78.5, for Agilent GCs, w/Quartz Wool, Premium Deactivation	5-pk.



Restek Electronic Leak Detector

Catalog No.	Product Name	Units
28500	Restek Electronic Leak Detector (includes carrying case; universal AC power adaptor [U.S., UK, Europe, Australia, Japan]; 6-ft USB charging cable)	ea.



Benzo[a]pyrene

1000 µg/mL, Acetone, 1 mL/ampul

Catalog No.	Contains	Units
31271	Benzo[a]pyrene (50-32-8)	ea.

Dibenz[a,h]anthracene

1000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31276	Dibenz[a,h]anthracene (53-70-3)	ea.

Indeno[1,2,3-cd]pyrene

1000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31279	Indeno[1,2,3-cd]pyrene (193-39-5)	ea.

Benzo[g,h,i]perylene

1000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31273	Benzo[g,h,i]perylene (191-24-2)	ea.




What Makes RMX Columns Better?

Highly Effective TriMax Deactivation Protects Analytes From Surface Interactions, Improving Peak Shape and Sensitivity for a Wide Range of Compound Chemistries



TriMax Deactivation




Inactive



Active


Non-TriMax Deactivation



Inactive



Active



- Inactive compounds: *alkanes, alkenes, alkynes, etc.*
- Active compounds: *acids, bases, alcohols, esters, ethers, etc.*
- Residual active site

Learn More About
RMX Columns Today!



For information on Restek patents and trademarks, visit www.restek.com/patents-trademarks. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2026 Restek Corporation. All rights reserved.

www.restek.com



Lit. Cat.# EVAN5233-UNV