


GC-MS半揮発性有機化合物分析：微量分析の限界に挑む

低濃度SVOC分析におけるGC-MSの感度とピーク形状の課題

By Jessi Collier, Erica Pack, Ramkumar Dhandapani, and Colton Myers

- 次世代TriMax不活性化処理により、より低い検出限界でも品質基準を満たすデータが取得可能に
- 高い不活性度により、分析が難しい成分を含む多様なSVOCにおいてシャープで対称性の高いピーク形状を実現
- 高感度分析 (SIM) と汎用的なスクリーニング (scan) の両方に適したGC-MS分析条件を確立



SVOCのGC-MS分析で検出限界の改善が求められる背景

半揮発性有機化合物の分析を実施するラボでは、顧客要求や溶媒削減（規制対応・コスト削減）のため、検出限界の引き下げが常に求められています。MSベースのメソッド（例：EPA 8270）は一般的なアプローチであり、MSを活用することでより低い検出限界を実現できます。

EPA Method 8270に代表されるMSを用いたSVOCの分析手法は広く活用されており、ラボは目的に応じて最適な検出限界を達成できる装置を柔軟に選択できます。例えば、GC-MS/MSのような高感度な装置を採用する場合もある一方で、GC-MSを用いてSIMによる感度向上を図る運用や、より汎用的なscanモードによる分析を選択する場合もあります。

SVOC分析におけるGC-MSのscan法とSIM法の使い分け

SIMとscanはどちらもGC-MSで実施可能ですが、scanが幅広い質量範囲の信号をまとめて取得するシンプルで汎用的な手法であるのに対し、SIMは特定イオンを特定の保持時間に絞って測定することができるため、より高い感度が得られます。SIMはメソッドの設定や維持が煩雑になる傾向がありますが、検出限界の引き下げを検討する場合、新たな装置投資を行わずに感度を高める選択肢として、GC-MSのSIM運用が有力候補になります。

低濃度SVOC分析を想定したGC-MS測定ワークフロー

本研究では、シンプルで堅牢な scanと、GC-MS/MSに代わり得る高感度なSIMという二つのワークフローを目的に応じて選択できるよう、scanとSIMの両モードにおけるGC-MS半揮発性有機化合物分析の条件を確立しました (Figure 1 & Figure 2)。

RMX-5Sil MSカラムによるSVOC分析のクロマトグラフィー性能評価

いずれのメソッドもRMX-5Sil MSカラムを採用しました。このカラムは標準的な5silの選択性を備えながらも、次世代のTriMax不活性化処理により、極めて高い不活性度を実現することが示されています[1,2]。Figure 2に示すように、ベンジジン、ペンタクロロフェノール、2,4-ジニトロフェノールといった主要化合物において、ピーク対称性は一貫して品質目標 (≤ 2) を満たしました。さらに、分離が難しいPAH類については、EPA 8270のシステム適合性で求められるvalley% 基準 ($> 50\%$) に対し、benzo[b]fluorantheneとbenzo[k]fluorantheneのピーク間valley%は94%、indeno[1,2,3-cd]pyreneとdibenz[a,h]anthraceneのvalley%は95%と、極めて良好な結果が得られました。

pgレベルにおけるSVOCの感度とピーク形状の評価結果

また、本稿で示した条件下で RMX-5Sil MS カラムを使用したところ、scanモードでは1 ng on-column、SIMモードでは200 pg on-columnという極めて低いオンカラム量においても、幅広い半揮発性有機化合物に対して良好なピーク形状と十分な分離性能が維持され、MSの検出限界近傍においても感度およびクロマトグラフィー性能の両立が可能であることが示されました。

メソッドのカスタマイズが必要ですか？

EZGC Chromatogram modelerに目的化合物を入力するだけで、最適化されたメソッド条件をすぐに作成できます。


Try it now!



GC-MSによる低濃度SVOC 165成分分析の実測データ

GC-MSによる低濃度SVOC 165成分分析におけるRMX-5Si1 MSカラムのクロマトグラフィー性能評価のクロマトグラムを以下に示します。

Figure 1: scan Mode (1 ng on-column)



Column RMX-5Si1 MS, 30 m, 0.25 mm ID, 0.25 μ m (cat.# [17323](#))
Standard/Sample

Methapyrilene (cat.# [32460](#))
Appendix IX mix #1, revised (cat.# [32459](#))
SVOC additions (cat.# [31909](#))
Benzoic acid (cat.# [31879](#))
8270 MegaMix (cat.# [31850](#))
Appendix IX mix #2 (cat.# [31806](#))
Acid surrogate mix (4/89 SOW) (cat.# [31025](#))
Base neutral surrogate mix (4/89 SOW) (cat.# [31024](#))
Revised SV internal standard mix (cat.# [31886](#))

Diluent: Methylene chloride
Conc.: 1 μ g/mL

Injection

Inj. Vol.: 1 μ L split (split ratio 10:1)
Liner: Topaz 4.0 mm ID single taper liner w/wool (cat.# [23303](#))
Inj. Temp.: 250 $^{\circ}$ C
Split Vent Flow Rate: 12 mL/min

Oven

Oven Temp.: 40 $^{\circ}$ C (hold 1 min) to 280 $^{\circ}$ C at 12.4 $^{\circ}$ C/min to 330 $^{\circ}$ C at 3.3 $^{\circ}$ C/min (hold 3.65 min)

Carrier Gas

He, constant flow

Flow Rate: 1.2 mL/min

Linear Velocity: 39.7 cm/sec

Detector

MS

Mode: Scan

Scan Program:

Group	Start Time (min)	Scan Range (amu)	Scan Rate (scans/sec)
1	2	35-550	5

Transfer Line Temp.: 280 $^{\circ}$ C
Analyzer Type: Quadrupole
Source Temp.: 330 $^{\circ}$ C
Quad Temp.: 180 $^{\circ}$ C
Solvent Delay Time: 2 min
Tune Type: PFTBA
Ionization Mode: EI

Instrument

Agilent 7890B GC & 5977B MSD



Notes To simplify ordering, the SVOC MegaMix 150 kit (cat.# [31907](#)) contains one ampul each of the following standards that were used in this experiment.

- 8270 MegaMix (cat.# [31850](#))
- SVOC Additions standard (cat.# [31909](#))
- Appendix IX mix #1, revised (cat.# [32459](#))
- Methapyrilene (cat.# [32460](#))
- Appendix IX mix #2 (cat.# [31806](#))
- Benzoic acid (cat.# [31879](#))

Figure 1: (続き)

Peaks	t _R (min)	Conc. (µg/mL)	Peaks	t _R (min)	Conc. (µg/mL)
1. 1,4-Dioxane-d8	2.400	1	84. 3-Nitroaniline	11.595	1
2. 1,4-Dioxane	2.428	1	85. Acenaphthene-d10	11.647	1
3. N-Nitrosodimethylamine	2.702	1	86. Acenaphthene	11.704	1
4. Pyridine	2.774	1	87. 2,4-Dinitrophenol	11.809	1
5. Ethyl methacrylate	3.175	1	88. 4-Nitrophenol	11.996	1
6. 2-Picoline	3.554	1	89. Pentachlorobenzene	11.996	1
7. N-Nitrosomethylethylamine	3.648	1	90. Dibenzofuran	12.040	1
8. Methyl methanesulfonate	4.021	1	91. 2,4-Dinitrotoluene	12.070	1
9. 2-Fluorophenol	4.258	1	92. 1-Naphthylamine (1-aminonaphthalene)	12.184	1
10. Acrylamide	4.543	1	93. 2,3,5,6-Tetrachlorophenol	12.236	1
11. N-Nitrosodiethylamine	4.875	1	94. 2,3,4,6-Tetrachlorophenol	12.319	1
12. Ethyl methanesulfonate	4.977	1	95. 2-Naphthylamine (2-aminonaphthalene)	12.341	1
13. Benzaldehyde	5.454	1	96. Diethylphthalate	12.602	1
14. Phenol-d6	5.612	1	97. n-Hexadecane (C16)	12.699	1
15. Phenol	5.629	1	98. Fluorene	12.699	1
16. Aniline	5.638	1	99. Zalophus (thionazine)	12.745	1
17. Pentachloroethane	5.705	1	100. 4-Chlorophenyl phenyl ether	12.745	1
18. Bis(2-chloroethyl)ether	5.743	1	101. 4-Nitroaniline	12.745	1
19. 2-Chlorophenol	5.814	1	102. 5-Nitro-o-tolidine	12.745	1
20. n-Decane (C10)	5.962	1	103. 4,6-Dinitro-2-methylphenol (Dinitro-o-cresol)	12.853	1
21. 1,3-Dichlorobenzene	6.066	1	104. Diphenylamine†	12.979	1
22. 1,4-Dichlorobenzene-d4	6.153	1	105. Azobenzene*	13.042	1
23. 1,4-Dichlorobenzene	6.180	1	106. 2,4,6-Tribromophenol	13.162	1
24. Benzyl alcohol	6.393	1	107. Sulfotepp	13.379	1
25. 1,2-Dichlorobenzene	6.426	1	108. 1,3,5-Trinitrobenzene	13.579	1
26. 2-Methylphenol (o-cresol)	6.612	1	109. Phorate	13.579	1
27. 2,2'-Oxybis(1-chloropropane)	6.634	1	110. Phenacetin	13.625	1
28. N-Nitrosopyrrolidine	6.765	1	111. 4-Bromophenyl phenyl ether	13.677	1
29. Acetophenone	6.825	1	112. Diallylate	13.734	1
30. N-Nitrosomorpholine	6.836	1	113. Hexachlorobenzene	13.785	1
31. N-Nitroso-di-n-propylamine	6.847	1	114. Dimethoate	13.865	1
32. 3-Methylphenol (m-cresol)	6.885	1	115. Atrazine	14.071	1
33. 4-Methylphenol (p-cresol)	6.885	1	116. 4-Aminobiphenyl	14.168	1
34. Indene	6.885	1	117. Pentachlorophenol	14.170	1
35. o-Tolidine	6.890	1	118. Pentachloronitrobenzene (Quintozene)	14.214	1
36. Hexachloroethane	6.989	1	119. Propyzamide	14.380	1
37. Nitrobenzene-d5	7.067	1	120. Phenanthrene-d10	14.477	1
38. Nitrobenzene	7.101	1	121. n-Octadecane (C18)	14.517	1
39. N-Nitrosopiperidine	7.347	1	122. Phenanthrene	14.517	1
40. Isophorone	7.536	1	123. Dinoseb	14.599	1
41. 2-Nitrophenol	7.673	1	124. Disulfoton	14.599	1
42. 2,4-Dimethylphenol	7.816	1	125. Anthracene	14.614	1
43. Benzoic acid	7.943	1	126. Carbazole	14.963	1
44. O,O,O-Triethyl phosphorothioate	7.943	1	127. Methyl parathion	15.295	1
45. Bis(2-chloroethoxy)methane	7.970	1	128. Di-n-butyl phthalate	15.775	1
46. 2,4-Dichlorophenol	8.113	1	129. 4-Nitroquinoline-N-oxide	16.027	1
47. α,α-Dimethylphenethylamine (phentermine)	8.268	1	130. Parathion (ethyl parathion)	16.078	1
48. 1,2,4-Trichlorobenzene	8.268	1	131. Methapyrilene hydrochloride	16.273	1
49. Naphthalene-d8	8.342	1	132. Isodrin	16.518	1
50. Naphthalene	8.376	1	133. Fluoranthene	16.804	1
51. α-Terpineol	8.451	1	134. Benzidine	17.127	1
52. 4-Chloroaniline	8.519	1	135. Pyrene	17.222	1
53. 2,6-Dichlorophenol	8.519	1	136. p-Terphenyl-d14	17.622	1
54. Hexachloropropene	8.571	1	137. Aramite isomer 2	17.668	1
55. Hexachlorobutadiene	8.674	1	138. Aramite isomer 1	17.822	1
56. Quinoline	8.988	1	139. p-Dimethylaminoazobenzene	17.908	1
57. ε-Caprolactam	9.085	1	140. Chlorobenzilate	18.034	1
58. 1,4-Phenylenediamine	9.160	1	141. Famphur	18.485	1
59. N-Nitrosodibutylamine	9.195	1	142. 3,3'-Dimethylbenzidine (o-tolidine)	18.548	1
60. 4-Chloro-3-methylphenol	9.480	1	143. Kepone	18.548	1
61. Isosafrole isomer 2	9.560	1	144. Benzyl butyl phthalate	18.640	1
62. 2-Methylnaphthalene	9.680	1	145. Bis(2-ethylhexyl)adipate	18.891	1
63. 1-Methylnaphthalene	9.857	1	146. 2-Acetylaminofluorene	19.017	1
64. 1,2,4,5-Tetrachlorobenzene	10.023	1	147. 4,4'-Methylene-bis(2-chloroaniline)	19.108	1
65. Hexachlorocyclopentadiene	10.023	1	148. 3,3'-Dichlorobenzidine	19.560	1
66. 2,3-Dichloroaniline	10.217	1	149. Benz[a]anthracene	19.560	1
67. Isosafrole isomer 1	10.217	1	150. Chrysene-d12	19.560	1
68. 2,4,6-Trichlorophenol	10.248	1	151. Chrysene	19.612	1
69. 2,4,5-Trichlorophenol	10.303	1	152. Bis(2-ethylhexyl)phthalate	19.823	1
70. 2-Fluorobiphenyl	10.412	1	153. 6-Methylchrysene	20.475	1
71. Safrole	10.547	1	154. Di-n-octyl phthalate	21.030	1
72. Biphenyl	10.585	1	155. Benzo[b]fluoranthene	21.590	1
73. 2-Chloronaphthalene	10.585	1	156. 7,12-Dimethylbenz[a]anthracene	21.590	1
74. 1-Chloronaphthalene	10.623	1	157. Benzo[k]fluoranthene	21.641	1
75. 2-Nitroaniline	10.801	1	158. Benzo[a]pyrene	22.248	1
76. Diphenyl ether	10.801	1	159. Perylene-d12	22.379	1
77. 1,4-Naphthoquinone	10.926	1	160. 3-Methylcholanthrene	23.082	1
78. 1,4-Dinitrobenzene	11.075	1	161. Dibenz(a,h)acridine	24.357	1
79. Dimethylphthalate	11.224	1	162. Dibenz[a,j]acridine	24.483	1
80. 1,3-Dinitrobenzene	11.304	1	163. Indeno[1,2,3-cd]pyrene	24.883	1
81. 2,6-Dinitrotoluene	11.304	1	164. Dibenz[a,h]anthracene	24.969	1
82. 1,2-Dinitrobenzene	11.361	1	165. Benzo[g,h,i]perylene	25.541	1
83. Acenaphthylene	11.361	1			

Figure 2: SIM Mode (200 pg on-column)



Column RMX-5Sil MS, 30 m, 0.25 mm ID, 0.25 µm (cat.# [17323](#))

Standard/Sample

Methapyrene (cat.# [32460](#))
 Appendix IX mix #1, revised (cat.# [32459](#))
 SVOC additions (cat.# [31909](#))
 Benzoic acid (cat.# [31879](#))
 8270 MegaMix (cat.# [31850](#))
 Appendix IX mix #2 (cat.# [31806](#))
 Acid surrogate mix (4/89 SOW) (cat.# [31025](#))
 Base neutral surrogate mix (4/89 SOW) (cat.# [31024](#))
 Revised SV internal standard mix (cat.# [31886](#))

Diluent:

Conc.:

Injection

Inj. Vol.:

Liner:

Inj. Temp.:

Split Vent Flow Rate:

Oven

Oven Temp.:

Carrier Gas

Flow Rate:

Linear Velocity:

Detector

Mode:

Methylene chloride

200 ng/mL

1 µL split (split ratio 10:1)

Topaz 4.0 mm ID single taper liner w/wool (cat.# [23303](#))

250 °C

12 mL/min

40 °C (hold 1 min) to 280 °C at 12.4 °C/min to 330 °C at 3.3 °C/min (hold 3.65 min)

He, constant flow

1.2 mL/min

39.7 cm/sec

MS

SIM

SIM Program:

Group	Start Time (min)	Ion(s) (m/z)	Dwell (ms)
1	2	see peak list	300
2	5.2	see peak list	286
3	7.2	see peak list	289
4	9.35	see peak list	285
5	10.7	see peak list	273
6	12.45	see peak list	290
7	13.29	see peak list	280
8	14.8	see peak list	300
9	17.4	see peak list	297
10	19.3	see peak list	285

Transfer Line Temp.:

Analyzer Type:

Source Temp.:

Quad Temp.:

Solvent Delay Time:

Tune Type:

Ionization Mode:

Instrument

Notes

280 °C

Quadrupole

330 °C

180 °C

2 min

PFTBA

EI

Agilent 7890B GC & 5977B MSD

To simplify ordering, the SVOC MegaMix 150 kit (cat.# [31907](#)) contains one ampul each of the following standards that were used in this experiment.

- 8270 MegaMix (cat.# [31850](#))
- SVOC Additions standard (cat.# [31909](#))
- Appendix IX mix #1, revised (cat.# [32459](#))
- Methapyrene (cat.# [32460](#))
- Appendix IX mix #2 (cat.# [31806](#))
- Benzoic acid (cat.# [31879](#))

Figure 2: (続き)

Peaks	t _r (min)	Conc. (ng/mL)	SIM Ion	Dwell Time (ms)	Ion Group	Peaks	t _r (min)	Conc. (ng/mL)	SIM Ion	Dwell Time (ms)	Ion Group
1. 1,4-Dioxane-d8	2.401	200	96	25	1	84. 3-Nitroaniline	11.591	200	65	13	5
2. 1,4-Dioxane	2.427	200	88	25	1	85. Acenaphthene-d10	11.649	200	162	13	5
3. N-Nitrosodimethylamine	2.702	200	74	25	1	86. Acenaphthene	11.706	200	153	13	5
4. Pyridine	2.774	200	79	25	1	87. 2,4-Dinitrophenol	11.809	200	184	13	5
5. Ethyl methacrylate	3.175	200	69	25	1	88. 4-Nitrophenol	11.996	200	139	13	5
6. 2-Picoline	3.554	200	93	25	1	89. Pentachlorobenzene	11.996	200	250	13	5
7. N-Nitrosomethylethylamine	3.648	200	42	25	1	90. Dibenzofuran	12.043	200	84	13	5
8. Methyl methanesulfonate	4.021	200	80	25	1	91. 2,4-Dinitrotoluene	12.183	200	89	13	5
9. 2-Fluorophenol	4.258	200	112	25	1	92. 1-Naphthylamine (1-aminonaphthalene)	12.235	200	143	13	5
10. Acrylamide	4.543	200	71	25	1	93. 2,3,5,6-Tetrachlorophenol	12.339	200	232	13	5
11. N-Nitrosodiethylamine	4.977	200	102	25	1	94. 2,3,4,6-Tetrachlorophenol	12.339	200	230	13	5
12. Ethyl methanesulfonate	4.977	200	109	25	1	95. 2-Naphthylamine (2-aminonaphthalene)	12.339	200	115	13	5
13. Benzaldehyde	5.454	200	77	11	2	96. Diethylphthalate	12.596	200	149	13	5
14. Phenol-d6	5.612	200	99	11	2	97. n-Hexadecane (C16)	12.703	200	57	29	6
15. Phenol	5.633	200	94	11	2	98. Fluorene	12.703	200	166	29	6
16. Aniline	5.633	200	93	11	2	99. Zalophus (thionazine)	12.746	200	107	29	6
17. Pentachloroethane	5.7	200	167	11	2	100. 4-Chlorophenyl phenyl ether	12.746	200	204	29	6
18. Bis(2-chloroethyl)ether	5.743	200	63	11	2	101. 4-Nitroaniline	12.746	200	65	29	6
19. 2-Chlorophenol	5.814	200	128	11	2	102. 5-Nitro-o-toluidine	12.746	200	152	29	6
20. n-Decane (C10)	5.962	200	57	11	2	103. 4,6-Dinitro-2-methylphenol (Dinitro-o-cresol)	12.852	200	198	29	6
21. 1,3-Dichlorobenzene	6.066	200	146	11	2	104. Diphenylamine†	12.98	200	169	29	6
22. 1,4-Dichlorobenzene-d4	6.154	200	150	11	2	105. Azobenzene*	13.044	200	77	29	6
23. 1,4-Dichlorobenzene	6.18	200	148	11	2	106. 2,4,6-Tribromophenol	13.162	200	330	29	6
24. Benzyl alcohol	6.393	200	79	11	2	107. Sulfotep	13.385	200	322	14	7
25. 1,2-Dichlorobenzene	6.426	200	111	11	2	108. 1,3,5-Trinitrobenzene	13.581	200	213	14	7
26. 2-Methylphenol (o-cresol)	6.611	200	108	11	2	109. Phenacetin	13.581	200	108	14	7
27. 2,2'-Oxybis(1-chloropropane)	6.634	200	45	11	2	110. Phorate	13.618	200	75	14	7
28. N-Nitrosopyrrolidine	6.765	200	100	11	2	111. 4-Bromophenyl phenyl ether	13.677	200	248	14	7
29. Acetophenone	6.825	200	105	11	2	112. Diallyl	13.735	200	43	14	7
30. N-Nitrosomorpholine	6.825	200	56	11	2	113. Hexachlorobenzene	13.788	200	284	14	7
31. N-Nitroso-di-n-propylamine	6.825	200	43	11	2	114. Dimethoate	13.857	200	87	14	7
32. 3-Methylphenol (m-cresol)	6.885	200	80	11	2	115. Atrazine	14.064	200	200	14	7
33. 4-Methylphenol (p-cresol)	6.885	200	107	11	2	116. 4-Aminobiphenyl	14.17	200	169	14	7
34. Indene	6.885	200	117	11	2	117. Pentachlorophenol	14.17	200	266	14	7
35. o-Toluidine	6.885	200	106	11	2	118. Pentachloronitrobenzene (Quintozone)	14.213	200	237	14	7
36. Hexachloroethane	6.989	200	119	11	2	119. Propylamide	14.377	200	173	14	7
37. Nitrobenzene-d5	7.068	200	82	11	2	120. Phenanthrene-d10	14.477	200	188	14	7
38. Nitrobenzene	7.1	200	123	11	2	121. n-Octadecane (C18)	14.52	200	57	14	7
39. N-Nitrosopiperidine	7.347	200	114	17	3	122. Phenanthrene	14.52	200	178	14	7
40. Isophorone	7.531	200	82	17	3	123. Dinoseb	14.616	200	211	14	7
41. 2-Nitrophenol	7.671	200	139	17	3	124. Disulfoton	14.616	200	88	14	7
42. 2,4-Dimethylphenol	7.811	200	122	17	3	125. Anthracene	14.616	200	179	14	7
43. Benzoic acid	7.898	200	105	17	3	126. Carbazole	14.961	200	167	30	8
44. O,O,O-Triethyl phosphorothioate	7.942	200	121	17	3	127. Methyl parathion	15.296	200	109	30	8
45. Bis(2-chloroethoxy)methane	7.969	200	93	17	3	128. Di-n-butyl phthalate	15.779	200	149	30	8
46. 2,4-Dichlorophenol	8.113	200	162	17	3	129. 4-Nitroquinoline-N-oxide	16.026	200	190	30	8
47. α,α-Dimethylphenethylamine (phentermine)	8.178	200	58	17	3	130. Parathion (ethyl parathion)	16.075	200	291	30	8
48. 1,2,4-Trichlorobenzene	8.264	200	180	17	3	131. Methapyrilene hydrochloride	16.273	200	58	30	8
49. Naphthalene-d8	8.34	200	136	17	3	132. Isodrin	16.52	200	193	30	8
50. Naphthalene	8.378	200	128	17	3	133. Fluoranthene	16.806	200	202	30	8
51. α-Terpineol	8.453	200	59	17	3	134. Benzidine	17.124	200	184	30	8
52. 4-Chloroaniline	8.518	200	127	17	3	135. Pyrene	17.223	200	203	30	8
53. 2,6-Dichlorophenol	8.518	200	164	17	3	136. p-Terphenyl-d14	17.627	200	244	27	9
54. Hexachloropropene	8.572	200	213	17	3	137. Aramite isomer 2	17.627	200	175	27	9
55. Hexachlorobutadiene	8.674	200	225	17	3	138. Aramite isomer 1	17.818	200	135	27	9
56. Quinoline	8.993	200	129	15	4	139. p-Dimethylaminoazobenzene	17.905	200	120	27	9
57. ε-Caprolactam	9.085	200	55	15	4	140. Chlorobenzilate	18.041	200	251	27	9
58. 1,4-Phenylenediamine	9.160	200	108	15	4	141. Famphur	18.487	200	218	27	9
59. N-Nitrosodibutylamine	9.195	200	84	15	4	142. 3,3'-Dimethylbenzidine (o-tolidine)	18.553	200	254	15	9
60. 4-Chloro-3-methylphenol	9.477	200	107	15	4	143. Kepone	18.553	200	272	27	9
61. Isosafrole isomer 2	9.563	200	181	15	4	144. Benzyl butyl phthalate	18.64	200	149	27	9
62. 2-Methylnaphthalene	9.681	200	142	15	4	145. Bis(2-ethylhexyl)adipate	18.89	200	129	27	9
63. 1-Methylnaphthalene	9.858	200	141	15	4	146. 2-Acetylaminofluorene	19.015	200	181	27	9
64. 1,2,4,5-Tetrachlorobenzene	10.024	200	216	15	4	147. 4,4'-Methylene-bis(2-chloroaniline)	19.561	200	231	15	9
65. Hexachlorocyclopentadiene	10.024	200	237	15	4	148. 3,3'-Dichlorobenzidine	19.561	200	212	27	9
66. 2,3-Dichloroaniline	10.221	200	161	15	4	149. Benz[a]anthracene	19.561	200	228	15	10
67. Isosafrole isomer 1	10.221	200	104	15	4	150. Chrysene-d12	19.561	200	240	15	10
68. 2,4,6-Trichlorophenol	10.25	200	196	15	4	151. Chrysene	19.61	200	226	15	10
69. 2,4,5-Trichlorophenol	10.303	200	198	15	4	152. Bis(2-ethylhexyl)phthalate	19.824	200	167	15	10
70. 2-Fluorobiphenyl	10.41	200	172	15	4	153. 6-Methylchrysene	20.473	200	242	15	10
71. Safrole	10.549	200	131	15	4	154. Di-n-octyl phthalate	21.031	200	149	15	10
72. Biphenyl	10.586	200	154	15	4	155. Benzo[b]fluoranthene	21.589	200	57	15	10
73. 2-Chloronaphthalene	10.586	200	162	15	4	156. 7,12-Dimethylbenz[a]anthracene	21.589	200	256	15	10
74. 1-Chloronaphthalene	10.624	200	127	15	4	157. Benzo[k]fluoranthene	21.637	200	252	15	10
75. 2-Nitroaniline	10.803	200	138	13	5	158. Benzo[a]pyrene	22.249	200	253	15	10
76. Diphenyl ether	10.803	200	170	15	5	159. Perylene-d12	22.377	200	264	15	10
77. 1,4-Naphthoquinone	10.927	200	158	13	5	160. 3-Methylcholanthrene	23.08	200	268	15	10
78. 1,4-Dinitrobenzene	11.073	200	168	13	5	161. Dibenz[a,h]acridine	24.357	200	279	15	10
79. Dimethylphthalate	11.223	200	163	13	5	162. Dibenz[a,j]acridine	24.48	200	280	15	10
80. 1,3-Dinitrobenzene	11.301	200	76	13	5	163. Indeno[1,2,3-cd]pyrene	24.882	200	277	15	10
81. 2,6-Dinitrotoluene	11.301	200	165	13	5	164. Dibenz[a,h]anthracene	24.968	200	278	15	10
82. 1,2-Dinitrobenzene	11.363	200	50	13	5	165. Benzo[g,h,i]perylene	25.537	200	276	15	10
83. Acenaphthylene	11.363	200	152	13	5						

まとめ：低濃度SVOCのGC-MS分析で感度とピーク形状を両立するために

低濃度の半揮発性有機化合物 (SVOC) をGC-MSで分析する際には、感度低下やピーク形状の劣化が大きな課題となります。本記事で示した結果から、カラムの不活性度を含むシステム全体の最適化により、scanおよびSIMのいずれにおいても高品質な分析が可能であることが示されました。

参考文献

1. E. Pack, J. Hoisington, C. English, R. Dhandapani, and C. Myers, [微量レベル半揮発性有機化合物のGC-MS/MS一斉分析法：データ品質維持と検出下限の改善](#), Application note, EVAN4920-JA, Restek Corporation, 2025
2. [RMX GCカラムカタログ](#), GNBR4923A-JA, Restek Corporation, 2026

注目製品

RMX-5Sil MS GCキャピラリーカラム

カタログ番号	製品名	入数単位
17323	RMX-5Sil MS GCキャピラリーカラム, 30 m, 0.25 mm ID, 0.25 μ m	1



Topaz Single Taper 注入口ライナー

カタログ番号	製品名	入数単位
23303	Topaz, Single Taper Inlet Liner, 4.0 mm x 6.5 x 78.5, for Agilent GCs, w/Quartz Wool, Premium Deactivation	5-pk.



Restekリークディテクタ

カタログ番号	製品名	入数単位
28500	Restek Electronic Leak Detector (持ち運び用ケース、ユニバーサル充電器セット[日本、米国、ヨーロッパ、英国、オーストラリア]、6-ftのUSBケーブルを含む)	1



アプリケーションや製品選定で
お困りですか？
テクニカルサービスまでお問い合わせください！



SVOC MegaMix 150 Kit, 1 mL/ampul; 6 ampuls/kit

1アンプルあたりに含まれる内容

カタログ番号 **31850: 8270 MegaMix**
 1000 µg/mL each in methylene chloride (3-methylphenol and
 4-methylphenol at 500 µg/mL), 1 mL/ampul
 Acenaphthene (83-32-9)
 Acenaphthylene (208-96-8)
 Aniline (62-53-3)
 Anthracene (120-12-7)
 Azobenzene (103-33-3)*
 Benz[a]anthracene (56-55-3)
 Benzo[a]pyrene (50-32-8)
 Benzo[b]fluoranthene (205-99-2)
 Benzo[g,h,i]perylene (191-24-2)
 Benzo[k]fluoranthene (207-08-9)
 Benzyl alcohol (100-51-6)
 Benzyl butyl phthalate (85-68-7)
 Bis(2-chloroethoxy)methane (111-91-1)
 Bis(2-chloroethyl)ether (111-44-4)
 Bis(2-ethylhexyl)adipate (103-23-1)
 Bis(2-ethylhexyl)phthalate (117-81-7)
 4-Bromophenyl phenyl ether (101-55-3)
 Carbazole (86-74-8)
 4-Chloroaniline (106-47-8)
 4-Chloro-3-methylphenol (59-50-7)
 2-Chloronaphthalene (91-58-2)
 2-Chlorophenol (95-57-8)
 4-Chlorophenyl phenyl ether (7005-72-3)
 Chrysene (218-01-9)
 Dibenz[a,h]anthracene (53-70-3)
 Dibenzofuran (132-64-9)
 1,2-Dichlorobenzene (95-50-1)
 1,3-Dichlorobenzene (541-73-1)
 1,4-Dichlorobenzene (106-46-7)
 2,4-Dichlorophenol (120-83-2)
 Diethylphthalate (84-66-2)
 2,4-Dimethylphenol (105-67-9)
 Dimethylphthalate (131-11-3)
 Di-n-butyl phthalate (84-74-2)
 1,2-Dinitrobenzene (528-29-0)
 1,3-Dinitrobenzene (99-65-0)
 1,4-Dinitrobenzene (100-25-4)
 4,6-Dinitro-2-methylphenol (Dinitro-o-cresol) (534-52-1)
 2,4-Dinitrophenol (51-28-5)
 2,4-Dinitrotoluene (121-14-2)
 2,6-Dinitrotoluene (606-20-2)
 Di-n-octyl phthalate (117-84-0)
 Diphenylamine (122-39-4)†
 Fluoranthene (206-44-0)
 Fluorene (86-73-7)
 Hexachlorobenzene (118-74-1)
 Hexachlorobutadiene (87-68-3)
 Hexachlorocyclopentadiene (77-47-4)
 Hexachloroethane (67-72-1)
 Indeno[1,2,3-cd]pyrene (193-39-5)
 Isophorone (78-59-1)
 1-Methylnaphthalene (90-12-0)
 2-Methylnaphthalene (91-57-6)
 2-Methylphenol (o-cresol) (95-48-7)
 3-Methylphenol (m-cresol) (108-39-4)
 4-Methylphenol (p-cresol) (106-44-5)
 Naphthalene (91-20-3)
 2-Nitroaniline (88-74-4)
 3-Nitroaniline (99-09-2)

4-Nitroaniline (100-01-6)
 Nitrobenzene (98-95-3)
 2-Nitrophenol (88-75-5)
 4-Nitrophenol (100-02-7)
 N-Nitrosodimethylamine (62-75-9)
 N-Nitroso-di-n-propylamine (621-64-7)
 2,2'-Oxybis(1-chloropropane) (108-60-1)
 Pentachlorophenol (87-86-5)
 Phenanthrene (85-01-8)
 Phenol (108-95-2)
 Pyrene (129-00-0)
 Pyridine (110-86-1)
 2,3,4,6-Tetrachlorophenol (58-90-2)
 2,3,5,6-Tetrachlorophenol (935-95-5)
 1,2,4-Trichlorobenzene (120-82-1)
 2,4,5-Trichlorophenol (95-95-4)
 2,4,6-Trichlorophenol (88-06-2)

カタログ番号 **31909: SVOC Additions**
 1000 µg/mL each in methylene chloride, 1 mL/ampul
 Acrylamide (79-06-1)
 Benzidine (92-87-5)
 n-Decane (C10) (124-18-5)
 Dibenz(a,h)acridine (226-36-8)
 2,3-Dichloroaniline (608-27-5)
 3,3'-Dichlorobenzidine (91-94-1)
 Dimethoate (60-51-5)
 Dinoseb (88-85-7)
 Disulfoton (298-04-4)
 Famphur (52-85-7)
 n-Hexadecane (C16) (544-76-3)
 Indene (95-13-6)
 Methyl parathion (298-00-0)
 6-Methylchrysene (1705-85-7)
 4,4'-Methylene-bis(2-chloroaniline) (101-14-4)
 n-Octadecane (C18) (593-45-3)
 Parathion (ethyl parathion) (56-38-2)
 Phorate (298-02-2)
 Quinoline (91-22-5)
 Sulfotepp (3689-24-5)
 α-Terpineol (98-55-5)
 o,o,o-Triethyl phosphorothioate (126-68-1)
 Zalophus (thionazine) (297-97-2)

カタログ番号 **32459: Appendix IX mix #1, revised**
 2000 µg/mL each in methylene chloride, 1 mL/ampul
 2-Acetylaminofluorene (53-96-3)
 4-Aminobiphenyl (92-67-1)
 p-Dimethylaminoazobenzene (60-11-7)
 3,3'-Dimethylbenzidine (o-tolidine) (119-93-7)
 α,α-Dimethylphenethylamine (phentermine) (122-09-8)
 1-Naphthylamine (1-aminonaphthalene) (134-32-7)
 2-Naphthylamine (2-aminonaphthalene) (91-59-8)
 N-Nitrosodibutylamine (924-16-3)
 N-Nitrosodiethylamine (55-18-5)
 N-Nitrosomethylethylamine (10595-95-6)
 N-Nitrosomorpholine (59-89-2)
 N-Nitrosopiperidine (100-75-4)
 N-Nitrosopyrrolidine (930-55-2)
 5-Nitro-o-toluidine (99-55-8)
 1,4-Phenylenediamine (106-50-3)
 2-Picoline (109-06-8)
 o-Toluidine (95-53-4)

カタログ番号 **32460: Methapyrilene**
 2000 µg/mL in methylene chloride, 1 mL/ampul
 Methapyrilene hydrochloride (135-23-9)

カタログ番号 **31806: Appendix IX mix #2**
 1000 µg/mL each in methylene chloride, 1 mL/ampul
 Acetophenone (98-86-2)
 Aramite (140-57-8)
 Atrazine (1912-24-9)
 Benzaldehyde (100-52-7)
 Biphenyl (92-52-4)
 ε-Caprolactam (105-60-2)
 Chlorobenzilate (510-15-6)
 1-Chloronaphthalene (90-13-1)
 Diallate (2303-16-4)
 Dibenz[a,j]acridine (224-42-0)
 2,6-Dichlorophenol (87-65-0)
 7,12-Dimethylbenz[a]anthracene (57-97-6)
 1,4-Dioxane (123-91-1)
 Diphenyl ether (101-84-8)
 Ethyl methacrylate (97-63-2)
 Ethyl methanesulfonate (62-50-0)
 Hexachloropropene (1888-71-7)
 Isodrin (465-73-6)
 Isosafrole (cis & trans) (120-58-1)
 Kepone (143-50-0)
 3-Methylcholanthrene (56-49-5)
 Methyl methanesulfonate (66-27-3)
 1,4-Naphthoquinone (130-15-4)
 4-Nitroquinoline-N-oxide (56-57-5)
 Pentachlorobenzene (608-93-5)
 Pentachloroethane (76-01-7)
 Pentachloronitrobenzene (Quintozone) (82-68-8)
 Phenacetin (62-44-2)
 Propylamide (23950-58-5)
 Safrole (94-59-7)
 1,2,4,5-Tetrachlorobenzene (95-94-3)
 1,3,5-Trinitrobenzene (99-35-4)

カタログ番号 **31879: Benzoic acid**
 2000 µg/mL in methylene chloride, 1 mL/ampul
 Benzoic acid (65-85-0)

*1,2-ジフェニルヒドラジン (8270収載分析種) は注入口内でアゾベンジン (混合成分) に分解されます。
 †N-ニトロソジフェニルアミンは反応性が高く、混合物中の他成分の早期分解を誘発する可能性があります。そのため、本混合物の調整にはジフェニルアミンが使用されています。



カタログ番号

31907

入数単位

kit

GC-MS Tuning Mix

1000 µg/mL, Methylene Chloride, 1 mL/ampul

カタログ番号	内容	入数単位
31615	Benzidine (92-87-5) 4,4'-DDT (50-29-3) DFTPP (decafluorotriphenylphosphine) (5074-71-5) Pentachlorophenol (87-86-5)	1

Acid Surrogate Mix (4/89 SOW)

2000 µg/mL, Methanol, 1 mL/ampul

カタログ番号	内容	入数単位
31025	2-Fluorophenol (367-12-4) Phenol-d6 (13127-88-3) 2,4,6-Tribromophenol (118-79-6)	1

Base Neutral Surrogate Mix (4/89 SOW)

1000 µg/mL, Methylene Chloride, 1 mL/ampul

カタログ番号	内容	入数単位
31024	2-Fluorobiphenyl (321-60-8) Nitrobenzene-d5 (4165-60-0) p-Terphenyl-d14 (1718-51-0)	1

Revised SV Internal Standard Mix

4000 µg/mL, Methylene Chloride, 1 mL/ampul

カタログ番号	内容	入数単位
31886	Acenaphthene-d10 (15067-26-2) Chrysene-d12 (1719-03-5) 1,4-Dichlorobenzene-d4 (3855-82-1) 1,4-Dioxane-d8 (17647-74-4) Naphthalene-d8 (1146-65-2) Perylene-d12 (1520-96-3) Phenanthrene-d10 (1517-22-2)	1



お見積りやご購入などのご相談は…

こちらからお問い合わせ下さい！







が優れている理由とは？

高い効果を発揮するTriMax不活性化処理が、不要な表面相互作用から化合物を保護し、幅広い化合物に対してピーク形状と感度を大きく向上させます。





TriMax 不活性化処理あり




不活性化化合物 活性化化合物

TriMax 不活性化処理なし



不活性化化合物 活性化化合物



- 不活性化化合物：
アルカン、アルケン、アルキンなど
- 活性化化合物：
酸、塩基、アルコール、エステル、エーテルなど
- 残存活性点

RMXカラムに関する
最新情報はこちらから！



Restekの特許および商標に関する情報については、www.restek.com/patents-trademarksをご覧ください。Restekからの情報配信や停止設定について変更をご希望の場合は、www.restek.com/subscribeより手続きが可能です。販売に関するお問合せやその他のご質問は、直接弊社までお気軽にご連絡ください。

© 2026 Restek Corporation. All rights reserved.

www.restek.com



Lit. Cat.# EVFA4921-JA