


Pushing the Boundaries of Low-Level GC-MS Semivolatiles Analysis

By Jessi Collier, Erica Pack, Ramkumar Dhandapani, and Colton Myers

- Leverage next-generation TriMax column deactivation to lower detection limits and still meet data quality requirements.
- Exceptionally inert sample flow path ensures sharp, symmetrical peaks for a wide range of challenging semivolatiles.
- Established GC-MS conditions for both a robust scan method and a sensitive SIM method.



Environmental laboratories analyzing semivolatiles are under continual pressure to achieve lower detection limits due to customer demands or the need to reduce solvent use for either regulatory reasons or to reduce costs and improve profitability. MS-based methods, such as EPA Method 8270, are a common approach for analysis, and offer scientists the flexibility of choosing instrumentation to achieve the lowest possible limits of detection via MS. Some labs may opt for sensitive instrumentation, such as GC-MS/MS, or use GC-MS with selected ion monitoring (SIM), while others may opt for more traditional scan methods. SIM and scan methods can typically be performed on the same instrument, where scan methods are simple and versatile, collecting signals for a range of compound masses, while SIM allows for greater sensi-

tivity for specific ions in specific retention time windows. SIM methods are typically more complicated to set up and maintain, but as labs explore strategies to lower detection limits, GC-MS in SIM mode may be considered as an alternative to GC-MS/MS without labs having to budget for a new instrument.

In this study, we established conditions for GC-MS semivolatiles analysis under both scan and SIM modes to provide labs with both a simple, robust scan workflow as well as a highly sensitive SIM alternative to GC-MS/MS (Figure 1). An RMX-5Sil MS column was selected for both methods because it provides standard 5sil selectivity with a next-generation TriMax deactivation that has been proven to provide an exceptionally inert sample flow path [1,2]. As shown in Figure 1, peak asymmetry consistently met quality objectives (≤ 2) for key analytes (benzidine, pentachlorophenol, and 2,4-dinitrophenol). In addition, resolution criteria ($>50\%$) were easily met for challenging PAHs with the valley between benzo[b]fluoranthene and benzo[k]fluoranthene $>94\%$, and the valley between indeno[1,2,3-cd]pyrene, and dibenz[a,h]anthracene $>95\%$. Using an RMX-5Sil MS column under the conditions shown here, excellent results were obtained for a wide range of semivolatiles at very low levels (1 ng on-column in scan mode and 200 pg on-column in SIM mode), allowing labs to work confidently at the lower end of MS detection limits.

Want to customize a method?

Enter your target analytes into Restek's free Pro EZGC chromatogram modeling software and instantly generate optimized method conditions for your specific compound list.

Try it now!




Pure Chromatography

www.restek.com

Figure 1: RMX-5Sil MS columns provide outstanding chromatographic results for 165 semivolatiles analyzed at low levels by GC-MS.

Scan Mode (1 ng on-column)



Column RMX-5Sil MS, 30 m, 0.25 mm ID, 0.25 µm (cat.# 17323)
Standard/Sample

Methapyrilene (cat.# 32460)
 Appendix IX mix #1, revised (cat.# 32459)
 SVOC additions (cat.# 31909)
 Benzoic acid (cat.# 31879)
 8270 MegaMix (cat.# 31850)
 Appendix IX mix #2 (cat.# 31806)
 Acid surrogate mix (4/89 SOW) (cat.# 31025)
 Base neutral surrogate mix (4/89 SOW) (cat.# 31024)
 Revised SV internal standard mix (cat.# 31886)

Diluent:
Conc.:

Methylene chloride
 1 µg/mL

Injection

Inj. Vol.:
Liner:
Inj. Temp.:
Split Vent Flow Rate:

1 µL split (split ratio 10:1)
 Topaz 4.0 mm ID single taper liner w/wool (cat.# 23303)
 250 °C
 12 mL/min

Oven

Oven Temp.:

40 °C (hold 1 min) to 280 °C at 12.4 °C/min to 330 °C at 3.3 °C/min (hold 3.65 min)

Carrier Gas

Flow Rate:
Linear Velocity:

He, constant flow
 1.2 mL/min
 39.7 cm/sec

Detector

Mode:
Scan Program:

MS
 Scan

Group	Start Time (min)	Scan Range (amu)	Scan Rate (scans/sec)
1	2	35-550	5

Transfer Line Temp.:

280 °C

Analyzer Type:

Quadrupole

Source Temp.:

330 °C

Quad Temp.:

180 °C

Solvent Delay Time:

2 min

Tune Type:

PFTBA

Ionization Mode:

EI

Instrument

Agilent 7890B GC & 5977B MSD

Notes

To simplify ordering, the SVOC MegaMix 150 kit (cat.# 31907) contains one ampul each of the following standards that were used in this experiment.



- 8270 MegaMix (cat.# 31850)
- SVOC Additions standard (cat.# 31909)
- Appendix IX mix #1, revised (cat.# 32459)
- Methapyrilene (cat.# 32460)
- Appendix IX mix #2 (cat.# 31806)
- Benzoic acid (cat.# 31879)

Figure 1: RMX-5Sil MS columns provide outstanding chromatographic results for 165 semivolatiles analyzed at low levels by GC-MS. (cont.)

Peaks	t _r (min)	Conc. (µg/mL)	Peaks	t _r (min)	Conc. (µg/mL)
1. 1,4-Dioxane-d8	2.400	1	84. 3-Nitroaniline	11.595	1
2. 1,4-Dioxane	2.428	1	85. Acenaphthene-d10	11.647	1
3. N-Nitrosodimethylamine	2.702	1	86. Acenaphthene	11.704	1
4. Pyridine	2.774	1	87. 2,4-Dinitrophenol	11.809	1
5. Ethyl methacrylate	3.175	1	88. 4-Nitrophenol	11.996	1
6. 2-Picoline	3.554	1	89. Pentachlorobenzene	11.996	1
7. N-Nitrosomethylethylamine	3.648	1	90. Dibenzofuran	12.040	1
8. Methyl methanesulfonate	4.021	1	91. 2,4-Dinitrotoluene	12.070	1
9. 2-Fluorophenol	4.258	1	92. 1-Naphthylamine (1-aminonaphthalene)	12.184	1
10. Acrylamide	4.543	1	93. 2,3,5,6-Tetrachlorophenol	12.236	1
11. N-Nitrosodiethylamine	4.875	1	94. 2,3,4,6-Tetrachlorophenol	12.319	1
12. Ethyl methanesulfonate	4.977	1	95. 2-Naphthylamine (2-aminonaphthalene)	12.341	1
13. Benzaldehyde	5.454	1	96. Diethylphthalate	12.602	1
14. Phenol-d6	5.612	1	97. n-Hexadecane (C16)	12.699	1
15. Phenol	5.629	1	98. Fluorene	12.699	1
16. Aniline	5.638	1	99. Zalophus (thionazine)	12.745	1
17. Pentachloroethane	5.705	1	100. 4-Chlorophenyl phenyl ether	12.745	1
18. Bis(2-chloroethyl)ether	5.743	1	101. 4-Nitroaniline	12.745	1
19. 2-Chlorophenol	5.814	1	102. 5-Nitro-o-toluidine	12.745	1
20. n-Decane (C10)	5.962	1	103. 4,6-Dinitro-2-methylphenol (Dinitro-o-cresol)	12.853	1
21. 1,3-Dichlorobenzene	6.066	1	104. Diphenylamine†	12.979	1
22. 1,4-Dichlorobenzene-d4	6.153	1	105. Azobenzene*	13.042	1
23. 1,4-Dichlorobenzene	6.180	1	106. 2,4,6-Tribromophenol	13.162	1
24. Benzyl alcohol	6.393	1	107. Sulfotep	13.379	1
25. 1,2-Dichlorobenzene	6.426	1	108. 1,3,5-Trinitrobenzene	13.579	1
26. 2-Methylphenol (o-cresol)	6.612	1	109. Phorate	13.579	1
27. 2,2'-Oxybis(1-chloropropane)	6.634	1	110. Phenacetin	13.625	1
28. N-Nitrosopyrrolidine	6.765	1	111. 4-Bromophenyl phenyl ether	13.677	1
29. Acetophenone	6.825	1	112. Diallylate	13.734	1
30. N-Nitrosomorpholine	6.836	1	113. Hexachlorobenzene	13.785	1
31. N-Nitroso-di-n-propylamine	6.847	1	114. Dimethoate	13.865	1
32. 3-Methylphenol (m-cresol)	6.885	1	115. Atrazine	14.071	1
33. 4-Methylphenol (p-cresol)	6.885	1	116. 4-Aminobiphenyl	14.168	1
34. Indene	6.885	1	117. Pentachlorophenol	14.170	1
35. o-Toluidine	6.890	1	118. Pentachloronitrobenzene (Quintozone)	14.214	1
36. Hexachloroethane	6.989	1	119. Propyzamide	14.380	1
37. Nitrobenzene-d5	7.067	1	120. Phenanthrene-d10	14.477	1
38. Nitrobenzene	7.101	1	121. n-Octadecane (C18)	14.517	1
39. N-Nitrosopiperidine	7.347	1	122. Phenanthrene	14.517	1
40. Isophorone	7.536	1	123. Dinoseb	14.599	1
41. 2-Nitrophenol	7.673	1	124. Disulfoton	14.599	1
42. 2,4-Dimethylphenol	7.816	1	125. Anthracene	14.614	1
43. Benzoic acid	7.943	1	126. Carbazole	14.963	1
44. O,O,O-Triethyl phosphorothioate	7.943	1	127. Methyl parathion	15.295	1
45. Bis(2-chloroethoxy)methane	7.970	1	128. Di-n-butyl phthalate	15.775	1
46. 2,4-Dichlorophenol	8.113	1	129. 4-Nitroquinoline-N-oxide	16.027	1
47. α,α-Dimethylphenethylamine (phentermine)	8.268	1	130. Parathion (ethyl parathion)	16.078	1
48. 1,2,4-Trichlorobenzene	8.268	1	131. Methapyrilene hydrochloride	16.273	1
49. Naphthalene-d8	8.342	1	132. Isodrin	16.518	1
50. Naphthalene	8.376	1	133. Fluoranthene	16.804	1
51. α-Terpineol	8.451	1	134. Benzidine	17.127	1
52. 4-Chloroaniline	8.519	1	135. Pyrene	17.222	1
53. 2,6-Dichlorophenol	8.519	1	136. p-Terphenyl-d14	17.622	1
54. Hexachloropropene	8.571	1	137. Aramite isomer 2	17.668	1
55. Hexachlorobutadiene	8.674	1	138. Aramite isomer 1	17.822	1
56. Quinoline	8.988	1	139. p-Dimethylaminoazobenzene	17.908	1
57. ε-Caprolactam	9.085	1	140. Chlorobenzilate	18.034	1
58. 1,4-Phenylenediamine	9.160	1	141. Famphur	18.485	1
59. N-Nitrosodibutylamine	9.195	1	142. 3,3'-Dimethylbenzidine (o-tolidine)	18.548	1
60. 4-Chloro-3-methylphenol	9.480	1	143. Kepone	18.548	1
61. Isosafrole isomer 2	9.560	1	144. Benzyl butyl phthalate	18.640	1
62. 2-Methylnaphthalene	9.680	1	145. Bis(2-ethylhexyl)adipate	18.891	1
63. 1-Methylnaphthalene	9.857	1	146. 2-Acetylaminofluorene	19.017	1
64. 1,2,4,5-Tetrachlorobenzene	10.023	1	147. 4,4'-Methylene-bis(2-chloroaniline)	19.108	1
65. Hexachlorocyclopentadiene	10.023	1	148. 3,3'-Dichlorobenzidine	19.560	1
66. 2,3-Dichloroaniline	10.217	1	149. Benz[a]anthracene	19.560	1
67. Isosafrole isomer 1	10.217	1	150. Chrysene-d12	19.560	1
68. 2,4,6-Trichlorophenol	10.248	1	151. Chrysene	19.612	1
69. 2,4,5-Trichlorophenol	10.303	1	152. Bis(2-ethylhexyl)phthalate	19.823	1
70. 2-Fluorobiphenyl	10.412	1	153. 6-Methylchrysene	20.475	1
71. Safrole	10.547	1	154. Di-n-octyl phthalate	21.030	1
72. Biphenyl	10.585	1	155. Benzo[b]fluoranthene	21.590	1
73. 2-Chloronaphthalene	10.585	1	156. 7,12-Dimethylbenz[a]anthracene	21.590	1
74. 1-Chloronaphthalene	10.623	1	157. Benzo[k]fluoranthene	21.641	1
75. 2-Nitroaniline	10.801	1	158. Benzo[a]pyrene	22.248	1
76. Diphenyl ether	10.801	1	159. Perylene-d12	22.379	1
77. 1,4-Naphthoquinone	10.926	1	160. 3-Methylcholanthrene	23.082	1
78. 1,4-Dinitrobenzene	11.075	1	161. Dibenz[a,h]acridine	24.357	1
79. Dimethylphthalate	11.224	1	162. Dibenz[a,j]acridine	24.483	1
80. 1,3-Dinitrobenzene	11.304	1	163. Indeno[1,2,3-cd]pyrene	24.883	1
81. 2,6-Dinitrotoluene	11.304	1	164. Dibenz[a,h]anthracene	24.969	1
82. 1,2-Dinitrobenzene	11.361	1	165. Benzo[g,h,i]perylene	25.541	1
83. Acenaphthylene	11.361	1			

Figure 1: RMX-5Sil MS columns provide outstanding chromatographic results for 165 semivolatiles analyzed at low levels by GC-MS. (cont.)

SIM Mode (200 pg on-column)



Column RMX-5Sil MS, 30 m, 0.25 mm ID, 0.25 μ m (cat.# 17323)
Standard/Sample

Methapyrilene (cat.# 32460)
Appendix IX mix #1, revised (cat.# 32459)
SVOC additions (cat.# 31909)
Benzoic acid (cat.# 31879)
8270 MegaMix (cat.# 31850)
Appendix IX mix #2 (cat.# 31806)
Acid surrogate mix (4/89 SOW) (cat.# 31025)
Base neutral surrogate mix (4/89 SOW) (cat.# 31024)
Revised SV internal standard mix (cat.# 31886)
Methylene chloride
200 ng/mL

Diluent:
Conc.:

Injection

Inj. Vol.: 1 μ L split (split ratio 10:1)
Liner: Topaz 4.0 mm ID single taper liner w/wool (cat.# 23303)
Inj. Temp.: 250 °C
Split Vent Flow Rate: 12 mL/min

Oven

Oven Temp.: 40 °C (hold 1 min) to 280 °C at 12.4 °C/min to 330 °C at 3.3 °C/min (hold 3.65 min)

Carrier Gas

Flow Rate: He, constant flow
1.2 mL/min

Linear Velocity: 39.7 cm/sec

Detector

MS
Mode: SIM

SIM Program:

Group	Start Time (min)	Ion(s) (m/z)	Dwell (ms)
1	2	see peak list	300
2	5.2	see peak list	286
3	7.2	see peak list	289
4	9.35	see peak list	285
5	10.7	see peak list	273
6	12.45	see peak list	290
7	13.29	see peak list	280
8	14.8	see peak list	300
9	17.4	see peak list	297
10	19.3	see peak list	285

Transfer Line Temp.: 280 °C

Analyzer Type: Quadrupole

Source Temp.: 330 °C

Quad Temp.: 180 °C

Solvent Delay Time: 2 min

Tune Type: PFTBA

Ionization Mode: EI

Instrument Notes

Agilent 7890B GC & 5977B MSD

To simplify ordering, the SVOC MegaMix 150 kit (cat.# 31907) contains one ampul each of the following standards that were used in this experiment.

- 8270 MegaMix (cat.# 31850)
- SVOC Additions standard (cat.# 31909)
- Appendix IX mix #1, revised (cat.# 32459)
- Methapyrilene (cat.# 32460)
- Appendix IX mix #2 (cat.# 31806)
- Benzoic acid (cat.# 31879)

Figure 1: RMX-5Sil MS columns provide outstanding chromatographic results for 165 semivolatiles analyzed at low levels by GC-MS. (cont.)

Peaks	t _r (min)	Conc. (ng/mL)	SIM Ion	Dwell Time (ms)	Ion Group	Peaks	t _r (min)	Conc. (ng/mL)	SIM Ion	Dwell Time (ms)	Ion Group
1. 1,4-Dioxane-d8	2.401	200	96	25	1	84. 3-Nitroaniline	11.591	200	65	13	5
2. 1,4-Dioxane	2.427	200	88	25	1	85. Acenaphthene-d10	11.649	200	162	13	5
3. N-Nitrosodimethylamine	2.702	200	74	25	1	86. Acenaphthene	11.706	200	153	13	5
4. Pyridine	2.774	200	79	25	1	87. 2,4-Dinitrophenol	11.809	200	184	13	5
5. Ethyl methacrylate	3.175	200	69	25	1	88. 4-Nitrophenol	11.996	200	139	13	5
6. 2-Picoline	3.554	200	93	25	1	89. Pentachlorobenzene	11.996	200	250	13	5
7. N-Nitrosomethylethylamine	3.648	200	42	25	1	90. Dibenzofuran	12.043	200	84	13	5
8. Methyl methanesulfonate	4.021	200	80	25	1	91. 2,4-Dinitrotoluene	12.183	200	89	13	5
9. 2-Fluorophenol	4.258	200	112	25	1	92. 1-Naphthylamine (1-aminonaphthalene)	12.235	200	143	13	5
10. Acrylamide	4.543	200	71	25	1	93. 2,3,5,6-Tetrachlorophenol	12.339	200	232	13	5
11. N-Nitrosodiethylamine	4.977	200	102	25	1	94. 2,3,4,6-Tetrachlorophenol	12.339	200	230	13	5
12. Ethyl methanesulfonate	4.977	200	109	25	1	95. 2-Naphthylamine (2-aminonaphthalene)	12.339	200	115	13	5
13. Benzaldehyde	5.454	200	77	11	2	96. Diethylphthalate	12.596	200	149	13	5
14. Phenol-d6	5.612	200	99	11	2	97. n-Hexadecane (C16)	12.703	200	57	29	6
15. Phenol	5.633	200	94	11	2	98. Fluorene	12.703	200	166	29	6
16. Aniline	5.633	200	93	11	2	99. Zolophus (thionazine)	12.746	200	107	29	6
17. Pentachloroethane	5.7	200	167	11	2	100. 4-Chlorophenyl phenyl ether	12.746	200	204	29	6
18. Bis(2-chloroethyl)ether	5.743	200	63	11	2	101. 4-Nitroaniline	12.746	200	65	29	6
19. 2-Chlorophenol	5.814	200	128	11	2	102. 5-Nitro-o-toluidine	12.746	200	152	29	6
20. n-Decane (C10)	5.962	200	57	11	2	103. 4,6-Dinitro-2-methylphenol (Dinitro-o-cresol)	12.852	200	198	29	6
21. 1,3-Dichlorobenzene	6.066	200	146	11	2	104. Diphenylamine†	12.98	200	169	29	6
22. 1,4-Dichlorobenzene-d4	6.154	200	150	11	2	105. Azobenzene*	13.044	200	77	29	6
23. 1,4-Dichlorobenzene	6.18	200	148	11	2	106. 2,4,6-Tribromophenol	13.162	200	330	29	6
24. Benzyl alcohol	6.393	200	79	11	2	107. Sulfotep	13.385	200	322	14	7
25. 1,2-Dichlorobenzene	6.426	200	111	11	2	108. 1,3,5-Trinitrobenzene	13.581	200	213	14	7
26. 2-Methylphenol (o-cresol)	6.611	200	108	11	2	109. Phenacetin	13.581	200	108	14	7
27. 2,2'-Oxybis(1-chloropropane)	6.634	200	45	11	2	110. Phorate	13.618	200	75	14	7
28. N-Nitrosopyrrolidine	6.765	200	100	11	2	111. 4-Bromophenyl phenyl ether	13.677	200	248	14	7
29. Acetophenone	6.825	200	105	11	2	112. Diallate	13.735	200	43	14	7
30. N-Nitrosomorpholine	6.825	200	56	11	2	113. Hexachlorobenzene	13.788	200	284	14	7
31. N-Nitroso-di-n-propylamine	6.825	200	43	11	2	114. Dimethoate	13.857	200	87	14	7
32. 3-Methylphenol (m-cresol)	6.885	200	80	11	2	115. Atrazine	14.064	200	200	14	7
33. 4-Methylphenol (p-cresol)	6.885	200	107	11	2	116. 4-Aminobiphenyl	14.17	200	169	14	7
34. Indene	6.885	200	117	11	2	117. Pentachlorophenol	14.17	200	266	14	7
35. o-Toluidine	6.885	200	106	11	2	118. Pentachloronitrobenzene (Quintozene)	14.213	200	237	14	7
36. Hexachloroethane	6.989	200	119	11	2	119. Propylamide	14.377	200	173	14	7
37. Nitrobenzene-d5	7.068	200	82	11	2	120. Phenanthrene-d10	14.477	200	188	14	7
38. Nitrobenzene	7.1	200	123	11	2	121. n-Octadecane (C18)	14.52	200	57	14	7
39. N-Nitrosopiperidine	7.347	200	114	17	3	122. Phenanthrene	14.52	200	178	14	7
40. Isophorone	7.531	200	82	17	3	123. Dinoseb	14.616	200	211	14	7
41. 2-Nitrophenol	7.671	200	139	17	3	124. Disulfoton	14.616	200	88	14	7
42. 2,4-Dimethylphenol	7.811	200	122	17	3	125. Anthracene	14.616	200	179	14	7
43. Benzoic acid	7.898	200	105	17	3	126. Carbazole	14.961	200	167	30	8
44. O,O,O-Triethyl phosphorothioate	7.942	200	121	17	3	127. Methyl parathion	15.296	200	109	30	8
45. Bis(2-chloroethoxy)methane	7.969	200	93	17	3	128. Di-n-butyl phthalate	15.779	200	149	30	8
46. 2,4-Dichlorophenol	8.113	200	162	17	3	129. 4-Nitroquinoline-N-oxide	16.026	200	190	30	8
47. α,α-Dimethylphenethylamine (phentermine)	8.178	200	58	17	3	130. Parathion (ethyl parathion)	16.075	200	291	30	8
48. 1,2,4-Trichlorobenzene	8.264	200	180	17	3	131. Methapyriline hydrochloride	16.273	200	58	30	8
49. Naphthalene-d8	8.34	200	136	17	3	132. Isodrin	16.52	200	193	30	8
50. Naphthalene	8.378	200	128	17	3	133. Fluoranthene	16.806	200	202	30	8
51. α-Terpineol	8.453	200	59	17	3	134. Benzidine	17.124	200	184	30	8
52. 4-Chloroaniline	8.518	200	127	17	3	135. Pyrene	17.223	200	203	30	8
53. 2,6-Dichlorophenol	8.518	200	164	17	3	136. p-Terphenyl-d14	17.627	200	244	27	9
54. Hexachloropropene	8.572	200	213	17	3	137. Aramite isomer 2	17.627	200	175	27	9
55. Hexachlorobutadiene	8.674	200	225	17	3	138. Aramite isomer 1	17.818	200	135	27	9
56. Quinoline	8.993	200	129	15	4	139. p-Dimethylaminoazobenzene	17.905	200	120	27	9
57. ε-Caprolactam	9.085	200	55	15	4	140. Chlorobenzilate	18.041	200	251	27	9
58. 1,4-Phenylenediamine	9.160	200	108	15	4	141. Famphur	18.487	200	218	27	9
59. N-Nitrosodibutylamine	9.195	200	84	15	4	142. 3,3'-Dimethylbenzidine (o-tolidine)	18.553	200	254	15	9
60. 4-Chloro-3-methylphenol	9.477	200	107	15	4	143. Kepone	18.553	200	272	27	9
61. Isosafrole isomer 2	9.563	200	181	15	4	144. Benzyl butyl phthalate	18.64	200	149	27	9
62. 2-Methylnaphthalene	9.681	200	142	15	4	145. Bis(2-ethylhexyl)adipate	18.89	200	129	27	9
63. 1-Methylnaphthalene	9.858	200	141	15	4	146. 2-Acetylaminofluorene	19.015	200	181	27	9
64. 1,2,4,5-Tetrachlorobenzene	10.024	200	216	15	4	147. 4,4'-Methylene-bis(2-chloroaniline)	19.561	200	231	15	9
65. Hexachlorocyclopentadiene	10.024	200	237	15	4	148. 3,3'-Dichlorobenzidine	19.561	200	212	27	9
66. 2,3-Dichloroaniline	10.221	200	161	15	4	149. Benz[a]anthracene	19.561	200	228	15	10
67. Isosafrole isomer 1	10.221	200	104	15	4	150. Chrysene-d12	19.561	200	240	15	10
68. 2,4,6-Trichlorophenol	10.25	200	196	15	4	151. Chrysene	19.61	200	226	15	10
69. 2,4,5-Trichlorophenol	10.303	200	198	15	4	152. Bis(2-ethylhexyl)phthalate	19.824	200	167	15	10
70. 2-Fluorobiphenyl	10.41	200	172	15	4	153. 6-Methylchrysene	20.473	200	242	15	10
71. Safrole	10.549	200	131	15	4	154. Di-n-octyl phthalate	21.031	200	149	15	10
72. Biphenyl	10.586	200	154	15	4	155. Benzo[b]fluoranthene	21.589	200	57	15	10
73. 2-Chloronaphthalene	10.586	200	162	15	4	156. 7,12-Dimethylbenz[a]anthracene	21.589	200	256	15	10
74. 1-Chloronaphthalene	10.624	200	127	15	4	157. Benzo[k]fluoranthene	21.637	200	252	15	10
75. 2-Nitroaniline	10.803	200	138	13	5	158. Benzo[a]pyrene	22.249	200	253	15	10
76. Diphenyl ether	10.803	200	170	15	5	159. Perylene-d12	22.377	200	264	15	10
77. 1,4-Naphthoquinone	10.927	200	158	13	5	160. 3-Methylcholanthrene	23.08	200	268	15	10
78. 1,4-Dinitrobenzene	11.073	200	168	13	5	161. Dibenz[a,h]acridine	24.357	200	279	15	10
79. Dimethylphthalate	11.223	200	163	13	5	162. Dibenz[a,j]acridine	24.48	200	280	15	10
80. 1,3-Dinitrobenzene	11.301	200	76	13	5	163. Indeno[1,2,3-cd]pyrene	24.882	200	277	15	10
81. 2,6-Dinitrotoluene	11.301	200	165	13	5	164. Dibenz[a,h]anthracene	24.968	200	278	15	10
82. 1,2-Dinitrobenzene	11.363	200	50	13	5	165. Benzo[g,h,i]perylene	25.537	200	276	15	10
83. Acenaphthylene	11.363	200	152	13	5						

References

- 1. E. Pack, J. Hoisington, C. English, R. Dhandapani, and C. Myers, Comprehensive trace-level semivolatiles analysis by GC-MS/MS (EPA Method 8270E), Application note, EVAN4919-US, Restek Corporation, 2025.
- 2. RMX GC columns brochure, GNBR4923-UNV, Restek Corporation, 2026.

Featured Products

RMX-5Sil MS GC Capillary Column

Catalog No.	Product Name	Units
17323	RMX-5Sil MS GC Capillary Column, 30 m, 0.25 mm ID, 0.25 µm	ea.



Topaz Single Taper Inlet Liner

Catalog No.	Product Name	Units
23303	Topaz, Single Taper Inlet Liner, 4.0 mm x 6.5 x 78.5, for Agilent GCs, w/Quartz Wool, Premium Deactivation	5-pk.



Restek Electronic Leak Detector

Catalog No.	Product Name	Units
28500	Restek Electronic Leak Detector (includes carrying case; universal AC power adaptor [U.S., UK, Europe, Australia, Japan]; 6-ft USB charging cable)	ea.



Need support with applications or product selection?
Contact Technical Service!



Want help with quoting, purchasing, and more?
Contact Us!



SVOC MegaMix 150 Kit, 1 mL/ampul; 6 ampuls/kit

Contains one ampul of each of the following.

Cat.# 31850: 8270 MegaMix

1000 µg/mL each in methylene chloride (3-methylphenol and 4-methylphenol at 500 µg/mL), 1 mL/ampul

Acenaphthene (83-32-9)
Acenaphthylene (208-96-8)
Aniline (62-53-3)
Anthracene (120-12-7)
Azobenzene (103-33-3)*
Benz[a]anthracene (56-55-3)
Benzo[a]pyrene (50-32-8)
Benzo[b]fluoranthene (205-99-2)
Benzo[g,h,i]perylene (191-24-2)
Benzo[k]fluoranthene (207-08-9)
Benzyl alcohol (100-51-6)
Benzyl butyl phthalate (85-68-7)
Bis(2-chloroethoxy)methane (111-91-1)
Bis(2-chloroethyl)ether (111-44-4)
Bis(2-ethylhexyl)adipate (103-23-1)
Bis(2-ethylhexyl)phthalate (117-81-7)
4-Bromophenyl phenyl ether (101-55-3)
Carbazole (86-74-8)
4-Chloroaniline (106-47-8)
4-Chloro-3-methylphenol (59-50-7)
2-Chloronaphthalene (91-58-7)
2-Chlorophenol (95-57-8)
4-Chlorophenyl phenyl ether (7005-72-3)
Chrysene (218-01-9)
Dibenz[a,h]anthracene (53-70-3)
Dibenzofuran (132-64-9)
1,2-Dichlorobenzene (95-50-1)
1,3-Dichlorobenzene (541-73-1)
1,4-Dichlorobenzene (106-46-7)
2,4-Dichlorophenol (120-83-2)
Diethylphthalate (84-66-2)
2,4-Dimethylphenol (105-67-9)
Dimethylphthalate (131-11-3)
Di-n-butyl phthalate (84-74-2)
1,2-Dinitrobenzene (528-29-0)
1,3-Dinitrobenzene (99-65-0)
1,4-Dinitrobenzene (100-25-4)
4,6-Dinitro-2-methylphenol (Dinitro-o-cresol) (534-52-1)
2,4-Dinitrophenol (51-28-5)
2,4-Dinitrotoluene (121-14-2)
2,6-Dinitrotoluene (606-20-2)
Di-n-octyl phthalate (117-84-0)
Diphenylamine (122-39-4)†
Fluoranthene (206-44-0)
Fluorene (86-73-7)
Hexachlorobenzene (118-74-1)
Hexachlorobutadiene (87-68-3)
Hexachlorocyclopentadiene (77-47-4)
Hexachloroethane (67-72-1)
Indeno[1,2,3-cd]pyrene (193-39-5)
Isophorone (78-59-1)
1-Methylnaphthalene (90-12-0)
2-Methylnaphthalene (91-57-6)
2-Methylphenol (o-cresol) (95-48-7)
3-Methylphenol (m-cresol) (108-39-4)
4-Methylphenol (p-cresol) (106-44-5)
Naphthalene (91-20-3)
2-Nitroaniline (88-74-4)
3-Nitroaniline (99-09-2)

4-Nitroaniline (100-01-6)
Nitrobenzene (98-95-3)
2-Nitrophenol (88-75-5)
4-Nitrophenol (100-02-7)
N-Nitrosodimethylamine (62-75-9)
N-Nitroso-di-n-propylamine (621-64-7)
2,2'-Oxybis(1-chloropropane) (108-60-1)
Pentachlorophenol (87-86-5)
Phenanthrene (85-01-8)
Phenol (108-95-2)
Pyrene (129-00-0)
Pyridine (110-86-1)
2,3,4,6-Tetrachlorophenol (58-90-2)
2,3,5,6-Tetrachlorophenol (935-95-5)
1,2,4-Trichlorobenzene (120-82-1)
2,4,5-Trichlorophenol (95-95-4)
2,4,6-Trichlorophenol (88-06-2)

Cat.# 31909: SVOC Additions

1000 µg/mL each in methylene chloride, 1 mL/ampul

Acrylamide (79-06-1)
Benzidine (92-87-5)
n-Decane (C10) (124-18-5)
Dibenz(a,h)acridine (226-36-8)
2,3-Dichloroaniline (608-27-5)
3,3'-Dichlorobenzidine (91-94-1)
Dimethoate (60-51-5)
Dinoseb (88-85-7)
Disulfoton (298-04-4)
Famphur (52-85-7)
n-Hexadecane (C16) (544-76-3)
Indene (95-13-6)
Methyl parathion (298-00-0)
6-Methylchrysene (1705-85-7)
4,4'-Methylene-bis(2-chloroaniline) (101-14-4)
n-Octadecane (C18) (593-45-3)
Parathion (ethyl parathion) (56-38-2)
Phorate (298-02-2)
Quinoline (91-22-5)
Sulfotepp (3689-24-5)
α-Terpineol (98-55-5)
o,o,o-Triethyl phosphorothioate (126-68-1)
Zalophus (thionazine) (297-97-2)

Cat.# 32459: Appendix IX mix #1, revised

2000 µg/mL each in methylene chloride, 1 mL/ampul

2-Acetylaminofluorene (53-96-3)
4-Aminobiphenyl (92-67-1)
p-Dimethylaminoazobenzene (60-11-7)
3,3'-Dimethylbenzidine (o-tolidine) (119-93-7)
α,α-Dimethylphenethylamine (phentermine) (122-09-8)
1-Naphthylamine (1-aminonaphthalene) (134-32-7)
2-Naphthylamine (2-aminonaphthalene) (91-59-8)
N-Nitrosodibutylamine (924-16-3)
N-Nitrosodiethylamine (55-18-5)
N-Nitrosomethylethylamine (10595-95-6)
N-Nitrosomorpholine (59-89-2)
N-Nitrosopiperidine (100-75-4)
N-Nitrosopyrrolidine (930-55-2)
5-Nitro-o-toluidine (99-55-8)
1,4-Phenylenediamine (106-50-3)
2-Picoline (109-06-8)
o-Toluidine (95-53-4)

Cat.# 32460: Methapyrilene

2000 µg/mL in methylene chloride, 1 mL/ampul
Methapyrilene hydrochloride (135-23-9)

Cat.# 31806: Appendix IX mix #2

1000 µg/mL each in methylene chloride, 1 mL/ampul

Acetophenone (98-86-2)
Aramite (140-57-8)
Atrazine (1912-24-9)
Benzaldehyde (100-52-7)
Biphenyl (92-52-4)
ε-Caprolactam (105-60-2)
Chlorobenzilate (510-15-6)
1-Chloronaphthalene (90-13-1)
Diallate (2303-16-4)
Dibenz[a,j]acridine (224-42-0)
2,6-Dichlorophenol (87-65-0)
7,12-Dimethylbenz[a]anthracene (57-97-6)
1,4-Dioxane (123-91-1)
Diphenyl ether (101-84-8)
Ethyl methacrylate (97-63-2)
Ethyl methanesulfonate (62-50-0)
Hexachloropropene (1888-71-7)
Isodrin (465-73-6)
Isosafrole (cis & trans) (120-58-1)
Kepone (143-50-0)
3-Methylcholanthrene (56-49-5)
Methyl methanesulfonate (66-27-3)
1,4-Naphthoquinone (130-15-4)
4-Nitroquinoline-N-oxide (56-57-5)
Pentachlorobenzene (608-93-5)
Pentachloroethane (76-01-7)
Pentachloronitrobenzene (Quintozone) (82-68-8)
Phenacetin (62-44-2)
Propylamide (23950-58-5)
Safrole (94-59-7)
1,2,4,5-Tetrachlorobenzene (95-94-3)
1,3,5-Trinitrobenzene (99-35-4)

Cat.# 31879: Benzoic acid

2000 µg/mL in methylene chloride, 1 mL/ampul
Benzoic acid (65-85-0)

*1,2-diphenylhydrazine (8270-listed analyte) decomposes to azobenzene (mix component) in the injector.

†N-Nitrosodiphenylamine is a reactive species that can initiate premature decomposition of other compounds in the mix. For these reasons, diphenylamine is used in the preparation of this mixture.



Catalog No.

31907

Units

kit

GC-MS Tuning Mix

1000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31615	Benzidine (92-87-5) 4,4'-DDT (50-29-3) DFTPP (decafluorotriphenylphosphine) (5074-71-5) Pentachlorophenol (87-86-5)	ea.

Acid Surrogate Mix (4/89 SOW)

2000 µg/mL, Methanol, 1 mL/ampul

Catalog No.	Contains	Units
31025	2-Fluorophenol (367-12-4) Phenol-d6 (13127-88-3) 2,4,6-Tribromophenol (118-79-6)	ea.

Base Neutral Surrogate Mix (4/89 SOW)

1000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31024	2-Fluorobiphenyl (321-60-8) Nitrobenzene-d5 (4165-60-0) p-Terphenyl-d14 (1718-51-0)	ea.

Revised SV Internal Standard Mix

4000 µg/mL, Methylene Chloride, 1 mL/ampul

Catalog No.	Contains	Units
31886	Acenaphthene-d10 (15067-26-2) Chrysene-d12 (1719-03-5) 1,4-Dichlorobenzene-d4 (3855-82-1) 1,4-Dioxane-d8 (17647-74-4) Naphthalene-d8 (1146-65-2) Perylene-d12 (1520-96-3) Phenanthrene-d10 (1517-22-2)	ea.





What Makes RMX Columns Better?

Highly Effective TriMax Deactivation Protects Analytes From Surface Interactions, Improving Peak Shape and Sensitivity for a Wide Range of Compound Chemistries



TriMax Deactivation




Inactive



Active


Non-TriMax Deactivation



Inactive



Active



- Inactive compounds: *alkanes, alkenes, alkynes, etc.*
- Active compounds: *acids, bases, alcohols, esters, ethers, etc.*
- Residual active site

Learn More About
RMX Columns Today!



For information on Restek patents and trademarks, visit www.restek.com/patents-trademarks. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2026 Restek Corporation. All rights reserved.

www.restek.com



Lit. Cat.# EVFA4921-UNV