

Reduce Helium Consumption by 68% Using Nitrogen Purge Gas for VOCs in Water

- Save 490 mL of helium per sample by switching to nitrogen purge gas.
- Spend less money on lab gases and reduce your dependence on helium.
- Easily resolve critical Method 524.4 compounds using an Rtx-VMS column.



Reduce Helium Consumption by 68% Using Nitrogen Purge Gas for VOCs in Water

By now, everyone has felt the impact of the helium supply problem—even regulatory agencies. In order to help laboratories reduce helium use, EPA has revised Method 524.3 to allow for the use of nitrogen as a purge gas. Although labs have been doing this for years with other methods, substantial reductions in helium consumption can now be obtained when analyzing purgeable organic compounds in water by GC-MS using the revised method. By switching to nitrogen purge gas using Method 524.4, you can save an impressive 490 mL of helium per sample, which translates into a 68% reduction in helium consumption (Table I). Saving nearly 0.5 L of helium per sample quickly adds up to considerable cost savings and also insulates labs from the impact of fluctuating helium availability.

In addition to reducing helium consumption by using nitrogen purge gas, employing an Rtx-VMS column for this analysis allows all Method 524.4 criteria to be easily met. The Rtx-VMS column is recommended for purge-and-trap GC-MS analysis of VOCs by Method 524.4 because the selectivity of this column provides ample separation between all critical compounds. As shown in Figure 1, good resolution is obtained for target analytes, including o-xylene and styrene, as well as 1,1,1-trichloroethane and carbon tetrachloride. The Rtx-VMS column is listed in Method 524.4 and was used by the EPA to establish method performance [1]. No interference from overlapping peaks was observed, and the small bore 0.25 mm column results in higher column efficiency and improved separations.

Labs analyzing purgeable organic compounds in water can save money and reduce helium dependence by using Method 524.4 with nitrogen purge gas and an Rtx-VMS column.

References


[1] U.S. Environmental Protection Agency, Method 524.4, Measurement of Purgeable Organic Compounds in Water by Gas Chromatography/Mass Spectrometry Using Nitrogen Purge Gas, May 2013. <http://water.epa.gov/scitech/drinkingwater/labcert/upload/815R13002.pdf> (accessed December 19, 2013).

Table I: Using nitrogen purge gas allows labs to reduce helium use by 68%.

Instrument	Analysis Step	He Volume	N ₂ Volume
P&T	Purge	—	440
	Dry Purge	—	50
	Desorb	30.9	—
GC-MS	Split Vent*	140	—
	Carrier	13.5	—
	Septum Purge	45	—
	Total Gas Volume		719.4
	Volume of Helium Saved		490
	Percent of Helium Saved		68%

*Gas saver is turned to 10 mL/min @ 1 min.

Figure 1: The Rtx-VMS column is specifically designed for separating purgeable organic compounds, which means it provides excellent resolution of the Method 524.4 VOCs analyzed here using nitrogen purge gas.



For EPA 524.4 analysis, Restek recommends...



EPA 524.3 Reference Standards

(5 separate mixes, including a MegaMix standard)

Each ampul sold separately.

Cat.# 30013: 524.3 VOA MegaMix (69 components)

Allyl chloride (3-chloropropene) (107-05-1)
tert-Amyl ethyl ether (TAE) (919-94-8)
tert-Amyl methyl ether (TAME) (994-05-8)
 Benzene (71-43-2)
 Bromobenzene (108-86-1)
 Bromochloromethane (74-97-5)
 Bromodichloromethane (75-27-4)
 Bromoform (75-25-2)
tert-Butanol (TBA) (75-65-0)
n-Butylbenzene (104-51-8)
sec-Butylbenzene (135-98-8)
tert-Butylbenzene (98-06-6)
 Carbon disulfide (75-15-0)
 Carbon tetrachloride (56-23-5)
 Chlorobenzene (108-90-7)
 Chloroform (67-66-3)
 1-Chlorobutane (butyl chloride) (109-69-3)
 2-Chlorotoluene (95-49-8)
 4-Chlorotoluene (106-43-4)
 Dibromochloromethane (124-48-1)
 1,2-Dibromo-3-chloropropane (96-12-8)
 Dibromomethane (74-95-3)
 1,2-Dibromoethane (EDB) (106-93-4)
 1,2-Dichlorobenzene (95-50-1)
 1,3-Dichlorobenzene (541-73-1)
 1,4-Dichlorobenzene (106-46-7)
 1,1-Dichloroethane (75-34-3)
 1,2-Dichloroethane (107-06-2)
 1,1-Dichloroethene (75-35-4)
cis-1,2-Dichloroethene (156-59-2)
trans-1,2-Dichloroethene (156-60-5)
 1,2-Dichloropropane (78-87-5)
 1,3-Dichloropropane (142-28-9)
 1,1-Dichloropropene (563-58-6)
cis-1,3-Dichloropropene (10061-01-5)
trans-1,3-Dichloropropene (10061-02-6)
 Diethyl ether (ethyl ether) (60-29-7)
 Diisopropyl ether (DIPE) (108-20-3)
 Ethylbenzene (100-41-4)
 Ethyl-*tert*-butyl ether (ETBE) (637-92-3)
 Ethyl methacrylate (97-63-2)
 Hexachlorobutadiene (87-68-3)
 Hexachloroethane (67-72-1)
 Iodomethane (methyl iodide) (74-88-4)
 Isopropylbenzene (cumene) (98-82-8)
 4-Isopropyltoluene (*p*-cymene) (99-87-6)
 Methyl acetate (79-20-9)
 Methyl-*tert*-butyl ether (MTBE) (1634-04-4)

Methylene chloride (dichloromethane) (75-09-2)
 Naphthalene (91-20-3)
 Pentachloroethane (76-01-7)
n-Propylbenzene (103-65-1)
 Styrene (100-42-5)
 Tetrachloroethene (127-18-4)
 1,1,1,2-Tetrachloroethane (630-20-6)
 1,1,2,2-Tetrachloroethane (79-34-5)
 Tetrahydrofuran (109-99-9)
 Toluene (108-88-3)
 1,2,3-Trichlorobenzene (87-61-6)
 1,2,4-Trichlorobenzene (120-82-1)
 1,1,1-Trichloroethane (71-55-6)
 1,1,2-Trichloroethane (79-00-5)
 Trichloroethene (79-01-6)
 1,2,3-Trichloropropane (96-18-4)
 1,2,4-Trimethylbenzene (95-63-6)
 1,3,5-Trimethylbenzene (108-67-8)
m-Xylene (108-38-3)
o-Xylene (95-47-6)
p-Xylene (106-42-3)

Cat.# 30014: 524.3 Gas Calibration Mix (7 components)

Bromomethane (methyl bromide) (74-83-9)
 1,3-Butadiene (106-99-0)
 Chlorodifluoromethane (CFC-22) (75-45-6)
 Chloromethane (methyl chloride) (74-87-3)
 Dichlorodifluoromethane (CFC-12) (75-71-8)
 Trichlorofluoromethane (CFC-11) (75-69-4)
 Vinyl chloride (75-01-4)

Cat.# 30015: 524.3 Internal Standard Mix (3 components)

Chlorobenzene-*d*5 (3114-55-4)
 1,4-Dichlorobenzene-*d*4 (3855-82-1)
 1,4-Difluorobenzene (540-36-3)

Cat.# 30016: 524.3 Surrogate Standard (3 components)

1-Bromo-4-fluorobenzene (BFB) (460-00-4)
 Methyl-*d*3-*tert*-butyl ether (29366-08-3)
 1,2-Dichlorobenzene-*d*4 (2199-69-1)

Cat.# 30017: 524.3 Internal Standard/Surrogate Mix (6 components)

1-Bromo-4-fluorobenzene (BFB) (460-00-4)
 Methyl-*d*3-*tert*-butyl ether (29366-08-3)
 Chlorobenzene-*d*5 (3114-55-4)
 1,2-Dichlorobenzene-*d*4 (2199-69-1)
 1,4-Dichlorobenzene-*d*4 (3855-82-1)
 1,4-Difluorobenzene (540-36-3)

Gas calibration and internal standards are also available at www.restek.com

Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
524.3 Gas Calibration Mix							
2000 µg/mL each in P&T methanol, 1 mL/ampul	Yes	6 months	36 months	Ambient	0 °C or colder	ea.	30014
524.3 Internal Standard Mix							
2000 µg/mL each in P&T methanol, 1 mL/ampul	Yes	6 months	60 months	Ambient	0 °C or colder	ea.	30015
524.3 Internal Standard/Surrogate Mix							
2000 µg/mL each in P&T methanol, 1 mL/ampul	Yes	6 months	60 months	Ambient	0 °C or colder	ea.	30017
524.3 Surrogate Standard							
2000 µg/mL each in P&T methanol, 1 mL/ampul	Yes	6 months	60 months	Ambient	0 °C or colder	ea.	30016
524.3 VOA MegaMix							
2000 µg/mL each in P&T methanol, 1 mL/ampul	Yes	6 months	36 months	Ambient	0 °C or colder	ea.	30013

Rtx-VMS Columns (fused silica)

proprietary Crossbond phase

- Application-specific columns for analyzing volatile organic pollutants by GC-MS including Methods TO-15, TMS, and EPA 8260.
- Complete separation of U.S. EPA Method 8260 compounds in less than 10 minutes.
- Stable to 260 °C.
- No known equivalent phases.

ID	df	Length	Temp. Limits	qty.	Similar to Part #	cat.#
0.25 mm	1.40 µm	30 m	-40 to 240/260 °C	ea.	Agilent 122-1534	19915
	1.40 µm	60 m	-40 to 240/260 °C	ea.	Agilent 122-1564	19916



Parker Nitrogen Gas Generators

- Produces ultra-pure nitrogen (up to 99.9999%).
- Requires only a compressed air source and 120 volt AC power.
- Typical applications include GC carrier gas, make-up gas, and low-flow sample concentrators.
- Maintenance kits include replacement filters.
- 12-month warranty from date of purchase.

Description	Model #	Certification/Compliance	qty.	cat.#
Nitrogen Generator	UHPN2-1100 (ultra-high purity zero grade)	CE	ea.	20697



20697

ordering notes

International power cords are available.
Contact Customer Service to order.

Parker Annual Maintenance Kits for Nitrogen Gas Generators

Description	Model #	Includes	Vendor cat.#	qty.	cat.#
Annual Maintenance Kit	for UHPN2-1100, 76-94, 76-96	1st, 2nd, and 3rd stage prefilters (1 ea.), and final filter (1)	MK7694	kit	21655
	for HPN2-1100, 76-92	1st and 2nd stage prefilters (1 ea.) and, final filter (1)	MK7692	kit	21649



21655



Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2022 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# EVSS1938A-UNV