



Avoid Mycotoxin Quantitation Errors When Using Stable Isotope Dilution Assay (SIDA)

There are significant benefits to using the stable isotope dilution assay (SIDA) calibration approach for mycotoxin analyses [1], and one of the biggest is being able to use a single calibration curve for multiple matrix types. Compared to making individual matrix-matched calibration standards for each sample type, the time savings offered by the SIDA technique can make up for the relatively high price of isotopically labeled internal standards (IS), especially for high-throughput food safety laboratories that frequently deal with a wide variety of samples. However, labeled standards are not always commercially available for every mycotoxin of interest. Due to these cost and availability issues, labs may be interested in using a practice that is common among many other types of methods that use internal standards: using the same internal standard to calibrate a group of analytes. That approach may seem logical, but in the case of mycotoxin SIDA methods, it can lead to significant errors.

Table I compares the observed concentrations of four different mycotoxins present in two maize flour reference materials to the concentrations reported by the reference material supplier. For the three mycotoxins that were quantified using their corresponding isotopically labeled internal standards, the agreement between the two values is excellent. However, in the case of the mycotoxin zearalenone, we did not have its analogous labeled internal standard, so we quantified it using $^{13}\text{C}_{17}$ -aflatoxin G1 because it elutes nearby. As the data show, the agreement between the observed results and those reported by the reference material supplier is very poor, which illustrates that similar chromatographic retention alone is not enough to predict the effects of sample preparation and/or matrix-related changes in ionization efficiency. A more detailed description of Restek's research on the LC-MS/MS analysis of mycotoxins in various foods comparing SIDA to matrix-matched calibration is published in a peer-reviewed study [2].

In light of these results, it is strongly recommended that only matching isotopically labeled internal standards be used for quantification in mycotoxin SIDA methods. If matching internal standards are not available, other calibration approaches, such as matrix-matched calibration, should be used.

Table I: Mycotoxin SIDA calibration should only be used for matching analyte and IS pairs.

Reference Material	Analyte	IS	Measured Concentration (ng/g), n=3	Assigned Concentration (ng/g)	Percent Accuracy (RSD, %)
TET017RM	Deoxynivalenol	$^{13}\text{C}_{15}$ -Deoxynivalenol	1867.9 \pm 37.36	1971 \pm 195	94.8 (2)
TET017RM	Aflatoxin B1	$^{13}\text{C}_{17}$ -Aflatoxin B1	8.68 \pm 0.434	9.49 \pm 0.85	91.4 (5)
TET017RM	Ochratoxin A	$^{13}\text{C}_{20}$ -Ochratoxin A	4.48 \pm 0.134	4.81 \pm 0.75	93.2 (3)
TET017RM	Zearalenone	$^{13}\text{C}_{17}$ -Aflatoxin G1	31.26 \pm 2.19 ^a	231 \pm 25	13.5 (7) ^a
T04301Q	Deoxynivalenol	$^{13}\text{C}_{15}$ -Deoxynivalenol	639.7 \pm 19.19	649 \pm 222	98.6 (3)
T04301Q	Aflatoxin B1	$^{13}\text{C}_{17}$ -Aflatoxin B1	8.69 \pm 0.348	9.21 \pm 4.05	94.4 (4)
T04301Q	Ochratoxin A	$^{13}\text{C}_{20}$ -Ochratoxin A	2.81 \pm 0.197	3.03 \pm 1.33	92.6 (7)
T04301Q	Zearalenone	$^{13}\text{C}_{17}$ -Aflatoxin G1	16.2 \pm 0.810 ^a	138.5 \pm 59.6	11.7 (5) ^a

* Results quantified using a nonmatched labeled IS.

References

1. K. Zhang, K., M.R. Schaab, G. Southwood, E.R. Tor, L.S. Aston, W. Song, B. Eitzer, S. Majumdar, T. Lapainis, H. Mai, K. Tran, A. El-Demerdash, V. Vega, Y. Cai, J.W. Wong, A.J. Krynsky, T.H. Begley, A collaborative study: determination of mycotoxins in corn, peanut butter, and wheat flour using stable isotope dilution assay (SIDA) and liquid chromatography-tandem mass spectrometry (LC-MS/MS), J. of Agric. Food Chem. 65 (33) (2017) 7138–7152. DOI: 10.1021/acs.jafc.6b04872 <https://pubs.acs.org/doi/10.1021/acs.jafc.6b04872>
2. D. Li, J.A. Steimling, J.D. Konschnik, S. Grossman, T. Kahler, Quantitation of mycotoxins in four food matrices comparing stable isotope dilution assay (SIDA) with matrix-matched calibration methods by LC-MS/MS, J. AOAC Int. (2019) DOI: 10.5740/jaoacint.19-0028 <https://doi.org/10.5740/jaoacint.19-0028>

Raptor Biphenyl LC Columns (USP L11)



ID	Length	qty.	cat.#
1.8 µm Particles			
2.1 mm	30 mm	ea.	9309232
	50 mm	ea.	9309252
	100 mm	ea.	9309212
	150 mm	ea.	9309262
3.0 mm	50 mm	ea.	930925E
	100 mm	ea.	930921E
2.7 µm Particles			
2.1 mm	30 mm	ea.	9309A32
	50 mm	ea.	9309A52
	100 mm	ea.	9309A12
	150 mm	ea.	9309A62
3.0 mm	30 mm	ea.	9309A3E
	50 mm	ea.	9309A5E
	100 mm	ea.	9309A1E
	150 mm	ea.	9309A6E
4.6 mm	30 mm	ea.	9309A35
	50 mm	ea.	9309A55
	100 mm	ea.	9309A15
	150 mm	ea.	9309A65
5 µm Particles			
2.1 mm	50 mm	ea.	9309552
	100 mm	ea.	9309512
	150 mm	ea.	9309562
	30 mm	ea.	930953E
3.0 mm	50 mm	ea.	930955E
	100 mm	ea.	930951E
	150 mm	ea.	930956E
	50 mm	ea.	9309555
4.6 mm	100 mm	ea.	9309515
	150 mm	ea.	9309565
	250 mm	ea.	9309575




Want even better performance when analyzing mycotoxins? Check out Inert LC columns at www.restek.com/inert



RESTEK
Pure Chromatography

Raptor EXP Guard Cartridges



Protect your investment and extend the life of our already-rugged LC columns and change guard column cartridges by hand without breaking fluid connections—no tools needed!

EXP Direct Connect Holder

Description	qty.	cat.#
EXP Direct Connect Holder for EXP Guard Cartridges (includes hex-head fitting & 2 ferrules)	ea.	25808

Intellectual Property: optimizetech.com/patents



Raptor EXP Guard Column Cartridges

Description	Particle Size	Size	qty.	cat.#
Raptor Biphenyl EXP Guard Column Cartridge	UHPLC	5 x 2.1 mm	3-pk.	9309U0252
	UHPLC	5 x 3.0 mm	3-pk.	9309U0253
	2.7 µm	5 x 2.1 mm	3-pk.	9309A0252
	2.7 µm	5 x 3.0 mm	3-pk.	9309A0253
	2.7 µm	5 x 4.6 mm	3-pk.	9309A0250
	5 µm	5 x 2.1 mm	3-pk.	930950252
	5 µm	5 x 3.0 mm	3-pk.	930950253
	5 µm	5 x 4.6 mm	3-pk.	930950250

Maximum cartridge pressure: 1034 bar/15,000 psi* (UHPLC), 600 bar/8700 psi (2.7 µm); 400 bar/5800 psi (5 µm)

* For maximum lifetime, recommended maximum pressure for UHPLC particles is 830 bar/12,000 psi.

Intellectual Property: optimizetech.com/patents

Raptor Inert Biphenyl HPLC Columns

- Inert LC column technology reduces nonspecific binding of chelating analytes, enabling sensitive analysis and smooth integration of peaks.
- Ideal for the analysis of metal-sensitive compounds, such as mycotoxins.
- Increased response and analyte recovery, allowing lower detection limits.
- Improved peak shape without additional passivation or mobile phase additives.
- Part of Restek's Raptor Biphenyl column line featuring 2.7 µm SPP core-shell silica.

ID	Length	Particle Size	Units	Cat.#
2.1 mm	100	2.7 µm	ea.	9309A12-T
3.0 mm	100	2.7 µm	ea.	9309A1E-T
2.1 mm	50	2.7 µm	ea.	9309A52-T
3.0 mm	50	2.7 µm	ea.	9309A5E-T

Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# FFAR3036B-UNV