



Analysis of Furan and Alkylfurans in Food Commodities by SPME Arrow GC-MS

By Dr. Nathaly Reyes-Garcés

Related Products

Restek PAL SPME Arrow; Restek PAL SPME Arrow GC-Specific Conversion Kits; Accessories for SPME Arrow; Rxi-624Sil MS; Topaz 1.8 mm ID Straight/SPME Inlet Liner; Headspace Crimp Vials; Magnetic Seals w/Septa; Headspace Screw-Thread Vials; Magnetic Screw-Thread Caps

Abstract

In this study, SPME Arrow sample preparation was tested for the analysis of furan and alkylfurans in baby formula and coffee samples. Optimized conditions were established for GC-MS analysis using a PAL CTC autosampler and an Rxi-624Sil MS column. Results demonstrated the method may be suitable for diverse sample types containing both high and low concentrations of furan and alkylfurans. This approach could allow food safety labs to improve productivity by using SPME Arrow as a single sample preparation approach compared to methods that use headspace alone or SPME for low-concentration samples and headspace for high-concentration samples.

Introduction

Analysis of furan and alkylfurans is emerging as a growth area in food safety testing due to increasing recognition of their potential health risks. These compounds form through a diverse array of chemical reactions that occur during heating processes, such as roasting or sterilization, so exposure can occur through a wide range of foods. Exposure has not been well characterized so, in early 2022, the European Commission adopted a recommendation calling for member states to monitor furan, 2-methylfuran, and 3-methylfuran in particular foods [1]. The agency also advised monitoring non-methylfurans, such as 2,5-dimethylfuran, 2-pentylfuran, and 2-ethylfuran, if methods are available that have been demonstrated to be reliable for this purpose.

Furan and alkylfurans are volatile compounds that can be analyzed by static headspace (HS) and/or solid phase microextraction (SPME) in combination with GC-MS. Static headspace GC-MS is typically used with samples that contain high levels of furan and alkylfurans, such as coffee, but it may not offer adequate sensitivity for low concentration samples. SPME is more sensitive and better for low-level samples, but traditional SPME fibers are very fragile and not well suited to high-throughput environments. To meet regulatory expectations efficiently, food safety labs need a simple, fast, cost-effective workflow that performs well for a wide range of analyte concentrations in different sample types.

For this reason, we decided to evaluate the performance of SPME Arrows to determine if a single approach could be used for the analysis of furan and alkylfurans. SPME Arrows are second-generation SPME devices that are mechanically enhanced with a sheath and inner rod to prevent breakage; they also have higher phase volumes to provide more loading capacity (Figure 1). In this study, based on an initial comparison of sample preparation techniques, we developed an optimized SPME Arrow GC-MS method for the analysis of furan and alkylfurans across a range of concentrations in both coffee and baby formula.

Figure 1: SPME Arrows have higher loading capacity than traditional fibers and are also sheathed in metal to prevent breakage.

SPME
Arrows



SPME
Fibers



Experimental

Standard Solutions

For each target analyte, 1 mg/mL solutions were prepared in methanol by spiking pure standards using a glass syringe. Then, 100 ppm stock solutions, one containing all the target analytes and one containing all the internal standards, were prepared in methanol. For the calibration curve covering a low concentration range, 0.1 and 0.5 µg/mL working solutions were prepared by spiking 10 and 50 µL of the 100 ppm stock, respectively, in 10 mL of LC-MS grade water. An internal standard solution of 1 µg/mL was prepared by adding 100 µL of the 100 ppm stock solution to 9.9 mL of water. For the analysis of coffee samples, a high-concentration calibration curve was used. For this purpose, aqueous working solutions at 10 and 25 µg/mL, for target analytes and internal standards, respectively, were made by spiking the original individual methanolic solutions directly into water. A 0.5 µg/mL solution was prepared by diluting 500 µL of the 10 µg/mL solution in water until reaching a volume of 10 mL.

Sample Preparation

Baby formula and instant coffee samples were purchased at a local grocery store. Preliminary experiments comparing headspace, traditional SPME fibers, and SPME Arrows were run to establish a sample preparation procedure that would be suitable for a wide range of analyte concentrations. Parameters tested included coating chemistry, incubation time and temperature, and extraction time [2,3,4].

Based on the results of the initial sample preparation experiments, SPME Arrow was selected as the best approach and was used under the optimized conditions shown in Table I for subsequent method evaluation experiments. All samples and calibration solutions were prepared in 20 mL glass vials and extracted and analyzed by automated headspace-SPME using a 120 µm wide carbon range/PDMS Restek PAL SPME Arrow (cat.# 27487) assembled in a PAL CTC autosampler.

Table I: Sample Preparation Conditions for Baby Formula and Coffee Samples

Sample	Baby Formula	Coffee
Amount		0.5 g
Volume of 30% NaCl solution added	10 mL	5 mL
Volume of internal standard solution added	50 µL (1 µg/mL solution)	40 µL (25 µg/mL solution)
Incubation time		10 min
Incubation and extraction temperature		50 °C
Agitation		250 rpm
HS-SPME extraction time	10 min	1 min
Desorption temperature		280 °C
Desorption time		1 min

Calibration Curves

Two calibration curves were prepared in 30% sodium chloride solution: one for low concentration analysis of furan and alkylfurans in baby formula, and one for the quantitation of high concentrations of these target analytes in coffee. For the low calibration range, vials were filled with 10 mL of sodium chloride solution, whereas for the high calibration range 5 mL of solution were added. Analytes were spiked by adding different volumes of working solutions to the calibration vials, as shown in Table II. Internal standard solutions were added to each calibration vial (50 µL of the 1 µg/mL solution for the low concentration calibration curve, and 40 µL of the 25 µg/mL solution for the high concentration calibration curve). To account for variation in recoveries due to matrix differences, it was essential to use isotopically labeled analogues for almost every analyte (except for 2- and 3-methylfuran because 2-methylfuran-d6 worked for both compounds). Calibration curves were constructed by plotting analyte area/internal standard average area ratios (n=2) versus spiked concentration.

Table II: Calibration Solutions

	Low Concentration Calibration Curve		High Concentration Calibration Curve	
	ng of analyte in vial	µL of working solution added (solution concentration)	ng of analyte in vial	µL of working solution added (solution concentration)
Level 1	1.25	12.5 (0.1 µg/mL)	25	50 (0.5 µg/mL)
Level 2	2.5	25 (0.1 µg/mL)	50	100 (0.5 µg/mL)
Level 3	5	10 (0.5 µg/mL)	100	10 (10 µg/mL)
Level 4	10	20 (0.5 µg/mL)	200	20 (10 µg/mL)
Level 5	20	40 (0.5 µg/mL)	500	50 (10 µg/mL)
Level 6	50	100 (0.5 µg/mL)	1000	100 (10 µg/mL)
Level 7	100	200 (0.5 µg/mL)	2000	200 (10 µg/mL)
Level 8	150	300 (0.5 µg/mL)	4000	400 (10 µg/mL)
Level 9			8000	800 (10 µg/mL)

Method Verification Samples

To evaluate method performance, baby formula and coffee samples were spiked at two concentration levels based on the standard method performance requirements listed by the AOAC [6]. Baby formula samples were spiked at 5 and 50 µg/kg (n=3), and coffee samples were spiked at 1000 and 4000 µg/kg (n=3). Sample blanks spiked only with internal standards were also analyzed (n=3).

Chromatographic Method

GC-MS analysis of furan and alkylfurans was performed on an Agilent 7890B GC paired with a 5977B MSD. GC-MS method details are shown in Tables III and IV. To enable the analysis of highly concentrated coffee samples with the wide carbon range SPME Arrows, three method parameters were modified: the GC split was set to 100:1 (with a split flow rate of 140 mL/min); the extraction time was set to 1 min; and the volume of sodium chloride solution was adjusted to 5 mL. Note that a 1.8 mm ID inlet liner was used because the wider bore diameter is necessary to accommodate the width of SPME Arrows.

Table III: GC-MS Instrument Conditions

Column	Rxi-624Sil MS, 30 m, 0.25 mm ID, 1.40 µm (cat.# 13868)
Injection Mode	Split (10:1 for baby formula and 100:1 for coffee)
Liner	Topaz 1.8 mm ID SPME/straight liner (cat.# 23280)
Inj. Temp.	280 °C
Split Flow	• Baby formula: 14.0 mL/min (10:1) • Coffee: 140 mL/min (100:1)
Purge Flow	5 mL/min
Oven	35 °C (hold 3 min) to 75 °C by 8 °C/min, then to 200 °C by 25°C/min (hold 1 min)
Carrier Gas	He, constant flow
Flow Rate	1.40 mL/min
Analyzer	MS (quadrupole)
Acquisition Type	SIM
Ionization Mode	El (70 eV)
Transfer Line Temp.	280 °C
Source Temp.	325 °C
Quadrupole Temp.	200 °C
Solvent delay	2.2 min

Table IV: MS Parameters (SIM mode)

Segment Starting Time (min)	Compound (t _r , min)	Ions	Dwell Time (ms)
2.2	Furan (2.447)	39	50
		68*	
	Furan-d4 (2.428)	42	
		72*	
4.2	2-Methylfuran (4.536)	53	30
		81	
		82*	
	3-Methylfuran (4.846)	53	30
		81	
		82*	
	2-Methylfuran-d6 (4.464)	58	30
		88*	
6.6	2-Ethylfuran (7.100)	53	30
		81*	
		96	
	2-Ethylfuran-d5 (7.001)	55	30
		101*	
	2,5-Dimethylfuran (7.243)	67	30
		95*	
	2,5-Dimethylfuran-d3 (7.179)	84	30
10.6	2-Pentylfuran (11.570)	81	30
		138*	30
	2-Pentylfuran-d11 (11.501)	83	30
		149*	30


*Quantifier ions

Results and Discussion

Comparison of Headspace, SPME Fiber, and SPME Arrow Sample Preparation

A comparison of various sampling strategies was conducted at the beginning of the study. SPME Arrow (wide carbon range coating); traditional SPME fibers (wide carbon range coating); and headspace sampling were evaluated in terms of area counts after sampling from 30% sodium chloride solution (10 mL) spiked with all analytes (20 ng of each analyte). As can be seen in Figure 2, the SPME Arrow provided significantly higher responses for all target analytes compared to headspace or traditional SPME fibers. Compared to SPME fibers, SPME Arrows have the additional advantage of a stainless-steel sheath and stabilizing inner rod, which prevent breakage during use. Based on the response data and rugged, steel construction, SPME Arrow was chosen for the development of a single sample preparation GC-MS method.

Figure 2: Analyte response was significantly higher with SPME Arrows than with headspace or traditional SPME fibers.




SPME conditions: 2 min extraction, 5 min incubation at 40 °C, 1 min desorption at 280 °C, agitation at 250 rpm.

Headspace conditions: 1000 uL sample volume, 20 min incubation, 120 sec purge time, 80 °C incubation temperature, 150 °C HS tool temperature.

Chromatographic Performance

An Rxi-624Sil MS column (30 m, 0.25 mm ID, 1.40 μm [cat.# 13868]) was used for the analysis of furan and alkylfurans because it is highly selective for volatile organics and works well for both SPME and headspace methods. The column and optimized conditions provided good chromatographic separation of all compounds, and example chromatograms for a solvent standard and spiked samples are provided in Figures 3-5. Of particular note, chromatographic separation between 2,5-dimethylfuran and 2-ethylfuran (isomeric compounds) and between 2,5-dimethylfuran-d3 and 2-ethylfuran was achieved, which allows accurate identification and quantitation (Figure 4). Chromatographic separation is important because 2,5-dimethylfuran-d3 and 2-ethylfuran share most of the ions of their mass spectra (except for ions m/z 95 and 96) and, if they are not separated, 2,5-dimethylfuran-d3 can contribute to the most abundant qualifier ions of 2-ethylfuran (m/z 81 and 53). All compounds were adequately separated in a 14-minute analysis.

Figure 3: Analysis of furan and alkylfurans with internal standards in a 10 ng/mL standard. Good separation of 2-ethylfuran, 2,5-dimethylfuran-d3, and 2,5-dimethylfuran was achieved on an Rxi-624Sil MS column using the optimized gradient conditions.



Peaks	t_r (min)	Ion (Quantifier)	Ion (Qualifier)	Dwell Time (ms)
1. Furan-d4	2.428	72	42	50
2. Furan	2.447	68	39	50
3. 2-Methylfuran-d6	4.464	88	58	30
4. 2-Methylfuran	4.536	82	53	30
5. 3-Methylfuran	4.846	82	53	30
6. 2-Ethylfuran-d5	7.001	101	55	30
7. 2-Ethylfuran	7.100	81	96	30
8. 2,5-Dimethylfuran-d3	7.179	99	84	30
9. 2,5-Dimethylfuran	7.243	95	67	30
10. 2-Pentylfuran-d11	11.501	149	83	30
11. 2-Pentylfuran	11.570	138	81	30

GC_GN1218

Column
Standard/Sample

Diluent:
N/A

Conc.:
10 ng/mL

Injection
split (split ratio 10:1)

Liner:
Topaz 1.8 mm ID straight/SPME inlet liner (cat.# 23280)

Inj. Temp.:
280 °C

Split Vent Flow Rate:
14 mL/min

Oven
Oven Temp.: 35 °C (hold 3 min) to 75 °C at 8 °C/min to 200 °C at 25 °C/min (hold 1 min)

Carrier Gas
He, constant flow

Flow Rate:
1.4 mL/min

Detector
MS

Transfer Line Temp.:
280 °C

Analyzer Type:
Quadrupole

Source Temp.:
325 °C

Quad Temp.:
200 °C

Electron Energy:
70 eV

Tune Type:
PFTBA

Ionization Mode:
EI

Instrument
Agilent 7890B GC & 5977B MSD

Sample Preparation
Data was collected by extracting via HS-SPME from a 20 mL vial (cat.# 23083) capped with a magnetic screw-thread cap (cat.# 23091). The vial contained 10 mL of sodium chloride solution spiked with 100 ng of each analyte and internal standard. A Restek PAL SPME Arrow (120 μm carbon wide range [WR]/PDMS; cat.# 27487) was used.

Rxi-624Sil MS, 30 m, 0.25 mm ID, 1.40 μm (cat.# 13868)

Oven Temp.: 35 °C (hold 3 min) to 75 °C at 8 °C/min to 200 °C at 25 °C/min (hold 1 min)

Carrier Gas

He, constant flow

Flow Rate:

1.4 mL/min

Detector

MS

Transfer Line Temp.:

280 °C

Analyzer Type:

Quadrupole

Source Temp.:

325 °C

Quad Temp.:

200 °C

Electron Energy:

70 eV

Tune Type:

PFTBA

Ionization Mode:

EI

Instrument

Agilent 7890B GC & 5977B MSD


Sample Preparation
Data was collected by extracting via HS-SPME from a 20 mL vial (cat.# 23083) capped with a magnetic screw-thread cap (cat.# 23091). The vial contained 10 mL of sodium chloride solution spiked with 100 ng of each analyte and internal standard. A Restek PAL SPME Arrow (120 μm carbon wide range [WR]/PDMS; cat.# 27487) was used.

SPME Arrow sampling conditions: 10 min extraction time, 10 min incubation at 50 °C, 1 min desorption at 280 °C, agitation at 250 rpm.

Notes

A furan/alkylfurans standard (cat.# 33334) is now available.

Figure 4: Analysis of Furan and Alkylfurans in a Baby Formula Sample Spiked at 5 µg/kg



Peaks	t _r (min)	Ion (Quantifier)	Ion (Qualifier)	Dwell Time (ms)
1. Furan-d4	2.428	72	42	50
2. Furan	2.447	68	39	50
3. 2-Methylfuran-d6	4.464	88	58	30
4. 2-Methylfuran	4.536	82	53	30
5. 3-Methylfuran	4.846	82	53	30
6. 2-Ethylfuran-d5	7.001	101	55	30
7. 2-Ethylfuran	7.100	81	96	30
8. 2,5-Dimethylfuran-d3	7.179	99	84	30
9. 2,5-Dimethylfuran	7.243	95	67	30
10. 2-Pentylfuran-d11	11.501	149	83	30
11. 2-Pentylfuran	11.570	138	81	30

Column Rxi-624Sil MS, 30 m, 0.25 mm ID, 1.40 µm (cat.# 13868)

Standard/Sample

Diluent:

Conc.: 5 µg/kg

Injection split (split ratio 10:1)

Liner: Topaz 1.8 mm ID straight/SPME inlet liner (cat.# 23280)

Inj. Temp.: 280 °C

Split Vent Flow Rate: 14 mL/min

Oven

Oven Temp.: 35 °C (hold 3 min) to 75 °C at 8 °C/min to 200 °C at 25 °C/min (hold 1 min)

Carrier Gas

Flow Rate: He, constant flow

Detector

Transfer Line Temp.: 280 °C

Analyzer Type: Quadrupole

Source Temp.: 325 °C

Quad Temp.: 200 °C

Electron Energy: 70 eV

Tune Type: PFTBA

Ionization Mode: EI

Instrument Agilent 7890B GC & 5977B MSD

Sample Preparation


Data was collected by extracting via HS-SPME from a 20 mL vial (cat.# 23083) capped with a magnetic screw-thread cap (cat.# 23091). The vial contained 0.5 g of baby formula and 10 mL of sodium chloride solution spiked with 2.5 ng of each analyte (5 ng/g final concentration). Internal standards were spiked at 100 ng/g. A Restek PAL SPME Arrow (120 µm carbon wide range [WR]/PDMS; cat.# 27487) was used.

SPME Arrow sampling conditions: 10 min extraction time, 10 min incubation at 50 °C, 1 min desorption at 280 °C, agitation at 250 rpm.

Notes

A furan/alkylfurans standard (cat.# 33334) is now available.

Figure 5: Analysis of Furan and Alkylfurans in an Instant Coffee Sample Spiked at 1000 µg/kg



Peaks	t _r (min)	Ion (Quantifier)	Ion (Qualifier)	Dwell Time (ms)
1. Furan-d4	2.428	72	42	50
2. Furan	2.447	68	39	50
3. 2-Methylfuran-d6	4.464	88	58	30
4. 2-Methylfuran	4.536	82	53	30
5. 3-Methylfuran	4.846	82	53	30
6. 2-Ethylfuran-d5	7.001	101	55	30
7. 2-Ethylfuran	7.100	81	96	30
8. 2,5-Dimethylfuran-d3	7.179	99	84	30
9. 2,5-Dimethylfuran	7.243	95	67	30
10. 2-Pentylfuran-d11	11.501	149	83	30
11. 2-Pentylfuran	11.570	138	81	30

Column Rxi-624Sil MS, 30 m, 0.25 mm ID, 1.40 µm (cat.# 13868)

Standard/Sample

Diluent:

N/A

Conc.:

1000 µg/kg

Injection

split (split ratio 100:1)

Liner:

Topaz 1.8 mm ID straight/SPME inlet liner (cat.# 23280)

Inj. Temp.:

280 °C

Split Vent Flow Rate:

140 mL/min

Oven

Oven Temp.:

35 °C (hold 3 min) to 75 °C at 8 °C/min to 200 °C at 25 °C/min (hold 1 min)

Carrier Gas

He, constant flow

Flow Rate:

1.4 mL/min

Detector

MS

Transfer Line Temp.:

280 °C

Analyzer Type:

Quadrupole

Source Temp.:

325 °C

Quad Temp.:

200 °C

Electron Energy:

70 eV

Tune Type:

PFTBA

Ionization Mode:

EI

Instrument

Agilent 7890B GC & 5977B MSD

Sample Preparation

Data was collected by extracting via HS-SPME from a 20 mL vial (cat.# 23083) capped with a magnetic screw-thread cap (cat.# 23091). The vial contained 0.5 g of instant coffee and 5 mL of sodium chloride solution spiked with 500 ng of each analyte (1000 ng/g final concentration). Internal standards were spiked at 2000 ng/g. A Restek PAL SPME Arrow (120 µm carbon wide range [WR]/PDMS; cat.# 27487) was used.

SPME Arrow sampling conditions: 1 min extraction time, 10 min incubation at 50 °C, 1 min desorption at 280 °C, agitation at 250 rpm.

Notes


A furan/alkylfurans standard (cat.# 33334) is now available.

Linearity


Calibration curves were constructed for all the target analytes in sodium chloride solution (30%) using appropriate deuterated analogues for each compound. The use of isotopically labeled internal standards was essential to account for differences in analyte recovery from sample matrices versus from spiked sodium chloride solution. The quantitation strategy followed in this work was a slightly modified version of the method described by Frank et al. in 2020 [5]. The low concentration calibration curve ranged from 1.25 to 150 ng of analytes in vial (8 calibration points) while the high concentration calibration curve ranged from 25 to 8000 ng of analyte in vial (9 calibration points). As can be seen in Figure 6, all analytes showed good linearity within the tested ranges, and all the coefficients of determination were above 0.99. For pentylfuran, the low concentration calibration curve ranged from 1.25 to 100 ng, and the high concentration calibration curve was linear from 25 to 2000 ng in vial.

Figure 6: Calibration Curves for the Analysis of Furans and Alkylfurans via SPME Arrow

Low Calibration Range



High Calibration Range



Method Verification in Baby Formula and Coffee Samples

Baby formula and coffee both had significant concentrations of some of the analytes of interest in the blank samples. Baby formula was analyzed using the low concentration calibration curve, whereas instant coffee was analyzed using the high concentration calibration curve. Table V and VI present a summary of the collected results. All results for baby formula passed the 80-110% recovery criteria, except for 3-methylfuran, which had a recovery of 113% in the low spike. For coffee, all low concentration spikes passed, but some high bias was observed for the high concentration spikes. While further investigation into sample size, sample handling, and dilution could be done in future experiments, recoveries were still fairly close to the 110% upper limit. This indicates that with some modification SPME Arrow is likely to be an effective approach for high concentration samples as well as low concentration samples.

Table V: Analysis of Furan and Alkylfurans in Baby Formula

Analyte	Blank (n=3)		Low concentration (n=3), 5 µg/kg*		High concentration (n=3), 50 µg/kg*	
	Concentration, µg/kg	RSD, %	Accuracy, %	RSD, %	Accuracy, %	RSD, %
Furan	16	1	110	2	94	4
2-Methylfuran	60	2	97	2	100	3
3-Methylfuran	-	-	113	3	107	6
2-Ethylfuran	67	2	93	3	105	5
2,5-Dimethylfuran	-	-	101	4	97	14
2-Pentylfuran	219	3	87	11	108	12

*Accuracy was determined as follows: ((measured concentration – concentration in blank)/spiked concentration)*100

Table VI: Analysis of Furan and Alkylfurans in Instant Coffee

Analyte	Blank (n=3)		Low concentration (n=3), 1000 µg/kg*		High concentration (n=3), 4000 µg/kg*	
	Concentration, µg/kg	RSD, %	Accuracy, %	RSD, %	Accuracy, %	RSD, %
Furan	394	10	87	2	125	3
2-Methylfuran	843	10	94	3	116	2
3-Methylfuran	96	10	87	2	119	3
2-Ethylfuran	29	11	93	2	115	5
2,5-Dimethylfuran	46	11	98	2	120	4
2-Pentylfuran	-	-	83	3	83	11

*Accuracy was determined as follows: ((measured concentration – concentration in blank)/spiked concentration)*100

Conclusion

Analysis of furan and alkylfurans is complicated by a diverse array of sample types that can contain the target compounds at a wide range of concentrations. Headspace and traditional SPME fibers can be used for high and low concentration samples respectively, but this study sought to establish a single approach that could be adapted for all analyte levels. The method developed here, based on SPME Arrow sample preparation, was evaluated in matrices spiked at both low and high concentration levels: 5 and 50 µg/kg for baby formula, and 1000 and 4000 µg/kg for coffee. Satisfactory results in terms of linearity, accuracy, and precision were obtained in most cases. It is noted that for the analysis of highly concentrated samples, different split (1:100); extraction time (1 min); and sodium chloride solution (30%) (5 mL) conditions should be used, and further optimization may be necessary. Use of this SPME Arrow GC-MS method may allow labs to improve productivity by using a single sample preparation technique that can be adjusted to accommodate both high and low concentration samples.

References

1. European Commission recommendation 2022/495, On monitoring the presence of furan and alkylfurans in food, 2022. https://eur-lex.europa.eu/legal-content/EN/TXT/HTML/?uri=CELEX-32022H0495&from=EN#ntr1-L_2022100EN.01006001-E0001
2. N. Reyes-Garcés, Analysis of furan and alkylfurans in food (part 1): choosing GC-MS conditions and a sample preparation technique, ChromaBLOGraphy, Restek Corporation, 2022. <https://www.restek.com/chromablography/chromablography/analysis-of-furans-and-alkylfurans-in-food-samples-part-1/>
3. N. Reyes-Garcés, Analysis of furan and alkylfurans in food (part 2): optimizing experimental conditions for the analysis of food matrices using the SPME Arrow, ChromaBLOGraphy, Restek Corporation, 2022. <https://www.restek.com/chromablography/chromablography/analysis-of-furan-and-alkylfurans-in-food-samples-part-2/>
4. N. Reyes-Garcés, Analysis of furan and alkylfurans in food samples (part 3): Adjusting the experimental conditions for the analysis of samples with high concentrations of analytes and preparing calibration curves, ChromaBLOGraphy, Restek Corporation, 2022.
5. N. Frank, M. Dubois, J.F. Huertas Pérez, Detection of furan and five alkylfurans, including 2-pentylfuran, in various food matrices, J. Chromatogr. A, 1622 (2020) 461119. <https://doi.org/10.1016/j.chroma.2020.461119>
6. AOAC SMPR 2019.004, Standard method performance requirements (SMPRs) for furan and alkyl furans in coffee, baby foods, infant formula, cereals, and fruit juices, final version: October 3, 2019, AOAC International, https://www.aoac.org/wp-content/uploads/2019/09/SMPR_Furans_v6-For-Comment.pdf

Restek PAL Smart SPME Arrow

- Rugged, stainless-steel construction ensures longer lifetimes.
- Faster extraction means higher sample throughput.
- Better sensitivity allows lower LODs.
- Now with Smart technology—each SPME Arrow is equipped with a unique Smart chip containing parameters, ranges, and usage history.

Recommended maximum GC inlet pressure is 50 psi or less.

All PAL SPME Smart Arrows have 20 mm of phase bonded onto stainless steel.

PAL system compatibility: Smart SPME Arrows are fully backwards compatible with any generation of PAL3 Systems. Additional hardware will be also required, including the SPME Arrow Tool, Heatex Stirrer, Arrow Conditioning Module, and Agitator Module with the perforated lid. The SPME Arrow is not supported on the HTS/Combi PAL systems.



Ordering Note: The SPME Arrow used in these analyses (cat.# 27487) has been replaced with a Smart SPME Arrow equivalent (cat.# 28903-1).

The new Smart technology offers the same application, performance, and life-time as the products you've come to rely on with the added feature of a smart chip that can keep track of parameters, ranges, and usage history.

Due to the relatively large diameter of PAL Smart SPME Arrows, you must modify the GC inlet using an instrument-specific conversion kit from Restek prior to use. Converted inlets are compatible with all standard injection techniques (SPME Arrow, liquid syringe, headspace, etc.); no need to switch inlets after installation.

Catalog No.	Product Name	Color	Diameter	Material	Recommended Analytes	Units
28903-1	PAL Smart SPME Arrow 1.10 mm: Carbon-WR/PDMS, Phase Thickness 120 µm, Phase Length 20 mm, Light Blue, ea.	Light Blue	1.1 mm	120 µm Carbon-WR/PDMS (Carbon Wide Range/ Polydimethylsiloxane)	Highly volatile, 30–225 g/mol*	ea.
28903-3	PAL Smart SPME Arrow 1.10 mm: Carbon-WR/PDMS, Phase Thickness 120 µm, Phase Length 20 mm, Light Blue, 3-pk.	Light Blue	1.1 mm	120 µm Carbon-WR/PDMS (Carbon Wide Range/ Polydimethylsiloxane)	Highly volatile, 30–225 g/mol*	3-pk.

*These molecular weight ranges are a reasonable approximation; however, end users should verify suitability for their specific application.

More Smart SPME Arrows are available for a wide variety of applications.

Visit www.restek.com/SPME

Accessories for SPME Arrows

Description	Instrument	qty.	cat.#
Injector adaptor cup	for Agilent GC 6890/7890/8890 Split/Splitless Injector	ea.	27496
	for Shimadzu GC 2010 and GC 2030 Split/Splitless Injector	ea.	27497
	for Thermo GC TRACE 1300/1310, 1600/1610 Split/Splitless Injector	ea.	27498
Injection port weldment	for Thermo GC TRACE Ultra Split/Splitless Injector	ea.	27499
	for Shimadzu GC 2010 Split/Splitless Injector	ea.	27500
	for Shimadzu GC 2030 Split/Splitless Injector	ea.	27893
Needle guide/septum nut	for Shimadzu GC 2010 Split/Splitless Injector	ea.	27501
Split/splitless weldment; large, canister-type filter	for Agilent GC 6890 Split/Splitless Injector	ea.	27502
Septum nut for split/splitless weldments	for Agilent GC 6890/7890 Split/Splitless Injector	ea.	27503
Split/splitless weldment and septum nut	for Agilent GC 7890 Split/Splitless Injector	ea.	27504
	for Agilent 8890 Split/Splitless Injector	ea.	28933
	for Thermo GC TRACE 1300/1310, 1600/1610 Split/Splitless Injector	ea.	27505
Septum cap	for Thermo GC TRACE 1300/1310, 1600/1610 Split/Splitless Injector	ea.	27506
Liner cap/septum holder	for Thermo GC TRACE 1300/1310, 1600/1610 Split/Splitless Injector	ea.	27507
Septum holder and support	for Thermo GC TRACE Ultra Split/Splitless Injector	ea.	27507
Liner cap	for Thermo GC TRACE Ultra Split/Splitless Injector	ea.	27508





Restek PAL SPME Arrow GC-Specific Conversion Kits

Conversion kits designed for use with 1.1 and 1.5 mm Restek PAL SPME Arrows. Compatible with Smart SPME Arrows and non-smart SPME Arrows.

Due to the relatively large diameter of Restek PAL SPME Arrows, you must modify the GC inlet using an instrument-specific conversion kit prior to use. Converted inlets are compatible with all standard injection techniques (SPME Arrow, liquid syringe, headspace, etc.); so there is no need to switch inlets after installation.

Restek PAL SPME Arrow GC-Specific Conversion Kits are easy to install. For the weldment, install according to your instrument's owner's manual.

Description	Includes	Instrument	qty.	cat.#
Restek PAL SPME Arrow Conversion Kit	Topaz 1.8 mm ID straight/SPME inlet liner, 5-pk. (cat.# 23280); Thermolite Plus septa, 3-pk. (cat.# 23864); Split/splitless weldment; large, canister-type filter (cat.# 27502); Septum nut for 6890 split/splitless weldments (cat.# 27503); Injector adaptor cup (cat.# 27496)	for Agilent 6890 Split/Splitless Injector (for canister-type filters)	kit	27492
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23280); 1.1 mm Microseal (cat.# 23232); Split/Splitless Weldment; Large Canister Type Filter (cat.# 27502); Adaptor Cup (cat.# 27496); Nut (cat.# 23228)	for Agilent 6890 Split/Splitless Injector (for canister-type filters)	kit	27356
Restek PAL SPME Arrow Conversion Kit with 1.5 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23280); 1.5 mm Microseal (cat.# 23233); Split/Splitless Weldment and Large Canister Type Filter (cat.# 27502); Adaptor Cup (cat.# 27496); Nut (cat.# 23228)	for Agilent 6890 Split/Splitless Injector (for canister-type filters)	kit	27361
Restek PAL SPME Arrow Conversion Kit	Topaz 1.8 mm ID straight/SPME inlet liner, 5-pk. (cat.# 23280); Thermolite Plus septa, 3-pk. (cat.# 23864); Agilent split/splitless weldment and septum nut (cat.# 27504); Injector adaptor cup (cat.# 27496)	for Agilent 7890 Split/Splitless Injector	kit	27493
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23280); 1.1 mm Microseal (cat.# 23232); Agilent Weldment (cat.# 27504); Adaptor Cup (cat.# 27496); Nut (cat.# 23228)	for Agilent 7890 Split/Splitless Injector	kit	27357
Restek PAL SPME Arrow Conversion Kit with 1.5 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23280); 1.5 mm Microseal (cat.# 23233); Agilent Weldment (cat.# 27504); Adaptor Cup (cat.# 27496); Nut (cat.# 23228)	for Agilent 7890 Split/Splitless Injector	kit	27362
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23280); 1.1 mm Microseal (cat.# 23232); Agilent Weldment (cat.# 28933); Adaptor Cup (cat.# 27496); Nut (cat.# 23228)	for Agilent 8890 Split/Splitless Injector	kit	28934
Restek PAL SPME Arrow Conversion Kit	Topaz 1.8 mm ID straight/SPME inlet liner, 5-pk. (cat.# 23279); Thermolite Plus septa, 3 pk. (cat.# 23872); Injection port weldment (cat.# 27500); Needle guide/septum nut (cat.# 27501); Injector adaptor cup (cat.# 27497)	for Shimadzu GC-2010 Split/Splitless Injector (not compatible with SE or Plus models)	kit	27491
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23279); 1.1 mm Microseal (cat.# 23232); Port Weldment (cat.# 27500); Adaptor Cup (cat.# 27497); Adaptor Kit (cat.# 23229)	for Shimadzu GC-2010 Split/Splitless Injector (not compatible with SE or Plus models)	kit	27355
Restek PAL SPME Arrow Conversion Kit with 1.5 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23279); 1.5 mm Microseal (cat.# 23233); Port Weldment (cat.# 27500); Adaptor Cup (cat.# 27497); Adaptor Kit (cat.# 23229)	for Shimadzu GC-2010 Split/Splitless Injector (not compatible with SE or Plus models)	kit	27360
Restek PAL SPME Arrow Conversion Kit	Topaz 1.8 mm ID straight/SPME inlet liner, 5-pk. (cat.# 23279); Thermolite Plus septa, 3 pk. (cat.# 23872); Injection port Shimadzu 2030 GC; Needle guide/septum nut (cat.# 27501); Injector adaptor cup (cat.# 27497)	for Shimadzu GC-2030 Split/Splitless Injector	kit	27886
Restek PAL SPME Arrow Conversion Kit	Topaz 1.8 mm ID straight/SPME inlet liner, 5-pk. (cat.# 23278); Premium non-stick BTO septa, 3-pk. (cat.# 27090); Septum cap (cat.# 27505); Liner cap/septum holder (cat.# 27506); Injector adaptor cup (cat.# 27498)	for Thermo TRACE 1300/1310, 1600/1610 Split/Splitless Injector	kit	27494
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23278); 1.1 mm Microseal (cat.# 23232); Liner Cap/Septum Holder (cat.# 27506); Adaptor Cup (cat.# 27498); Nut (cat.# 23230)	for Thermo TRACE 1300/1310, 1600/1610 Split/Splitless Injector	kit	27358
Restek PAL SPME Arrow Conversion Kit with 1.5 mm Merlin Microseal	Topaz 1.8 mm ID Straight/SPME inlet liner, 5-pk. (cat.# 23278); 1.5 mm Microseal (cat.# 23233); Liner Cap/Septum Holder (cat.# 27506); Adaptor Cup (cat.# 27498); Adaptor Kit (cat.# 23230)	for Thermo TRACE 1300/1310, 1600/1610 Split/Splitless Injector	kit	27363
Restek PAL SPME Arrow Conversion Kit	2.0 mm ID straight inlet liner, 5-pk. (cat.# 22267); Premium non-stick BTO septa, 3-pk. (cat.# 27096); Septum holder and support (cat.# 27507); Liner cap (cat.# 27508); Injector adaptor cup (cat.# 27499)	for Thermo TRACE Ultra Split/Splitless Injector	kit	27495
Restek PAL SPME Arrow Conversion Kit with 1.1 mm Merlin Microseal	2.0 mm ID straight inlet liner, 5-pk. (cat.# 22267); 1.1 mm Microseal (cat.# 23232); Liner Cap (cat.# 27508); Adaptor Cup (cat.# 27499); Adaptor Kit (cat.# 23231)	for Thermo TRACE Ultra Split/Splitless Injector	kit	27359
Restek PAL SPME Arrow Conversion Kit with 1.5 mm Merlin Microseal	2.0 mm ID straight inlet liner, 5-pk. (cat.# 22267); 1.5 mm Microseal (cat.# 23233); Liner Cap (cat.# 27508); Adaptor Cup (cat.# 27499); Adaptor Kit (cat.# 23231)	for Thermo TRACE Ultra Split/Splitless Injector	kit	27364

Rxi-624Sil MS Columns (fused silica)

midpolarity Crossbond phase

- Low-bleed, high-thermal stability column—maximum temperatures up to 300–320 °C.
- Inert—excellent peak shape for a wide range of compounds.
- Selective—G43 phase highly selective for volatile organics and residual solvents, great choice for USP<467>.
- Manufactured for column-to-column reproducibility—well suited for validated methods.




ID	df	Length	Temp. Limits	qty.	Similar to Part #	cat.#
0.25 mm	1.40 µm	30 m	-20 to 300/320 °C	ea.	Agilent 122-1334; Phenomenex 7HG-G005-27	13868

ordering notes

Custom lengths and film thicknesses available. Contact Technical Service or your local Restek representative.

SAVE MONEY! Get six columns for the price of five. Contact Customer Service or your local Restek representative for details!





Headspace Crimp Vials, 20 mm

Description	Type	Volume	Color	Size	qty.	Similar to Part #	cat.#
Headspace Vial, Flat Bottom	20 mm Crimp-Top	6 mL	Clear	22 x 38 mm	100-pk.		21166
	20 mm Crimp-Top	6 mL	Clear	22 x 38 mm	1000-pk.		21167
	20 mm Crimp-Top	10 mL	Clear	23 x 46 mm	100-pk.	Agilent 8010-0029	24683
	20 mm Crimp-Top	10 mL	Clear	23 x 46 mm	1000-pk.		24684
Headspace Vial, Rounded Bottom	20 mm Crimp-Top	10 mL	Clear	23 x 46 mm	100-pk.		21164
	20 mm Crimp-Top	10 mL	Clear	23 x 46 mm	1000-pk.		21165
Headspace Vial, Rounded Bottom, Deactivated	20 mm Crimp-Top	10 mL	Clear	23 x 46 mm	1000-pk.		21165-221
Headspace Vial, Flat Bottom	20 mm Crimp-Top	20 mL	Clear	23 x 75 mm	100-pk.	Agilent 8010-0033	24685
	20 mm Crimp-Top	20 mL	Clear	23 x 75 mm	1000-pk.		24686
Headspace Vial, Rounded Bottom	20 mm Crimp-Top	20 mL	Clear	23 x 75 mm	100-pk.		21162
	20 mm Crimp-Top	20 mL	Clear	23 x 75 mm	1000-pk.		21163
Headspace Vial, Flat Bottom	20 mm Crimp-Top	27 mL	Clear	30 x 60 mm	100-pk.		21160
	20 mm Crimp-Top	27 mL	Clear	30 x 60 mm	1000-pk.		21161



Headspace Screw-Thread Vials, 18 mm

Description	Type	Volume	Color	Size	qty.	cat.#
Headspace Vial, Rounded Bottom	18-425 Screw-Thread	20 mL	Clear	22 x 75 mm	100-pk.	23082
	18-425 Screw-Thread	20 mL	Clear	22 x 75 mm	1000-pk.	23083
	18-425 Screw-Thread	20 mL	Amber	22 x 75 mm	100-pk.	23086
	18-425 Screw-Thread	20 mL	Amber	22 x 75 mm	1000-pk.	23087
	18-425 Screw-Thread	10 mL	Clear	22 x 45 mm	100-pk.	23084
	18-425 Screw-Thread	10 mL	Clear	22 x 45 mm	1000-pk.	23085
	18-425 Screw-Thread	10 mL	Amber	22 x 45 mm	100-pk.	23088
	18-425 Screw-Thread	10 mL	Amber	22 x 45 mm	1000-pk.	23089

Caps not included (sold separately).

Magnetic Screw-Thread Caps, 18 mm

Description	Type	Cap Size	Septa Material	qty.	cat.#
Magnetic Caps and Septa for SPME,	Screw-Thread	18-425	Blue PTFE/Silicone, 1.5 mm thick	100-pk.	23090
	Screw-Thread	18-425	Blue PTFE/Silicone, 1.5 mm thick	1000-pk.	23091
	Screw-Thread	18-425	Red PTFE/Silicone, 1.9 mm thick	100-pk.	23092
Magnetic Caps and Septa,	Screw-Thread	18-425	Red PTFE/Silicone, 1.9 mm thick	1000-pk.	23093
	Screw-Thread	18-425	PTFE/Red Chlorobutyl	100-pk.	23094
	Screw-Thread	18-425	PTFE/Red Chlorobutyl	1000-pk.	23095



Magnetic Seals w/Septa, 20 mm, preassembled

Description	Type	Cap Size	Color	Material	Septa Material	qty.	cat.#
Magnetic Seals w/Septa, with 5 mm Hole	Crimp-Top	20 mm	Gold	steel	PTFE/Silicone	100-pk.	22833
	Crimp-Top	20 mm	Gold	steel	PTFE/Silicone	1000-pk.	22834
	Crimp-Top	20 mm	Gold	steel	PTFE/Silicone	100-pk.	22831
Magnetic Seals w/Septa, with 8 mm Hole	Crimp-Top	20 mm	Gold	steel	PTFE/Silicone	1000-pk.	22832
	Crimp-Top	20 mm	Gold	steel	PTFE/Butyl	100-pk.	22835
	Crimp-Top	20 mm	Gold	steel	PTFE/Butyl	1000-pk.	22836
	Crimp-Top	20 mm	Red/Silver	aluminum/steel	PTFE/Silicone	100-pk.	22441
BiMetal Magnetic Seals w/Septa, with 8 mm Hole	Crimp-Top	20 mm	Red/Silver	aluminum/steel	PTFE/Silicone	1000-pk.	22442
	Crimp-Top	20 mm	Blue/Silver	aluminum/steel	PTFE/Silicone	100-pk.	22443
	Crimp-Top	20 mm	Blue/Silver	aluminum/steel	PTFE/Silicone	1000-pk.	22444



The 5 mm hole is compatible with the following systems: Carlo Erba HS500/HS800, CTC 500, Fisons HS500/HS800, and Bruker/Varian/Chrompack 9020/25. The 8 mm hole is compatible with the Combi Pal. BiMetal and standard magnetic crimp cap seals are more rigid than aluminum crimp cap seals and require more force to seal and decap. To prevent user fatigue, Restek recommends using high-power electronic crimpers/decappers instead of manual or battery-powered crimpers/decappers when using BiMetal and standard magnetic crimp cap seals.

Topaz GC Inlet Liners

Topaz GC inlet liners feature revolutionary technology and inertness to deliver you the next level of True Blue Performance:

- **Deactivation**—unbelievably low breakdown for accurate and precise low-level GC analyses.
- **Reproducibility**—unbeatable manufacturing controls and QC testing for superior reliability across compound classes.
- **Productivity**—unparalleled cleanliness for maximized GC uptime and lab throughput.
- **100% Satisfaction**—if a liner doesn't perform to your expectations, we will replace it or credit your account.*

Patented

* 100% SATISFACTION GUARANTEE: If your Topaz inlet liner does not perform to your expectations for any reason, simply contact Restek Technical Service or your local Restek representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty

Topaz 1.8 mm ID Straight/SPME Inlet Liner

for Agilent GCs equipped with split/splitless inlets

Description	Length	ID	OD	Deactivation	Material	qty	cat.#
1.8 mm ID Straight/SPME Liner	78.5 mm	1.8 mm	6.5 mm	Premium	Borosilicate Glass	5-pk.	23280

Looking for furan/alkylfuran standards?

A furan/alkylfurans multicomponent standard (cat.# 33334) is now available.

