




Featured Application: Polar Pesticides in Spinach on Raptor Polar X by LC-MS/MS





Fast and Rugged Direct Analysis of Polar Pesticides in Spinach

- No time-consuming derivatization, ion-pairing steps, or long waits for column equilibration.
- Unique stationary phase provides good retention and efficient elution, resulting in better peak shape.
- Robust column performance ensures consistent chromatography and reliable results.

The direct analysis of polar pesticides, such as glyphosate, its main metabolite AMPA, and glufosinate, is difficult because they are not retained well on reversed phase columns, and chelation with metal surfaces in the LC-MS/MS can significantly reduce response. Derivatization and ion-pairing reagents are often used to improve performance, but they require additional sample preparation time and can produce unfavorable side effects (precipitates, long equilibration times, system contamination, etc.).

Using a Raptor Polar X column is a better solution because the novel stationary phase offers stronger interaction with polar anionic compounds. It features hybrid retention modes (HILIC and ion exchange) and provides strong retention yet allows analytes to be quickly eluted as symmetrical peaks through simple mobile phase changes. In this analysis of polar pesticides in spinach, excellent chromatographic performance is seen for all three compounds with a total analysis time of just eight minutes. Further, peak shapes and retention times remain stable over hundreds of matrix injections, and a simple guard column change restores full performance. This fast, simple, and effective workflow is a good alternative for food labs interested in the direct analysis of polar pesticides without the need for complex and time-consuming procedures.



Reference Standards	Sample Handling	Analytical Column	Maintenance/Accessories
			
Glyphosate cat.# 32426 & 32427	50 mL centrifuge tubes cat.# 25846 Syringe filter 30 mm x 0.22 µm cat.# 23984	2.7 µm, 30 x 2.1 mm Raptor Polar X column cat.# 9311A32	Passivation solution cat.# 32475 Raptor Polar X guard column cat.# 9311A0252 EXP direct connect holder cat.# 25808



Glyphosate Standard

Glyphosate (N-(phosphonomethyl)glycine) (1071-83-6)

Description	CAS #	Conc. in Solvent	cat.#
Glyphosate	1071-83-6	1,000 µg/mL in DI water, 1 mL/ampul	32426 (ea.)
	1071-83-6	1,000 µg/mL in DI water, 5 mL/ampul	32427 (ea.)



25846

LC Passivation Solution

Methylenediphosphonic Acid (Medronic Acid) (1984-15-2)

Description	CAS #	Conc. in Solvent	cat.#
Methylenediphosphonic Acid (Medronic Acid)	1984-15-2	1,760 µg/mL, Methanol (HPLC grade)/Water (50:50), 1mL/ampul	32475 (ea.)

Empty Centrifuge Tubes, Polypropylene

Description	qty.	cat.#
Empty 50 mL Centrifuge Tube, Polypropylene w/Blue Cap	50-pk.	25846



23984

Syringe Filters with Luer Lock Inlet

- Luer lock inlet offers leak-tight syringe connection.
- Autoclavable to 121 °C for 15 minutes.
- Quantity break pricing for greater savings.

Description	Color	Diameter	Porosity	qty.	cat.#
PTFE (polytetrafluoroethylene)	White	30 mm	0.22 µm	100-pk.	23984

PTFE—hydrophobic applications. Syringe filters are for laboratory use only.

Raptor Polar X LC Columns

- Reliably analyze a wide variety of polar analytes (acidic, basic, and neutral) without time-consuming derivatization or complex ion pairing.
- Switch between HILIC and ion-exchange retention modes with simple mobile phase changes and short equilibration times.
- 2.7 μm Raptor core-shell particles provide UHPLC-like speed and efficiency on all makes and models of LC systems.
- Ideal for increasing sensitivity and selectivity in LC-MS analyses.

ID	Length	qty.	cat.#
2.7 μm Particles			
2.1 mm	30 mm	ea.	9311A32
	50 mm	ea.	9311A52
	100 mm	ea.	9311A12



Raptor EXP Guard Column Cartridges

- Free-Turn architecture lets you change cartridges by hand without breaking inlet/outlet fluid connections—no tools needed.
- Patented titanium hybrid ferrules can be installed repeatedly without compromising high-pressure seal.
- Auto-adjusting design provides ZDV (zero dead volume) connection to any 10-32 female port.
- Guard column cartridges require EXP direct connect holder (cat.# 25808).

Description	Particle Size	Size	qty.	cat.#
Raptor Polar X EXP Guard Column Cartridge	2.7 μm	5 x 2.1 mm	3-pk.	9311A0252

Maximum cartridge pressure: 600 bar/8,700 psi (2.7 μm)

Hybrid Ferrule U.S. Patent No. 8201854, EXP Holders U.S. Patent No. 8696902, EXP2 Wrench U.S. Patent No. D766055. Other U.S. and Foreign Patents Pending. The EXP, Free-Turn, and the Opti- prefix are registered trademarks of Optimize Technologies, Inc.



9311A0252

EXP Direct Connect Holder

Description	qty.	cat.#
EXP Direct Connect Holder for EXP Guard Cartridges (includes hex-head fitting & 2 ferrules)	ea.	25808

Maximum holder pressure: 20,000 psi (1,400 bar)

Hybrid Ferrule U.S. Patent No. 8201854, EXP Holders U.S. Patent No. 8696902, EXP2 Wrench U.S. Patent No. D766055. Other U.S. and Foreign Patents Pending. The EXP, Free-Turn, and the Opti- prefix are registered trademarks of Optimize Technologies, Inc.



25808

Related Products



26431

Bluestem Glass Solvent Filter

Prevent the particulates and microbial growth in your LC solvents from entering your instrument with the new Restek Bluestem glass solvent filter.

Description	qty.	cat.#
Frit Adaptor, PTFE	4-pk.	26392
Glass Solvent Filter, 15 µm frit	ea.	26431



25322

Survival Kit for HPLC, PEEK

For start-up and maintenance in all HPLC systems.

The PEEK Survival Kit is an invaluable parts kit that contains tubing, fittings, and tools essential for setting up and maintaining your HPLC system: PEEK tubing, connectors, and elbows; PTFE tubing; a tubing cutter and extra blades; a ValvTool wrench; open-end wrenches; and more.

Description	qty.	cat.#
Survival Kit for HPLC	kit	25322

*Kit contains 1 wrench, replacement (cat.# 20110) is a 2-pk.



23243

Limited-Volume 2.0 mL, 9 mm Screw-Thread Polypropylene Vials

- Fit all 2.0 mL, 12 x 32 mm, vial-based autosamplers.
- Compatible with all 9 mm screw-thread caps.
- PTFE-free—ideal for PFAS analysis (e.g., EPA 537) and other PFAS-sensitive methods.

Note: Polypropylene vials and caps prevent sample contamination from PTFE coated septa. However, since polypropylene caps do not reseal, evaporation occurs after injection. Multiple injections from the same vial are therefore not possible.

Description	Type	Volume	Color	Size	qty.	cat.#
Limited-Volume 2.0 mL,	9 mm Screw-Thread	700 µL	Clear	12 x 32 mm	100-pk.	23243
9 mm Screw-Thread Polypropylene Vials	9 mm Screw-Thread	700 µL	Clear	12 x 32 mm	1,000-pk.	23246



24669

2.0 mL, 9 mm Short-Cap, Screw-Vial Closures (Polypropylene, preassembled)

Description	Type	Cap Size	Color	Septa Material	qty.	cat.#
Short Screw Caps	Screw-Thread	9-425	Blue	PTFE/Silicone	100-pk.	24485
	Screw-Thread	9-425	Blue	PTFE/Silicone	1,000-pk.	24486
	Screw-Thread	9-425	Green	PTFE/Silicone	100-pk.	24487
	Screw-Thread	9-425	Green	PTFE/Silicone	1,000-pk.	24488
	Screw-Thread	9-425	Yellow	PTFE/Silicone	100-pk.	24493
	Screw-Thread	9-425	Yellow	PTFE/Silicone	1,000-pk.	24494
	Screw-Thread	9-425	Mixed	PTFE/Silicone	500-pk.	24669