



Applicazione descritta: Antibiotici veterinari multiclasse in LC-MS/MS sulla Raptor C18

Una sola analisi, una sola colonna, meno di 9 minuti per oltre 60 antibiotici multiclasse


- Separazione dei picchi altamente efficiente e tempi di analisi rapidi.
- Applicabilità universale per l'analisi residuale degli antibiotici, con possibilità di ottimizzare l'analisi di singole classi di composti ai fini della quantificazione:
 - Macrolidi, lincosamidi, e streptogramine (Figura 1)
 - Amfenicoli e tetracicline (Figura 2)
 - Chinoloni (Figura 3)
 - Penicilline, cefalosporine e tetracicline (Figura 4)
 - Sulfamidici (Figura 5)
 (Per gli ionofori è stata utilizzata la Raptor bifenilica. Figura 6)

L'uso di antibiotici sugli animali destinati alla produzione alimentare è un problema di salute e sicurezza pubblica in quanto potrebbero svilupparsi batteri farmacosensibili. Molti paesi dell'Unione europea e il Canada hanno vietato l'uso di questi farmaci per scopi non terapeutici, e gli Stati Uniti stanno implementando una politica volta a ridurre l'uso di antibiotici di rilievo clinico come promotori della crescita. Per regolamentare l'uso corretto degli antibiotici veterinari, la statunitense FDA ha fissato i limiti massimi dei residui (LMR) per diversi tessuti animali e prodotti alimentari (21 CFR Part 556). Per operare in conformità a queste disposizioni è necessario impiegare un metodo analitico sensibile, efficiente e affidabile per diverse classi di antibiotici, e la colonna LC Raptor C18 è la scelta ideale.

Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion	Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion
1. Desacetyl cephalirin	0.70	150	382.03	111.92	124.21	31. Sarafloxacin	2.59	10	386.20	342.20	368.15
2. Sulfanilamide	0.85	200	172.98	93.07	75.23	32. Difloxacin	2.65	10	400.23	356.17	299.13
3. Amoxicillin	0.92	100	366.24	349.10	208.07	33. Cefazolin	2.75	100	455.10	323.06	295.09
4. Cephalirin	1.36	50	424.17	292.08	124.14	34. Spiramycin	2.96	200	843.64	540.36	699.48
5. Tildipirosin	1.38	200	734.59	561.45	204.15	35. Pirlimycin	3.05	20	411.32	363.18	327.21
6. Desfuroyl ceftiofur cysteine disulfide	1.40	300	549.16	183.02	126.00	36. Chlortetracycline	3.08	25	479.27	154.07	371.06
7. Lincomycin	1.50	50	407.32	359.23	389.28	37. Sulfachlorpyridazine	3.12	20	285.05	156.03	108.09
8. Sulfadiazine	1.57	20	251.18	156.04	92.08	38. Gamithromycin	3.28	100	777.63	619.52	601.45
9. Cefquinome	1.73	200	529.19	134.10	125.12	39. Sulfadoxine	3.39	10	311.17	156.03	108.09
10. Ampicillin	1.78	50	350.19	106.07	160.06	40. Sulfamethoxazole	3.46	20	254.18	155.98	147.06
11. Sulfathiazole	1.79	10	256.16	156.03	92.08	41. Cefoperazone	3.52	100	646.26	143.07	148.02
12. Marbofloxacin	1.81	10	363.20	72.11	320.10	42. Florfenicol*	3.55	200	356.10	336.02	184.98
13. Cefalexin	1.82	100	348.10	158.05	174.05	43. Sulfathoxypridazine	3.56	20	295.17	267.07	156.03
14. Sulfapyridine	1.86	10	250.13	156.10	92.08	44. Tilmicosin	3.64	100	869.72	696.50	522.42
15. Norfloxacin	1.96	20	320.23	276.20	233.13	45. Sulfisoxazole	3.76	20	268.14	156.03	113.10
16. Ofloxacin	1.98	10	362.21	318.20	261.15	46. Oxolinic acid	3.94	5	262.10	244.06	215.96
17. Sulfamerazine	2.00	20	265.08	156.03	92.08	47. Chloramphenicol*	4.00	200	321.16	151.99	257.04
18. Cefalonium	2.01	100	459.16	337.03	123.10	48. Ceftiofur	4.11	50	524.14	241.08	125.24
19. Oxytetracycline	2.02	25	461.27	426.15	443.32	49. Erythromycin	4.31	25	734.64	576.40	558.38
20. Ciprofloxacin	2.04	20	332.18	288.22	245.15	50. Sulfadimethoxine	4.36	10	311.17	156.09	108.09
21. Cefacetile	2.09	300	362.07	258.08	178.01	51. Sulfaguanoxaline	4.42	20	301.18	156.04	108.02
22. Tulathromycin A	2.11	100	806.65	577.42	420.31	52. Tylosin**	4.67	100	916.62	772.49	598.36
23. Tetracycline	2.21	25	445.28	154.07	427.32	53. Penicillin G	5.07	100	335.18	176.07	160.07
24. Danofloxacin	2.23	20	358.22	340.16	314.21	54. Flumequine	5.34	5	262.15	244.11	202.03
25. Enrofloxacin	2.32	10	360.29	316.22	245.13	55. Penicillin V	5.56	100	351.10	160.06	114.07
26. Orbifloxacin	2.35	10	396.22	352.17	226.12	56. Oxacillin	5.99	100	402.15	160.05	114.06
27. Thiamphenicol*	2.38	200	354.16	290.04	184.98	57. Virginiamycin M1	6.13	50	526.43	508.31	355.10
28. Sulfamethazine	2.39	10	279.23	186.08	124.08	58. Tylvalosin	6.19	50	1042.71	814.46	640.39
29. Sulfamethizole	2.52	10	271.17	156.02	108.02	59. Cloxacillin	6.39	100	436.15	277.06	160.05
30. Sulfamethoxypridazine	2.56	10	281.14	156.03	126.07	60. Nafcillin	6.60	25	415.19	199.09	171.06
						61. Dicloxacillin	6.96	100	470.11	160.05	311.02

*Acquired in negative ion mode

Consulta il cromatogramma e le condizioni a pagina 2



LC_FS0506


Column Raptor C18 (cat.# 9304A12)
Dimensions: 100 mm x 2.1 mm ID
Particle Size: 2.7 µm
Pore Size: 90 Å
Guard Column: Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)
Temp.: 35 °C
Sample
Diluent: Water
Conc.: 5–300 ng/mL
Inj. Vol.: 2 µL
Mobile Phase
A: 0.1% Formic acid in water
B: 0.1% Formic acid in acetonitrile

Time (min)	Flow (mL/min)	%A	%B
0.00	0.4	90	10
4.50	0.4	65	35
7.00	0.4	45	55
7.01	0.4	90	10
9.00	0.4	90	10

Detector MS/MS
Ion Mode: ESI+/ESI-
Mode: Scheduled MRM
Instrument UHPLC
Notes


1. Positive and negative polarity data were collected simultaneously from a single injection.
 2. Amphenicol compounds (chloramphenicol, thiamphenicol, and flornfenicol) were detected with negative polarity.
 3. The MRM was scheduled at +/- 20 to 30 seconds for each analyte.
 4. Multiclass antibiotics include penicillin, cephalosporin, tetracycline, sulfonamide, macrolide, lincosamide, streptogramin, amphenicol, and quinolone.
- **The retention time for Tylosin is noted in the peak list; however, it was not included in the chromatogram.

Figura 1: Macrolidi, lincosamidi e streptogramine in LC-MS/MS sulla Raptor C18



Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion	Column
1. Tildipirosin	1.01	200	734.59	561.45	204.15	Raptor C18 (cat.# 9304A12)
2. Lincomycin	1.11	50	407.32	359.23	389.28	Dimensions: 100 mm x 2.1 mm ID
3. Tulathromycin A	1.19	100	806.65	577.42	420.31	Particle Size: 2.7 µm
4. Spiramycin	1.41	200	843.64	540.36	699.48	Pore Size: 90 Å
5. Pirlmycin	1.49	20	411.32	363.18	327.21	Guard Column: Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)
6. Gamithromycin	1.50	100	777.63	619.52	601.45	Temp.: 40 °C
7. Tilmicosin	1.58	100	869.72	696.50	522.42	Sample Diluent: Water
8. Erythromycin	1.80	25	734.64	576.40	558.38	Conc.: 20–200 ng/mL
9. Tylosin	1.87	100	916.62	772.49	598.36	Inj. Vol.: 2 µL
10. Tylvalosin	2.30	50	1042.71	814.46	640.39	Mobile Phase
11. Virginiamycin M1	2.45	50	526.43	508.31	355.10	A: 0.1% Formic acid in water
						B: 0.1% Formic acid in acetonitrile
						Time (min) Flow (mL/min) %A %B
						0.00 0.4 90 10
						3.00 0.4 20 80
						3.01 0.4 90 10
						5.00 0.4 90 10
						Detector MS/MS
						Ion Mode: ESI+
						Mode: MRM
						Instrument UHPLC

Figura 2: Amfenicoli e tetracicline in LC-MS/MS sulla Raptor C18



Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion	Column
1. Oxytetracycline	1.28	25	461.27	426.15	443.32	Raptor C18 (cat.# 9304A12)
2. Tetracycline	1.34	25	445.28	154.07	427.32	Dimensions: 100 mm x 2.1 mm ID
3. Thiamphenicol*	1.48	200	354.16	290.04	184.98	Particle Size: 2.7 µm
4. Chlortetracycline	1.56	25	479.27	154.07	371.06	Pore Size: 90 Å
5. Florfenicol*	1.86	200	356.10	336.02	184.98	Guard Column: Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)
6. Chloramphenicol*	1.95	200	321.16	151.99	257.04	Temp.: 40 °C

*Acquired in negative ion mode.

Sample
Diluent: Water
Conc.: 25–200 ng/mL
Inj. Vol.: 2 µL


Mobile Phase
A: 0.1% Formic acid in water
B: 0.1% Formic acid in acetonitrile

Time (min)	Flow (mL/min)	%A	%B
0.00	0.4	90	10
3.00	0.4	20	80
3.01	0.4	90	10
5.00	0.4	90	10

Detector
Ion Mode: MS/MS
Mode: ESI+/ESI-
MRM
UHPLC

Notes
Tetracyclines and amphenicols were analyzed with ESI+ and ESI- mode, respectively.

Figura 3: Chinoloni in LC-MS/MS sulla Raptor C18




Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion
1. Marbofloxacin	1.28	10	363.20	72.11	320.10
2. Norfloxacin	1.32	20	320.23	276.20	233.13
3. Ofloxacin	1.32	10	362.21	318.20	261.15
4. Ciprofloxacin	1.35	20	332.18	288.22	245.15
5. Danofloxacin	1.40	20	358.22	340.16	314.21
6. Enrofloxacin	1.44	10	360.29	316.22	245.13
7. Orbifloxacin	1.47	10	396.22	352.17	226.12
8. Sarafloxacin	1.55	10	386.20	342.20	368.15
9. Difloxacin	1.57	10	400.23	356.17	299.13
10. Oxolinic acid	2.23	5	262.10	244.06	215.96
11. Flumequine	2.78	5	262.15	244.11	202.03

Column	Raptor C18 (cat.# 9304A12)
Dimensions:	100 mm x 2.1 mm ID
Particle Size:	2.7 µm
Pore Size:	90 Å
Guard Column:	Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)
Temp.:	40 °C
Sample	
Diluent:	Water
Conc.:	5–20 ng/mL
Inj. Vol.:	2 µL
Mobile Phase	
A:	0.1% Formic acid in water
B:	0.1% Formic acid in acetonitrile

Time (min)	Flow (mL/min)	%A	%B
0.00	0.4	90	10
3.00	0.4	40	60
3.01	0.4	90	10
5.00	0.4	90	10

Detector	MS/MS
Ion Mode:	ESI+
Mode:	MRM
Instrument	UHPLC

Figura 4: Penicilline, cefalosporine, e tetraciclina in LC-MS/MS sulla Raptor C18




LC_FS0500

Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion
1. Desacetyl cephalirin	0.68	150	382.03	111.92	124.21
2. Amoxicillin	0.89	100	366.24	349.10	208.07
3. Cephalirin	1.05	50	424.17	292.08	124.14
4. Desfuroryl ceftiofur cysteine disulfide	1.06	300	549.16	183.02	126.00
5. Cefquinome	1.16	200	529.19	134.10	125.12
6. Ampicillin	1.18	50	350.19	106.07	160.06
7. Cefalexin	1.19	100	348.10	158.05	174.05
8. Oxytetracycline	1.26	50	461.27	426.15	443.32
9. Cefalonium	1.29	100	459.16	337.03	123.10
10. Tetracycline	1.32	50	445.28	154.07	427.32
11. Cefacetrile	1.37	300	362.07	258.08	178.01
12. Cefazolin	1.50	100	455.10	323.06	295.09
13. Chlortetracycline	1.54	50	479.27	154.07	371.06
14. Cefoperazone	1.69	100	646.26	143.07	148.02
15. Ceftiofur	1.85	50	524.14	241.08	125.24
16. Penicillin G	2.18	100	335.18	176.07	160.07
17. Penicillin V	2.33	100	351.10	160.06	114.07
18. Oxacillin	2.44	100	402.15	160.05	114.06
19. Cloxacillin	2.56	100	436.15	277.06	160.05
20. Nafcillin	2.63	25	415.19	199.09	171.06
21. Dicloxacillin	2.76	100	470.11	160.05	311.02

Column	Raptor C18 (cat.# 9304A12)		
Dimensions:	100 mm x 2.1 mm ID		
Particle Size:	2.7 µm		
Pore Size:	90 Å		
Guard Column:	Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)		
Temp.:	40 °C		
Sample			
Diluent:	Water		
Conc.:	25–300 ng/mL		
Inj. Vol.:	2 µL		
Mobile Phase			
A:	0.1% Formic acid in water		
B:	0.1% Formic acid in acetonitrile		
Time (min)	Flow (mL/min)	%A	%B
0.00	0.4	90	10
3.00	0.4	20	80
3.01	0.4	90	10
5.00	0.4	90	10
Detector	MS/MS		
Ion Mode:	ESI+		
Mode:	MRM		
Instrument	UHPLC		

Figura 5: Sulfamidici in LC-MS/MS sulla Raptor C18




LC_FS0501

Peaks	tr (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion	Column
1. Sulfanilamide	0.83	200	172.98	93.07	75.23	Raptor C18 (cat.# 9304A12)
2. Sulfadiazine	1.45	20	251.18	156.04	92.08	Dimensions: 100 mm x 2.1 mm ID
3. Sulfathiazole	1.60	10	256.16	156.03	92.08	Particle Size: 2.7 µm
4. Sulfapyridine	1.67	10	250.13	156.10	92.08	Pore Size: 90 Å
5. Sulfamerazine	1.79	20	265.08	156.03	92.08	Guard Column: Raptor C18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9304A0252)
6. Sulfamethazine	2.07	10	279.23	186.08	124.08	Temp.: 40 °C
7. Sulfamethizole	2.11	10	271.17	156.02	108.02	Sample Diluent: Water
8. Sulfamethoxypyridazine	2.16	10	281.14	156.03	126.07	Conc.: 10–200 ng/mL
9. Sulfachlorpyridazine	2.55	20	285.05	156.03	108.09	Inj. Vol.: 2 µL
10. Sulfadoxine	2.75	10	311.17	156.03	108.09	Mobile Phase
11. Sulfamethoxazole	2.78	20	254.18	155.98	147.06	A: 0.1% Formic acid in water
12. Sulfathoxypyridazine	2.84	20	295.17	267.07	156.03	B: 0.1% Formic acid in acetonitrile
13. Sulfisoxazole	2.98	20	268.14	156.03	113.10	
14. Sulfadimethoxine	3.37	10	311.17	156.09	108.09	
15. Sulfaquinoxaline	3.40	20	301.18	156.04	108.02	

Time (min)	Flow (mL/min)	%A	%B
0.00	0.4	90	10
3.50	0.4	60	40
3.51	0.4	90	10
5.00	0.4	90	10

Detector	MS/MS
Ion Mode:	ESI+
Mode:	MRM
Instrument	UHPLC

Figura 6: Ionofori in LC-MS/MS sulla Raptor C18



Peaks	tg (min)	Conc. (ng/mL)	Precursor Ion	Product Ion	Product Ion
1. Lasalocid A	1.92	100	613.42	377.28	595.40
2. Monensin	2.12	100	693.50	675.44	461.30
3. Salinomycin	2.19	100	773.57	431.24	531.39
4. Maduramicin	2.30	100	939.65	877.58	473.31
5. Narasin	2.58	100	787.59	431.27	531.35

Column Raptor Biphenyl (cat.# 9309A12)
Dimensions: 100 mm x 2.1 mm ID
Particle Size: 2.7 µm
Pore Size: 90 Å
Guard Column: Raptor Biphenyl EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9309A0252)
Temp.: 40 °C
Sample
Diluent: Water:methanol (10:90)
Conc.: 100 ng/mL
Inj. Vol.: 2 µL
Mobile Phase
A: 0.5% Formic acid in water
B: 0.5% Formic acid in acetonitrile

Time (min)	Flow (mL/min)	%A	%B
0.00	0.5	35	65
3.00	0.5	25	75
3.01	0.5	35	65
5.00	0.5	35	65

Detector MS/MS
Ion Mode: ESI+
Mode: MRM
Instrument UHPLC