

# Streamline Your MCPD and Glycidol Analysis with Our New Certified Reference Materials



# Streamline Your MCPD and Glycidol Analysis with Our New Certified Reference Materials

Our new MCPD esters and glycidyl stearate certified reference materials (CRMs) help labs streamline their analysis of these food processing contaminants. By relying on our manufacturing expertise, labs can minimize errors, save time, and reduce cost compared to sourcing neat materials.

A high concentration of up to 100 µg/mL adds additional flexibility in creating working standards, and an optimized GC-MS method helps labs simplify method development.

- Compounds meet DIN EN ISO 18363-1; AOCS Cd29c-18; and DGF C-VI 18 (10) method requirements.
- Lower concentration, higher volume ampul option of PP-3-MCPD-d5 (10 µg/mL, 10 mL/ampul) for automated systems analyzing 3-MCPD and glycidol.
- Certified Reference Material (CRM) manufactured and QC-tested in our ISO-accredited labs.
- Formulated for optimal stability and usability to ensure your calibration standards are accurate.



## 3-MCPD Standard

3-Chloro-1,2-propanediol (96-24-2)/2-Methylfuran (534-22-5)

Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
100 µg/mL, P&T Methanol, 1 mL/ampul	Yes	6 months	36 months	Ambient	10 °C or colder	ea.	31016

## PP-3-MCPD-d5 Standard

- Lower concentration, higher volume ampul option available for automated systems (cat.# 31017).

3-Chloro-1,2-propanediol Dipalmitate-d5 (1185057-55-9)

Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
10 µg/mL, Toluene, 10 mL/ampul	Yes	6 months	24 months	Ambient	10 °C or colder	ea.	31017
100 µg/mL, Toluene, 1 mL/ampul	Yes	6 months	24 months	Ambient	10 °C or colder	ea.	31018

## PP-3-MCPD Standard

3-Chloro-1,2-propanediol Dipalmitate (51930-97-3)


Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
100 µg/mL, Toluene, 1 mL/ampul	Yes	6 months	24 months	Ambient	10 °C or colder	ea.	31019

## Glycidyl Stearate Standard

Glycidyl Stearate (7460-84-6)

Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
100 µg/mL, Toluene, 1 mL/ampul	Yes	6 months	24 months	Ambient	10 °C or colder	ea.	31020

**Figure 1:** Restek's Rxi-17Sil MS GC columns provide baseline separation of 3-MCPD and 3-MCPD-d5 on a 30 m column.



Peaks	tr (min)	Conc. (µg/mL)
1. 3-MCPD-d5	5.954	4
2. 3-MCPD	5.994	4

<b>Column</b>	Rxi-17Sil MS GC, 30 m, 0.32 mm ID, 0.25 µm (cat.# 14124)
<b>Standard/Sample</b>	3-MCPD (cat.# 31016) PP-MCPD-d5 (cat.# 31018) PP-MCPD (cat.# 31019) Glycidyl stearate (cat.# 31020)
<b>Diluent:</b>	isooctane
<b>Conc.:</b>	4 µg/mL
<b>Injection</b>	
Inj. Vol.:	2 µL split (split ratio 10:1)
Liner:	Topaz, straight split inlet liner, 3.5 mm x 5.0 x 95, for Shimadzu 17A, 2010, 2014, and 2030 GCs, w/ quartz wool (cat.# 23319)
Inj. Temp.:	280 °C
Split Vent Flow Rate:	18.4 mL/min
<b>Oven</b>	
Oven Temp.:	100 °C (hold 0.5 min) to 180 °C at 12 °C/min (hold 0 min) to 320 °C at 25 °C/min (hold 4 min)
<b>Carrier Gas</b>	He, constant flow
Flow Rate:	1.4 mL/min
Linear Velocity:	44 cm/sec @ 280 °C
<b>Detector</b>	Shimadzu GCMS-QP 2020 NX
SIM Program:	147, 196, 150, 201 m/z, 0.3 ms dwell
Transfer Line Temp.:	280 °C
Analyzer Type:	Quadrupole
Source Temp.:	200 °C
Solvent Delay Time:	3 min
Tune Type:	PFTBA
Ionization Mode:	EI
<b>Instrument</b>	Shimadzu GC-2030
<b>Sample Preparation</b>	Two samples were prepared in parallel. 20 µL of PP-MCPD-d5 (100 µg/mL cat.# 31018), 20 µL of PP-MCPD (100 µg/mL cat.# 31019), and 100 mg of virgin olive oil were added to each sample. 20 µL of glycidyl stearate (100 µg/mL cat.# 31020) was added to the second sample. One-hundred microliters of MTBE was added to both samples, and the solutions were vortexed for 30 s. Two-hundred microliters of a NaOH solution (0.5 M in methanol) was added to each. Both samples were vortexed for 30 s and then allowed to sit for 5 min. Six-hundred microliters of acidified NaCl (3.4 M NaCl, 0.875% sulfuric acid) was added to quench the reaction. The solution was washed twice with 600 µL of n-hexane with the organic layer being discarded. The aqueous solutions were extracted into 600 µL of 6:4 ether:ethyl acetate solution three times. The combined organic layers were stirred over sodium sulfate for 5 mins. Twenty microliters of a saturated phenyl boronic acid was added to the each and then was evaporated to dryness under a stream of nitrogen. The samples were re-extracted into 500 µL of isooctane for final analysis.
<b>Notes</b>	Unreacted phenyl boronic acid appears at 11.32 min. Lower concentration, higher volume ampul option of PP-3-MCPD-d5 (10 µg/mL, 10 mL/ampul) available for automated systems (cat.# 31017).

## Rxi-17Sil MS GC Columns

df (Film Thickness)	Inner Diameter (ID)	Length	cat.#
0.25	0.32 mm	30 m	14124



## Topaz Straight Split Inlet Liner

for Shimadzu 17A, 2010, 2014, and 2030 GCs

Internal Diameter (ID)	Outer Diameter (OD)	Length	Packing Materials	Units	Cat.#
3.5 mm	5.0 mm	95 mm	Quartz Wool	5-pk.	23319



## Custom Reference Standards

Do you need specific compounds to meet your method requirements?

Visit [www.restek.com/standards](http://www.restek.com/standards) and start a request today!



Additional product configurations available at [www.restek.com](http://www.restek.com)

# RESTEK

For information on Restek patents and trademarks, visit [www.restek.com/patents-trademarks](http://www.restek.com/patents-trademarks). To unsubscribe from future Restek communications or to update your preferences, visit [www.restek.com/subscribe](http://www.restek.com/subscribe). To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.  
© 2025 Restek Corporation. All rights reserved.

[www.restek.com](http://www.restek.com)



Lit. Cat.# FSSS4250B-UNV