GC Columns Installation Guide

The Restek Applications and Technical Services specialists have found this to be a reliable sequence for avoiding problems when installing a fused silica or MXT capillary GC column.

Instrument Preparation & Column Installation

1. Cool all heated zones.
2. Visually inspect indicating oxygen and moisture traps. Replace saturated traps.
3. Examine the inlet and the detector. Clean or replace all dirty or corroded parts.
4. Replace the inlet liner, septum, and the injector seals (O-rings, inlet seals, ferrules, etc.).
5. Uncoil the column ends to make sure they are long enough to reach the injector and detector. Cut a few centimeters from each end of the column.

Note: To correctly cut a capillary column, use a smooth-edged ceramic scoring wafer (cat.# 20116) and follow the directions provided in its instruction sheet (200-00-001).

6. For MXT columns, use your hands to gently straighten the ends of the column to remove the curve in the tubing.
7. Mount the column in the oven onto a support. For fused silica columns, ensure the support protects it from scratches. Center the column in the oven. This ensures uniform heat exposure generating consistent retention times.

Note: Restek has several types of cages for fused silica columns. If you have a cage or column with high-temperature string, do not remove it from the column.

8. Slide the nut and appropriate size ferrule onto the inlet end of the column. Cut an additional 2 cm from the end of the column to remove any residual ferrule material from the end of the column.
9. Install the column the appropriate distance into the injection port as indicated in your instrument manual.
10. Set the carrier gas to the flow rate or inlet pressure recommended for the column or to your method flow rate/pressure. Confirm presence of column flow by immersing the column outlet in a vial of solvent and looking for bubbles.
11. Install the column into the detector as described in Steps 8–9. Set the detector gases to proper settings.
12. Set the injector and detector temperatures. Do not exceed the column's maximum isothermal operating temperature (listed on the column box). Check column connections for leaks using a Restek Electronic Leak Detector (cat.# 28500).
13. Purge the column at ambient temperature with carrier gas for at least 15 minutes.
14. Load your method to set the appropriate flow rates and temperatures.
15. Inject an unretained compound to verify the column is installed correctly and to determine the dead time. A symmetric peak indicates the column is installed correctly. Adjust the carrier gas flow as necessary to obtain proper dead time.
16. Condition the column 20 °C above the final analysis temperature of your method, but do not exceed the column's maximum isothermal operating temperature. For most applications, one hour of conditioning is sufficient. For sensitive detectors or low-level analysis, longer conditioning times or conditioning the column at the maximum isothermal operating temperature may be beneficial. Extended time at high temperatures will not damage the column if precautions are taken to make sure the carrier gas is clean and is filtered for oxygen and water. Additionally, it is important to ensure the connections are leak free.
17. To confirm proper column and system performance, analyze a column test mix or a known standard.
18. Your GC system is now ready to be calibrated and acquire data.

Note: When you intend to condition thick-film columns (film thickness > 1 µm) at temperatures near the maximum isothermal operating temperature, it is recommended to do the initial one to two hours conditioning without a connection to the detector, making sure the column end exits the GC oven.

Standby Conditions

Short-Term: Leave the column in the GC with the carrier gas flow on at an oven temperature of 100-150 °C.

Long-Term: Remove the column from the GC and seal the ends by carefully pushing each end into the curved edge of a septum. For fused silica columns, store the column in the original box away from strong lighting.

If you have any questions or problems installing a Restek column, visit www.restek.com/gcinstall or contact your Restek representative.

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Questions about this or any other Restek product?
Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com

GNOT1098E-UNV

Rev. date: 06/24