

Certified Reference Materials

Accurately Quantify 380 Pesticides

Using LC-MS/MS and GC-MS/MS

- CRM kits include comprehensive lists of pesticides of global concern.
- Optimized chromatography methods available using Raptor LC and Rxi GC columns by mass spec.
- Complete your full workflow solution with additional Restek standards, sample preparation products, and applications.



RESTEK

Pure Chromatography

www.restek.com



Accurately Quantify 380 Pesticides Using LC-MS/MS and GC-MS/MS with Restek's Multiresidue Pesticide CRM Kits

- Accurately detect and quantify pesticides of global environmental and food safety concern in a wide range of matrices.
- Formulated and grouped for maximum long-term stability.*
- Quantitatively tested to confirm composition; detailed support documentation provided.
- Optimized multiresidue pesticide method is offered free of charge; downloadable XLS file includes conditions and transition tables.
- Certified reference material (CRM) manufactured and QC-tested in Restek's ISO-accredited labs satisfies your ISO requirements.

* NOTE: When combining a large number of compounds with different chemical functionalities, mix stability can be an issue. In formulating these standards, we extensively studied all the compounds involved, then grouped them into as few mixes as possible while still ensuring maximum long-term stability and reliability. For quantitative analysis, we recommend analyzing each mix separately to ensure accurate results for every compound.

For full compound lists, chromatograms, and analysis conditions, visit
www.restek.com/multiresidue



Custom Solutions

Don't see exactly what you need? Restek can blend custom formulations to match your specific requirements

Build your mix today at
www.restek.com/solutions

LC Multiresidue Pesticide Kit

Full kit contains 204 compounds of interest, covering many LC-determined pesticides listed by government agencies; individual ampuls also sold separately. For chromatograms and optimized multiresidue pesticide method including conditions and transition tables, visit www.restek.com/multiresidue-LC

Cat. # 31972: LC Multiresidue Pesticide Standard #1 (13 components)

Organophosphorus Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Acephate (30560-19-1)
Carbaryl (Sevin) (63-25-2)
Dicrotophos (141-66-2)
Dimethoate (60-51-5)
Dimethomorph (110488-70-5)
Isocarbophos (24353-61-5)
Methamidophos (10265-92-6)
Mevinphos (7786-34-7)
Monocrotophos (6923-22-4)
Omethoate (1113-02-6)
Temephos (Abate) (3383-96-8)
Trichlorfon (Dylox) (52-68-6)
Vamidothion (Vamidoate) (2275-23-2)

Cat. # 31973: LC Multiresidue Pesticide Standard #2 (16 components)

Carbamate/Uron Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Alanycarb (83130-01-2)
Aldicarb (116-06-3)
Aldicarb sulfone (1646-88-4)
Aldicarb sulfoxide (1646-87-3)
Benfuracarb (82560-54-1)
Butocarboxim (34681-10-2)
Butoxycarboxim (34681-23-7)
Ethiofencarb (29973-13-5)
Furathiocarb (65907-30-4)
Methabenzthiazuron (18691-97-9)
Methiocarb (2032-65-7)
Methomyl (16752-77-5)
Oxamyl (23135-22-0)
Tebuthiuron (34014-18-1)
Thidiazuron (51707-55-2)
Thiophanate-methyl (23564-05-8)

Cat. # 31974: LC Multiresidue Pesticide Standard #3 (38 components)

Carbamate/Uron Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Bendiocarb (22781-23-3)
Bifenazate (149877-41-8)
Carbofuran (1563-66-2)
Chlorfluazuron (71422-67-8)
Chloroxuron (1982-47-4)
Chlortoluron (15545-48-9)
Cycluron (2163-69-1)
Diethofencarb (87130-20-9)
Diflubenuron (35367-38-5)
Dioxacarb (6988-21-2)
Diuron (330-54-1)
Fenobucarb (BPMC) (3766-81-2)
Fenoxycarb (72490-01-8)
Fenuron (101-42-8)
Flufenoxuron (101463-69-8)
Fluometuron (2164-17-2)
Forchlorfenuron (68157-60-8)
Hexaflumuron (86479-06-3)
3-Hydroxycarbofuran (16655-82-6)
Indoxacarb (173584-44-6)
Iprovalicarb (140923-17-7)

Isoprocab (2631-40-5)
Isoproturon (34123-59-6)
Linuron (330-55-2)
Lufenuron (103055-07-8)
Metobromuron (3060-89-7)
Monolinuron (1746-81-2)
Neburon (555-37-3)
Novaluron (116714-46-6)
Pirimicarb (23103-98-2)
Promecarb (2631-37-0)
Propham (122-42-9)
Propoxur (Baygon) (114-26-1)
Pyraclostrobin (175013-18-0)
Siduron (1982-49-6)
Teflubenzuron (83121-18-0)
Thiobencarb (28249-77-6)
Triflururon (64628-44-0)

Cat. # 31975: LC Multiresidue Pesticide Standard #4 (63 components)

Organonitrogen Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Abamectin (71751-41-2)
Acetamiprid (135410-20-7)
Ametryn (834-12-8)
Amiratz (33089-61-1)
Azoxytrobin (131860-33-8)
Benalaxyl (71626-11-4)
Benzoximate (29104-30-1)
Boscalid (188425-85-6)
Butafenacil (134605-64-4)
Carbetamide (16118-49-3)
Carfentrazone ethyl (128639-02-1)
Chlorantraniliprole (500008-45-7)
Clofentezine (74115-24-5)
Cymoxanil (57966-95-7)
Cyprodinil (121552-61-2)
Cyromazine (66215-27-8)
Dimoxystrobin (149961-52-4)
Dinotefuran (165252-70-0)
Doramectin (117704-25-3)
Eprinomectin (123997-26-2)
Famoxadon (131807-57-3)
Fenazaquin (120928-09-8)
Fenhexamid (126833-17-8)
(E)-Fenpyroximate (134098-61-6)
Fonicamid (158062-67-0)
Fluzinam** (79622-59-6)
Fludioxonil (131341-86-1)
Fluoxastrobin (361377-29-9)
Flutolanil (66332-96-5)
Furalaxyl (57646-30-7)
Halofenozide (112226-61-6)
Imazalil (35554-44-0)
Imidacloprid (138261-41-3)
Ivermectin (70288-86-7)
Kresoxim methyl (143390-89-0)
Mandipropamid (374726-62-2)
Mepanipyrim (110235-47-7)
Meprothion (55814-41-0)
Metaflumizone (139968-49-3)
Metalaxyl (57837-19-1)
Methoxyfenozide (161050-58-4)
Moxidectin (113507-06-5)

Myclobutanil (88671-89-0)
Nitenpyram (120738-89-8)
Oxadixyl (77732-09-3)
Picoxystrobin (117428-22-5)
Piperonyl butoxide (51-03-6)
Prochloraz (67747-09-5)
Prometon (1610-18-0)
Pymetrozine (123312-89-0)
Pyracarbolid (24691-76-7)
Pyrimethanil (53112-28-0)
Pyriproxyfen (95737-68-1)
Quinoxifen (124495-18-7)
Rotenone (83-79-4)
Sebumeeton (26259-45-0)
Spiroxamine (118134-30-8)
Tebufenozide (112410-23-8)
Tebufenpyrad (119168-77-3)
Terbumeton (33693-04-8)
Triadimefon (43121-43-3)
Trifloxystrobin (141517-21-7)
Zoxamide (156052-68-5)

Cat. # 31976: LC Multiresidue Pesticide Standard #5 (30 components)

Organonitrogen Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Acibenzolar-S-methyl (135158-54-2)
Bupirimate (41483-43-6)
Buprofezin (69327-76-0)
Carboxin (5234-68-4)
Clethodim (99129-21-2)
Clothianidin (210880-92-5)
Cyazofamid (120116-88-3)
Ethiprole (181587-01-9)
Ethofumesate (26225-79-6)
Fenamidone (161326-34-7)
Fipronil (120068-37-3)
Flubendiamide (272451-65-7)
Flufenacet (Fluthiamide) (142459-58-3)
Hexythiazox (78587-05-0)
Mefenacet (73250-68-7)
Mesotrione (104206-82-8)
Methoprotrene (841-06-5)
Metribuzin (21087-64-9)
Prometryne (7287-19-6)
Propargite (2312-35-8)
Prothioconazole (178928-70-6)
Pyridaben (96489-71-3)
Simetryn (1014-70-6)
Sulfentrazon (122836-35-5)
Terbutryn (886-50-0)
Thiabendazole (148-79-8)
Thiacloprid (111988-49-9)
Thiamethoxam (153719-23-4)
Thiofanox (39196-18-4)
Tricyclazole (Beam) (41814-78-2)

Cat. # 31977: LC Multiresidue Pesticide Standard #6 (28 components)

Organonitrogen Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Baycor (Bitteranol) (55179-31-2)
Bromuconazole (116255-48-2)
Cyproconazole (94361-06-5)

Diclobutrazol (75736-33-3)
Difenoconazole (119446-68-3)
Diniconazole (83657-24-3)
Epoxiconazole (133855-98-8)
Etaconazole (60207-93-4)
Ethinilol (23947-60-6)
Etozazole (153233-91-1)
Fenarimol (60168-88-9)
Fenbuconazole (114369-43-6)
Fluquinconazole (136426-54-5)
Flusilazole (85509-19-9)
Flutriafol (76674-21-0)
Fuberidazole (3878-19-1)
Hexaconazole (79983-71-4)
Ipconazole (125225-28-7)
Metconazole (125116-23-6)
Nuarimol (63284-71-9)
Paclobutrazol (76738-62-0)
Penconazole (66246-88-6)
Propiconazole (Tilt) (60207-90-1)
Tebuconazole (107534-96-3)
Tetraconazole (112281-77-3)
Triadimenol (55219-65-3)
Triflumizole (68694-11-1)
Triticonazole (131983-72-7)

Cat. # 31978: LC Multiresidue Pesticide Standard #7 (7 components)

Organonitrogen Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Emamectin-benzoate (155569-91-8)
Fenpropimorph (67564-91-4)
Spinetoram (J&L) (935545-74-7)†
Spinosad (A&D) (168316-95-8)††
Spirodiclofen (148477-71-8)
Spiromesifen (283594-90-1)
Spirotetramat (203313-25-1)

Cat. # 31979: LC Multiresidue Pesticide Standard #8

Organonitrogen Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Hydramethylnon (67485-29-4)

Cat. # 31980: LC Multiresidue Pesticide Standard #9 (7 components)

Carbamate/Uron Compounds

100 µg/mL each in acetonitrile, 1 mL/ampul
Aminocarb (2032-59-9)
Desmedipham (13684-56-5)
Formetanate HCL (23422-53-9)
Mexacarbate (Zectran) (315-18-4)
Monceren (Pencycuron) (66063-05-6)
Phenmedipham (13684-63-4)
Propamocarb free base (24579-73-5)

Cat. # 31981: LC Multiresidue Pesticide Standard #10

Carbamate/Uron Compounds

100 µg/mL each in methanol, 1 mL/ampul
Carbendazim (10605-21-7)

Description	Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
LC Multiresidue Pesticide Kit	Contains 1 mL each of these mixtures.	Yes	6 months	On Ice	-20 °C or colder	kit	31971

* Note: When combining a large number of compounds with different chemical functionalities, mix stability can be an issue. In formulating these standards, we extensively studied the 204 compounds involved, then grouped them into a few mixes as possible while still ensuring maximum long-term stability and reliability. For quantitative analysis, we recommend analyzing each mix separately to ensure accurate results for every compound.

** Note: In this standard, fluzinam should only be used for qualitative analysis. A single-component standard (cat.# 31982) is available for quantitative analysis.

†This reference material contains both the Spinetoram J and L isomers (187166-40-1 and 187166-15-0); however, CAS# 935545-74-7 is displayed on the certificate as this is the neat material dissolved in the solution. CAS# 935545-74-7 is a blend of Spinetoram J and L and the ratio of each material isomer are displayed on your certificate of analysis.

††This reference material contains Spinosad A and D isomers (131929-60-7 and 131929-63-0); however, CAS# 168316-95-8 is displayed on the certificate as this is the neat material dissolved in the solution. CAS# 168316-95-8 is a blend of Spinosad A and D and the ratio of each material isomer are displayed your certificate of analysis.

GC Multiresidue Pesticide Kit

Comprehensive 203-compound kit covers food safety lists by the FDA, USDA, and other global governmental agencies; individual ampuls also sold separately. For chromatograms and optimized multiresidue pesticide method including conditions and transition tables, visit www.restek.com/multiresidue-GC

Cat. # 32563: GC Multiresidue Pesticide Standard #1 (16 components)

Organophosphorus Compounds

100 µg/mL each in toluene, 1 mL/ampul
Azinphos ethyl (2642-71-9)
Azinphos methyl (86-50-0)
Chlorpyrifos (2921-88-2)
Chlorpyrifos methyl (5598-13-0)
Diazinon (333-41-5)
EPN (2104-64-5)
Fenitrothion (122-14-5)
Isazophos (42509-80-8)
Phosalone (2310-17-0)
Phosmet (732-11-6)
Pirimiphos ethyl (23505-41-1)
Pirimiphos methyl (29232-93-7)
Pyraclofos (89784-60-1)
Pyraclofos (13457-18-6)
Pyridaphenthion (119-12-0)
Quinalphos (13593-03-8)

Cat. # 32564: GC Multiresidue Pesticide Standard #2 (40 components)

Organochlorine Compounds

100 µg/mL each in toluene, 1 mL/ampul
Aldrin (309-00-2)
α-BHC (319-84-6)
β-BHC (319-85-7)
δ-BHC (319-86-8)
γ-BHC (Lindane) (58-89-9)
Chlorbenside (103-17-3)
cis-Chlordane (5103-71-9)
trans-Chlordane (5103-74-2)
Chlorfenson (Ovex) (80-33-1)
Chloroneb (2675-77-6)
2,4'-DDD (53-19-0)
4,4'-DDD (72-54-8)
2,4'-DDE (3424-82-6)
4,4'-DDE (72-55-9)
2,4'-DDT (789-02-6)
4,4'-DDT (50-29-3)
4,4'-Dichlorobenzophenone (90-98-2)
Dieldrin (60-57-1)
Endosulfan I (959-98-8)
Endosulfan II (33213-65-9)
Endosulfan ether (3369-52-6)
Endosulfan sulfate (1031-07-8)
Endrin (72-20-8)
Endrin aldehyde (7421-93-4)
Endrin ketone (53494-70-5)
Ethylan (Perthane) (72-56-0)
Fenson (80-38-6)
Heptachlor (76-44-8)
Heptachlor epoxide (isomer B) (1024-57-3)
Hexachlorobenzene (118-74-1)
Isodrin (465-73-6)
2,4'-Methoxychlor (30667-99-3)
4,4'-Methoxychlor olefin (2132-70-9)
Mirex (2385-85-5)
cis-Nonachlor (5103-73-1)
trans-Nonachlor (39765-80-5)
Pentachloroanisole (1825-21-4)
Pentachlorobenzene (608-93-5)
Pentachlorothioanisole (1825-19-0)
Tetradifon (116-29-0)

Cat. # 32565: GC Multiresidue Pesticide Standard #3 (25 components)

Organonitrogen Compounds

100 µg/mL each in toluene:acetonitrile (99:1), 1 mL/ampul
Benfluralin (1861-40-1)
Biphenyl (92-52-4)
Chlorothalonil (1897-45-6)
Dichlofluanid (1085-98-9)
Dichloran (99-30-9)
3,4-Dichloroaniline (95-76-1)
2,6-Dichlorobenzonitrile (Dichlobenil) (1194-65-6)
Diphenylamine (122-39-4)
Ethalfuralin (55283-68-6)
Fluchloralin (33245-39-5)
Isopropalin (33820-53-0)
Nitrilin (4726-14-1)
Nitrofen (1836-75-5)
Oxyfluorfen (42874-03-3)
Pendimethalin (40487-42-1)
Pentachloroaniline (527-20-8)
Pentachlorobenzonitrile (20925-85-3)
Pentachloronitrobenzene (Quintozone) (82-68-8)
Prodiamine (29091-21-2)
Profluralin (26399-36-0)
2,3,5,6-Tetrachloroaniline (3481-20-7)
Tetrachloronitrobenzene (Tecnazene) (117-18-0)
THPI (Tetrahydrophthalimide) (1469-48-3)
Tolylfluaniid (731-27-1)
Trifluralin (1582-09-8)

Cat. # 32566: GC Multiresidue Pesticide Standard #4 (28 components)

Organonitrogen Compounds

100 µg/mL each in toluene, 1 mL/ampul
Acetochlor (34256-82-1)
Alachlor (15972-60-8)
Allidochlor (93-71-0)
Clomazone (Command) (81777-89-1)
Cycloate (1134-23-2)
Diallate (cis & trans) (2303-16-4)
Dimethachlor (50563-36-5)
Diphenamid (957-51-7)
Fenpropathrin (39515-41-8)
Fluquinconazole (136426-54-5)
Flutolanil (66332-96-5)
Linuron (330-55-2)
Metazachlor (67129-08-2)
Methoxychlor (72-43-5)
Metolachlor (51218-45-2)
N-(2,4-Dimethylphenyl)formamide (60397-77-5)
Norflurazon (27314-13-2)
Oxadiazon (19666-30-9)
Pebulate (1114-71-2)
Pretilachlor (51218-49-6)
Prochloraz (67747-09-5)
Propachlor (1918-16-7)
Propanil (709-98-8)
Propisochlor (86763-47-5)
Propyzamide (23950-58-5)
Pyridaben (96489-71-3)
Tebufenpyrad (119168-77-3)
Triallate (2303-17-5)

Cat. # 32567: GC Multiresidue Pesticide Standard #5 (34 components)

Organonitrogen Compounds

100 µg/mL each in toluene, 1 mL/ampul
Atrazine (1912-24-9)
Bupirimate (41483-43-6)
Captafol (2425-06-1)
Captan (133-06-2)
Chlorfenapyr (122453-73-0)
Cyprodinil (121552-61-2)
Etofenprox (80844-07-1)
Etridiazole (2593-15-9)
Fenarimol (60168-88-9)
Fipronil (120068-37-3)
Fludioxonil (131341-86-1)
Fluridone (Sonar) (59756-60-4)
Flusilazole (85509-19-9)
Flutriafol (76674-21-0)
Folpet (133-07-3)
Hexazinone (Velpar) (51235-04-2)
Iprodione (36734-19-7)
Lenacil (2164-08-1)
MGK-264 (113-48-4)
Myclobutanil (88671-89-0)
Paclobutrazol (76738-62-0)
Penconazole (66246-88-6)
Procymidone (32809-16-8)
Propargite (2312-35-8)
Pyrimethanil (53112-28-0)
Pyriproxyfen (95737-68-1)
Tebuconazole (107534-96-3)
Terbacil (5902-51-2)
Terbutylazine (5915-41-3)
Triadimefon (43121-43-3)
Triadimenol (55219-65-3)
Tricyclazole (Beam) (41814-78-2)
Triflumizole (68694-11-1)
Vinclozolin (50471-44-8)

Cat. # 32568: GC Multiresidue Pesticide Standard #6 (18 components)

Synthetic Pyrethroid Compounds

100 µg/mL each in toluene, 1 mL/ampul
Acrinathrin (101007-06-1)
Anthraquinone (84-65-1)
Bifenthrin (82657-04-3)
Bioallethrin (584-79-2)
Cyfluthrin (68359-37-5)
lambda-Cyhalothrin (91465-08-6)
Cypermethrin (52315-07-8)
Deltamethrin (52918-63-5)
Fenvalerate (51630-58-1)
Flucythrinate (70124-77-5)
tau-Fluvalinate (102851-06-9)
cis-Permethrin (61949-76-6)
trans-Permethrin (61949-77-7)
Phenothrin (cis & trans) (26002-80-2)
Resmethrin (10453-86-8)
Tefluthrin (79538-32-2)
Tetramethrin (7696-12-0)
Transfluthrin (118712-89-3)

Cat. # 32569: GC Multiresidue Pesticide Standard #7 (10 components)

Herbicide Methyl Esters

100 µg/mL each in toluene, 1 mL/ampul
Acephinocyl (57960-19-7)
Bromopropylate (18181-80-1)
Carfentrazone ethyl (128639-02-1)
Chlorobenzilate (510-15-6)
Chlorpropham (101-21-3)
Chlozolate (84332-86-5)
DCPA methyl ester (Chlorthal-dimethyl) (1861-32-1)
Fluazifop-*p*-butyl (79241-46-6)
Metalaxyl (57837-19-1)
2-Phenylphenol (90-43-7)

Cat. # 32570: GC Multiresidue Pesticide Standard #8 (24 components)

Organophosphorus Compounds

100 µg/mL each in toluene, 1 mL/ampul
Bromfenvinfos-methyl (13104-21-7)
Bromfenvinfos (33399-00-7)
Bromophos ethyl (4824-78-6)
Bromophos methyl (2104-96-3)
Carbophenothion (786-19-6)
Chlorfenvinfos (470-90-6)
Chlorthiophos (60238-56-4)
Coumaphos (56-72-4)
Edifenphos (17109-49-8)
Ethion (563-12-2)
Fenamiphos (22224-92-6)
Fenchlorphos (Ronnell) (299-84-3)
Fenthion (55-38-9)
Iodofenphos (18181-70-9)
Leptophos (21609-90-5)
Malathion (121-75-5)
Methacrifos (62610-77-9)
Profenofos (41198-08-7)
Prothiofos (34643-46-4)
Sulfotepp (3689-24-5)
Sulprofos (35400-43-2)
Terbufos (13071-79-9)
Tetrachlorvinphos (22248-79-9)
Tolclofos-methyl (57018-04-9)

Cat. # 32571: GC Multiresidue Pesticide Standard #9 (8 components)

Organophosphorus Compounds

100 µg/mL each in toluene, 1 mL/ampul
Disulfoton (298-04-4)
Fonofos (944-22-9)
Methyl parathion (298-00-0)
Mevinphos (7786-34-7)
Parathion (ethyl parathion) (56-38-2)
Phorate (298-02-2)
Piperonyl butoxide (51-03-6)
Triazophos (24017-47-8)

Description	Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat.#
GC Multiresidue Pesticide Kit	Contains 1 mL each of these mixtures.	Yes	6 months	Ambient	10 °C or colder	kit	32562

Restek is Your Complete Supplier for World-Class LC-MS/MS and GC-MS/MS Multiresidue Pesticide Analysis.

We carry certified reference materials, analytical columns, sample preparation products, and accessories.

Labeled Pesticide Residue Internal Standards for Food Analysis

Isotopically labeled to provide the best approach for pesticide residue quantification.

Description	CAS #	Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat. #
Atrazine-d5	163165-75-1	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	36 months	Ambient	10 °C or colder	ea.	31984
Carbaryl-d7	362049-56-7	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	31 months	Ambient	10 °C or colder	ea.	31985
Diazinon-d10 (diethyl-d10)	100155-47-3	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	36 months	Ambient	10 °C or colder	ea.	31986
Dichlorvos-d6	203645-53-8	100 µg/mL in acetone, 1 mL/ampul	Yes	3 months	12 months	Ambient	10 °C or colder	ea.	31987
Dimethoate-d6	1219794-81-6	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	36 months	Ambient	10 °C or colder	ea.	31988
Diuron-d6	1007536-67-5	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	31 months	Ambient	10 °C or colder	ea.	31989
Linuron-d6	1219804-76-8	100 µg/mL in acetonitrile, 1 mL/ampul	Yes	6 months	31 months	Ambient	10 °C or colder	ea.	31990

QuEChERS Internal Standard Mix for GC-MS Analysis

(6 components)

PCB 18 (37680-65-2), 50 µg/mL
PCB 28 (7012-37-5), 50 µg/mL
PCB 52 (35693-99-3), 50 µg/mL

Triphenylmethane (519-73-3), 10 µg/mL
Triphenylphosphate (115-86-6), 20 µg/mL
Tris(1,3-dichloroisopropyl)phosphate (13674-87-8), 50 µg/mL

Conc. in Solvent	CRM?	Min Shelf Life on Ship Date	Max Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat. #
In acetonitrile, 5 mL/ampul	Yes	6 months	75 months	Ambient	10 °C or colder	ea.	33267



QuEChERS Performance Standards Kit

Contains 1 mL each of these mixtures. 31153: QuEChERS Performance Standard A 31154: QuEChERS Performance Standard B 31155: QuEChERS Performance Standard C 300 µg/mL each in acetonitrile/acetic acid (99.9:0.1), 1 mL/ampul. Blend equal volumes of all three ampuls for a 100 µg/mL final solution.

Cat. # 31153: QuEChERS Performance Standard A (16 components)
Acephate (30560-19-1)
Azinphos methyl (86-50-0)
Chlorpyrifos (2921-88-2)
Coumaphos (56-72-4)
Diazinon (333-41-5)
Dichlofluanid (1085-98-9)
Dichlorvos (DDVP) (62-73-7)
Dimethoate (60-51-5)
Fenthion (55-38-9)

Malathion (121-75-5)
Methamidophos (10265-92-6)
Mevinphos (7786-34-7)
Omethoate (1113-02-6)
Phosalone (2310-17-0)
Pirimiphos methyl (29232-93-7)
Propargite (2312-35-8)

Cat. # 31154: QuEChERS Performance Standard B (7 components)
γ-BHC (Lindane) (58-89-9)

Chlorothalonil (1897-45-6)
4,4'-DDT (50-29-3)
Dicofol (Kelthane) (115-32-2)
Endosulfan sulfate (1031-07-8)
Endrin (72-20-8)
2-Phenylphenol (90-43-7)

Cat. # 31155: QuEChERS Performance Standard C (17 components)
Bifenthrin (82657-04-3)
Captan (133-06-2)

Carbaryl (Sevin) (63-25-2)
Cyprodinil (121552-61-2)
Deltamethrin (52918-63-5)
Fenhexamid (126833-17-8)
Fenpropathrin (39515-41-8)
Folpet (133-07-3)
Imazalil (35554-44-0)
Iprodione (36734-19-7)
Metalaxyl (57837-19-1)
Methiocarb (2032-65-7)
Myclobutanil (88671-89-0)
cis-Permethrin (61949-76-6)

trans-Permethrin (61949-77-7)
Thiabendazole (148-79-8)
Vinclozolin (50471-44-8)

Description	CRM?	Min Shelf Life on Ship Date	Shipping Conditions	Storage Temp.	qty.	cat. #
QuEChERS Performance Standards Kit	Yes	3 months	Ambient	10 °C or colder	kit	31152



Raptor ARC-18 LC Columns (USP L1)

- Ideal for high-throughput LC-MS/MS applications with minimal sample preparation.
- Well-balanced retention profile for better detection and integration of large, multiclass analyte lists.
- Sterically protected to endure low-pH mobile phases without sacrificing retention or peak quality.
- Part of Restek's Raptor LC column line featuring 1.8, 2.7, and 5 μm SPP core-shell silica.

ID	Length	qty.	cat.#
2.7 μm Particles Raptor ARC-18			
2.1 mm	100 mm	ea.	9314A12

For more column dimensions and Raptor phases, visit www.restek.com/raptor



Rxi-5ms Columns (fused silica)

low-polarity phase; Crossbond diphenyl dimethyl polysiloxane

- Ideal for pesticides in food.
- General-purpose columns that can be used for phenols, residual solvents, drugs of abuse, pesticides, semivolatiles, PCB congeners (e.g., Aroclor mixes), and solvent impurities.
- Tested and guaranteed for ultra-low bleed; improved signal-to-noise ratio for better sensitivity and mass spectral integrity.
- Temperature range: -60 $^{\circ}\text{C}$ to 330/350 $^{\circ}\text{C}$.
- Equivalent to USP G27 and G36 phases.

ID	df	Length	Temp. Limits	qty.	Similar to Part #	cat.#
0.25 mm	0.25 μm	30 m	-60 to 330/350 $^{\circ}\text{C}$	ea.	Agilent 190915-433UI; Phenomenex 7HG-G032-11	13423

For more column dimensions and Rxi phases, visit www.restek.com/rxi



Low-Pressure GC (LPGC) Column Kit

Leverage Your MS Vacuum to Significantly Speed Up Separations

- 3x faster multiresidue pesticides analysis in foods.
- Factory-coupled, leak-free kit makes setting up LPGC as simple as a column change.
- Ideal for speeding up GC-MS and GC-MS/MS methods.
- Integrated transfer line reduces background and stabilization time.

Restek's low-pressure GC column kit has been specifically designed to easily install into your GC-MS or GC-MS/MS system, making it simpler to take advantage of the speed boost that is possible with low-pressure GC-MS (LPGC-MS). This kit is comprised of two factory-coupled columns:

Temp. Limits	Includes	qty.	cat.#
-60 to 340/340 $^{\circ}\text{C}$	Factory-coupled restrictor column (5 m x 0.18 mm ID) and Rtx-5ms analytical column (15 m, 0.53 mm ID, 1 μm plus 1 m integrated transfer line on the outlet end)	kit	11800

Learn more about Restek LPGC at www.restek.com/lpgc



Q-sep QuEChERS Extraction Salts

Free-flowing salts transfer easily and completely. Convenient slim packets fit perfectly into tubes to prevent spills.

Description	Material	Method	qty.	cat.#
Q-sep QuEChERS Extraction Kit	4 g MgSO ₄ , 1 g NaCl with 50 mL Centrifuge Tube	Original unbuffered	50 packets & 50 tubes	25848
Q-sep QuEChERS Extraction Salt Packets Only	4 g MgSO ₄ , 1 g NaCl	Original unbuffered	50 packets	25847
Q-sep QuEChERS Extraction Kit	4 g MgSO ₄ , 1 g NaCl, 1 g TSCD, 0.5 g DHS with 50 mL Centrifuge Tube	European EN 15662	50 packets & 50 tubes	25850
Q-sep QuEChERS Extraction Salt Packets Only	4 g MgSO ₄ , 1 g NaCl, 1 g TSCD, 0.5 g DHS	European EN 15662	50 packets	25849
Q-sep QuEChERS Extraction Kit	6 g MgSO ₄ , 1.5 g NaOAc with 50 mL Centrifuge Tube	AOAC 2007.01	50 packets & 50 tubes	25852
Q-sep QuEChERS Extraction Salt Packets Only	6 g MgSO ₄ , 1.5 g NaOAc	AOAC 2007.01	50 packets	25851

DHS – disodium hydrogen citrate sesquihydrate; MgSO₄ – magnesium sulfate; NaCl – sodium chloride; NaOAc – sodium acetate; TSCD – trisodium citrate dihydrate

For our full line of extraction salts, dSPE cleanup tubes, and more for QuEChERS analysis, **visit www.restek.com/quechers**



Thomson SINGLE StEP Standard Filter Vials

Easy-to-use vials offer fast sample filtration and require only a squeeze of your fingers.

Description	Color	Porosity	qty.	cat.#
Nylon				
Thomson SINGLE StEP Standard Filter Vial	black preslit cap	0.2 µm	100-pk.	25891
	pink preslit cap	0.45 µm	100-pk.	25892
PES (polyethersulfone)				
Thomson SINGLE StEP Standard Filter Vial	grey preslit cap	0.2 µm	100-pk.	25897
PTFE (polytetrafluoroethylene)				
Thomson SINGLE StEP Standard Filter Vials with Screw-Top Caps	green preslit cap	0.2 µm	100-pk.	28307 NEW!
	blue preslit cap	0.45 µm	100-pk.	28306 NEW!
PVDF (polyvinylidene fluoride)				
Thomson SINGLE StEP Standard Filter Vial	red preslit cap	0.2 µm	100-pk.	25895
	yellow preslit cap	0.45 µm	100-pk.	25896

Patent No. 7,790,117

Also available in eXtreme, low-evaporation, and nano versions at **www.restek.com/singlestep**





Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2023 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS2295D-UNV