

Improved Baseline Stability and Sensitivity in GC-MS Using a Novel Column with an Uncoated Interface Section

Laurent Cerveau², Javier Arrebola¹, Ramkumar Dhandapani², Jaap de Zeeuw³,

¹ University of Almeria, Almeria, Spain

² Restek Corporation, Bellefonte, PA, USA

³ CreaVisions, Boxtel, Netherlands

Background

Persistent challenges in GC-MS workflows - particularly baseline instability, elevated background noise, and sensitivity loss - are often linked to stationary-phase degradation and thermal stress near the interface region. These issues become more pronounced at higher transfer-line temperatures and can negatively impact detection of high-boiling-point or thermally sensitive analytes. To address these limitations, a new GC-MS column design incorporating an uncoated interface section was developed to reduce thermally driven degradation processes and stabilize baseline behavior.




Figure 1: Diagram of the GC-MS/MS system with a 5-meter guard column connected to the analytical column. The RMX-5SiI MS column has an integrated transfer line which is deactivated but not coated with stationary phase which reduces conditioning time and bleed.

Methods

The column incorporates a thermally robust, uncoated segment positioned at the GC-MS interface to limit stationary-phase exposure to high temperatures. Performance was evaluated using a diverse set of pesticide analytes covering a range of volatilities and thermal characteristics. Measurements included baseline stabilization time, background signal level, signal-to-noise ratios, and detection behavior in both full-scan and selected-ion monitoring (SIM) modes. Comparative assessments were performed under identical instrumental conditions to isolate the effect of the interface modification.

Integra-Transfer Line




Figure 2: RMX-5SiI MS column with an Integra-Transfer Line, the column is clearly marked indicating the detector end.

Dual

(Integra-Guard and Integra-Transfer Line)




Figure 3: The RMX-5SiI MS column with an integrated guard column and an integrated transfer line, reducing bleed and protecting the column from contaminated samples.




Figure 4: Use of the integrated guard and transferline at different MS transferline temperatures allows for the analysis of high molecular weight compounds at low levels.




Figure 5: Using a higher transferline temperature with the integrated transferline allows for better detection of high molecular weight compounds (Azoxystrobin) that can suffer from low sensitivity as a function of background noise (bleed).




Figure 6: Signal-to-noise ratio (S/N) obtained for extracted Ion Chromatogram (TIC) of Azoxystrobin at different transferline temperatures (m/z 344).




Figure 7: Evaluation of total bleed with after initial conditioning by changing the MS transfer line temperature. Every temperature evaluated resulted in a lower bleed and faster conditioning time using the integrated transfer line.

Observations

Results show that the modified column configuration provides substantially faster baseline stabilization and a measurable reduction in background noise. Improvements in signal-to-noise ratio were consistently observed for both scanning modes, with the greatest gains occurring in high-boiling-point and late-eluting analytes—compounds typically most affected by thermal stress and column bleed. The reduced bleed and improved chromatographic stability suggest lower analyte loss and enhanced reliability in detecting thermally labile species.



Figure 8: Data Courtesy of Research Group Analytical Chemistry of Contaminants, Department of Chemistry and Physics, University of Almería, Almería, Spain.