

Simplifying PFAS analyses with an improved dual bed solid-phase extraction method

Emanuele Ceccon, Sandra Ruiz Perez, Jason Hoisington, Alexis Shelow, Jason Thomas, Colton Myers, Diego A. Lopez
Restek Corporation, Bellefonte, PA, USA

PFAS and EPA 1633

Per- and Polyfluoroalkyl Substances (PFAS) are synthetic organofluorine chemicals that have been linked to a range of health effects, including decrease immune response, vaccine effectiveness, cancer risk, among others.

EPA Method 1633 was finalized early last year, and it represents a comprehensive analysis of 40 target PFAS compounds in different aqueous, solids, and tissue matrices. This method uses a weak anion exchange (WAX) solid phase extraction (SPE), coupled with a dispersive graphitized carbon black (GCB) clean-up step prior to LC-MS/MS analysis on a C18 column.

Statement of Problem

1. Dispersive GCB can be a cumbersome and laborious step that requires further filtration before LC-M/MS analysis.
2. EPA Method 1633 requires analysis of large volume (up to 500 mL) samples of complex aqueous matrices with a varying degree of suspended solids that tend to clog the SPE cartridge, even with addition of glass wool or settling and centrifugation steps. A secondary SPE tube is necessary if clogging occurs.

	#1	#2	#3	#4
Weight of wool used (n=3)	39 mg	150 mg	51 mg	116 mg
Wool Dev.	2 mg	13 mg	7 mg	49 mg
RSD	5%	8%	14%	42%

Figure 1. Inter-cartridge variation from the amount of wool used to fill up to "half the height" of the SPE tube. Manual packing carried out in triplicates by four different lab technicians.

Solution Approach


A dual bed WAX/GCB cartridge that leverages an integrated Filter Aid; therefore, eliminating the need for dispersive carbon clean up, as well as the requirement to hand-pack deactivated glass wool in each tube.



Sample Preparation Workflow

Samples for precision, accuracy, and method detection limit (MDL) studies were prepared in polypropylene bottles using 500 mL deionized water spiked with 25 µL of extracted internal standards (EIS, Wellington Laboratories, Cat. # MPFAC-HIF-ES) as per EPA Method 1633, Section 11.2.4. Four samples for precision and recovery analysis were spiked with 200 µL of native PFAS standards (Wellington Laboratories, Cat. # EPA-1633STK), giving pre-extraction concentrations of 100-2500 ng/L. Seven MDL samples were spiked with 20 µL of a 20:1 dilution of the native standard, giving pre-extraction concentrations of 0.5-12.5 ng/L. Seven blank replicates were prepared for the MDL study as well. The MDL samples were prepared and analyzed over three days.


Figure 2. Sample preparation procedure per EPA 1633 using Thermo Scientific AutoTrace 280 PFAS System.



LC-MS/MS Parameters

Figure 3. Analytical conditions (Waters Xevo TQ-S with ACQUITY Premier UPLC) and chromatogram showing all target PFAS and isotopically labeled compounds.

Analytical Column	1.8 µm Force C18 50x2.1 mm id	
PFAS Delay Column	5 µm PFAS Delay 50x2.1 mm id	
Mobile Phase A	Water, 5mM Ammonium Acetate	
Mobile Phase B	Methanol	
Gradient	Time (min)	%B
	0.00	20
	6.00	95
	6.60	95
	6.61	20
7.50	20	
Flow Rate	0.4 mL/min	
Injection Volume	3 µL	
Column Temp.	40°C	
Ion Mode	Scheduled MRM in negative ESI	



Results and Discussion

Table 1. Results from MDL, Blank, Accuracy, and Precision experiments for native PFAS.

Compound	Abbreviation	MDL (ng/L)	Blank (ng/L)	Accuracy (%)	%RSD
Perfluorobutanoic acid	PFBA	0.34	ND	112	5
Perfluoro-3-methylpropanoic acid	PF3MPA	0.20	ND	109	6
3-Perfluoropropyl propanoic acid	3-PTFA	0.31	ND	95	6
Perfluorobenzoic acid	PFBA	0.26	ND	111	6
Perfluorooctanoic acid	PF8OA	0.19	ND	101	5
Perfluoro-4-methylbutanoic acid	PF4MBA	0.18	ND	111	6
Perfluoro-2-methoxyethanesulfonic acid	PF2ESA	0.15	ND	94	5
Nonafluoro-3,6-dioxoheptanoic acid	NFDHA	0.46	ND	121	5
1H,1H,2H,2H-perfluorohexane sulfonate	4:2 FT6S	0.36	ND	104	5
2H,2H,3H,3H-perfluorooctanoic acid	5:3 FTCA	1.48	ND	102	6
Perfluorodecanoic acid	PF10DA	0.08	ND	113	5
Perfluoropentanoic acid	PF5PA	0.07	ND	122	4
Hexafluoroisopropyl ether diester acid	HFPICDA	0.40	ND	114	12
Perfluorooctanoic acid	PF8OA	0.08	ND	109	5
Perfluorodecanoic acid	PF10DA	0.07	ND	104	9
4,8-Dioxo-3H-perfluorooctanoic acid	ADOMA	0.39	ND	99	9
1H,1H,2H,2H-perfluorooctane sulfonate	6:2 FT8S	0.37	ND	107	4
3-Perfluoropropyl propanoic acid	3-PTFA	2.32	ND	100	6
2H,2H,3H,3H-perfluorooctanoic acid	5:3 FTCA	1.48	ND	102	6
Perfluorodecanoic acid	PF10DA	0.32	ND	104	2
Perfluorooctanoic acid	PF8OA	0.40	ND	114	12
Perfluorooctanesulfonamide	PF8OSA	0.21	ND	81	11
Perfluorodecanesulfonamide	PF10DSA	0.38	ND	104	13
N-methyl perfluorooctanesulfonamide	NMFO8A	0.07	ND	104	13
6-Chloroheptafluoro-3-oxooxane-1-sulfonic acid	6C-PPF3ONS	0.61	ND	85	11
Perfluorooctanesulfonic acid	PF8NS	0.25	ND	79	7
Perfluorodecanesulfonic acid	PF10NS	0.18	ND	117	8
Perfluorooctanesulfonamide	PF8OSA	0.78	ND	112	15
N-methyl perfluorooctanesulfonamide	8:2 FT8S	0.91	ND	113	11
1H,1H,2H,2H-perfluorodecane sulfonate	8:2 FT10S	0.26	ND	105	12
Perfluorodecanoic acid	PF10DA	0.26	ND	105	12
N-methyl perfluorodecanesulfonamide	NMFO10A	0.15	ND	95	12
perfluorooctanesulfonamide	PF8NSA	0.23	ND	93	12
Perfluorooctanesulfonic acid	PF8NS	0.25	ND	94	23
11-Chlorooctadecafluoro-3-oxooxane-1-sulfonic acid	11C-PPF18ONS	0.96	ND	77	25
Perfluorodecanoic acid	PF10DA	0.14	ND	112	15
Perfluorooctanesulfonamide	PF8OSA	0.27	ND	118	15
N-methyl perfluorooctanesulfonamide	NMFO8A	1.11	ND	106	9
N-methyl perfluorodecanesulfonamide	NMFO10A	1.10	ND	111	9
Perfluorodecanoic acid	PF10DA	0.33	ND	108	15
Perfluorodecanesulfonic acid	PF10NS	0.13	ND	89	12
Perfluorodecanesulfonamide	PF10NSA	0.32	ND	102	23

Table 2. Results from Precision and Accuracy experiments for isotope dilution standards.

Compound	Abbreviation	Accuracy (%)	%RSD
Perfluoro-n(1)-2,3,4,6-13C9butanoic acid	13C4-PFBA	88	4
Perfluoro-n(1)-2,3,4,5,13C5pentanoic acid	13C5-PFPA	88	5
Sodium perfluoro-n(1)-2,3,4,5,6,13C6hexanoic acid	13C6-PFHA	97	3
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C2)heptane sulfonate	13C2-4-FT7S	95	8
Perfluoro-n(1)-2,3,4,6-13C9heptanoic acid	13C9-PFHA	122	2
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3,6,6,6,6,13C12)undecanoic acid	13C12-HFPO-DA	92	4
Perfluoro-n(1)-2,3,4,6-13C9heptanoic acid	13C9-PFHA	107	4
Sodium perfluoro-n(1)-2,3,4,6,13C10octanoic acid	13C10-PFHA	102	4
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C2)decane sulfonate	13C2-6-FT8S	83	5
Perfluoro-n(1)-13C10decanoic acid	13C10-PF10A	106	5
Sodium perfluoro-n(1)-13C10decanoic acid	13C10-PFOS	77	7
Perfluoro-n(1)-13C10decanesulfonamide	13C10-PFOSA	94	11
N-methyl-4-3-perfluoro-1-octanesulfonamide	D3-NMFO8A	85	9
N-methyl-4-3-perfluoro-1-octanesulfonamide	D5-NMFO8A	87	11
perfluorooctanesulfonamide	PF8NSA	83	7
Perfluoro-n(1)-2,3,4,5,6,13C7heptanoic acid	13C7-PFPA	80	7
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C2)decane sulfonate	13C2-8-FT8S	84	5
N-methyl-D-7-perfluorooctanesulfonamide	D7-NMFO8E	108	13
N-methyl-D-9-perfluorodecanesulfonamide	D9-NMFO10E	88	17
N-methyl-D-3-perfluoro-1-octanesulfonamide	D3-NMFO8A	84	8
N-methyl-D-5-perfluoro-1-octanesulfonamide	D5-NMFO8A	111	10
Perfluoro-n(1)-2,3,4,5,6,7-13C7heptanoic acid	13C7-PFPA	78	13
Perfluoro-n(1)-2-13C2podecanoic acid	13C2-PF10A	45	24
Perfluoro-n(1)-2-13C2tridecanoic acid	13C2-PF13A	36	21

Figure 4. Accuracy (%) for all target PFAS compounds from four replicate spikes. Accuracy ranged from 77% to 125% of the spiked value, meeting the requirements in Table 5 of EPA Method 1633.




Figure 5. Average volume of 500 mL wastewater extracted before SPE clogging. Comparison between Filter Aid, glass wool, and control. ASTM substitute wastewater (ASTM International, D5905-98[2018]) diluted to ~100 mg/mL suspended solids was used.

Summary

Nondrinking water matrices present unique challenges to labs using SPE sample preparation. In this study, a dual bed WAX/GCB cartridge format exhibited consistency with the acceptance criteria in EPA method 1633. Furthermore, the *Filter Aid* greatly mitigates the clogging risk by this challenging matrix, allowing for 100% extraction of the recommended sample volume under method guidelines.