



Restek Refined

Restek's PLOT Column Family — The Benchmark for Performance!

- Innovative bonding process minimizes particle release.
- More consistent flow means stable retention times.
- Outstanding peak symmetry improves impurity analysis.

Visit us at www.restek.com/petro



RESTEK

Pure Chromatography

www.restek.com

Next Generation of Porous Layer Open Tubular (PLOT) Columns

- Stabilized particle layers improve robustness and reproducibility of retention and flow.
- Compatible with valve switching and Deans switching systems.
- Highly efficient, reproducible analyses; ideal for permanent gases, solvents, and hydrocarbons.
- Innovative manufacturing procedure reduces particle generation and improves performance of PLOT columns.
- Wound on a 7"-diameter, 11-pin cage unless otherwise noted.

Model your PLOT column using our **Pro EZGC Chromatogram Modeler**

Porous layer open tubular (PLOT) columns are very beneficial for solving application problems, especially for the analysis of volatile compounds. PLOT columns have a unique selectivity, allowing for the separation of volatile compounds at ambient temperature. Due to the adsorption mechanism of the stationary phases used in PLOT columns, permanent gases and light hydrocarbons can be resolved at ambient temperature; columns can then be programmed to higher temperatures to elute higher boiling compounds.

Traditional PLOT Columns Offer Poor Stability

The traditional PLOT column is built with a 5–50 μm layer of particles adhered to the tubing walls. Because this layer of particles generally lacks stability, PLOT columns must be used very carefully as particle release is common and can cause unpredictable changes in retention time and flow behavior. Traditional PLOT columns also must generally be used in conjunction with particle traps to prevent the contamination of valves, injectors, and GC detectors. Detectors contaminated with particles typically generate electronic noise, which shows up chromatographically as a spike in the baseline. In extreme cases, detector flow can be obstructed by particle buildup. Particles can also affect valves by becoming lodged in the valve and causing leaks or restricting flow. Figure 1 shows an example of blockage caused by particle accumulation inside a Press-Tight connector.

Restek PLOT Columns Offer Improved Stability to Minimize Particle Release

Restek has developed technology and procedures to manufacture PLOT columns with concentric stabilized adsorption layers. These next-generation PLOT columns show a constant flow behavior (permeability) and have significantly improved mechanical stability, resulting in easier operation, better chromatography, and reduced particle release. Greater particle stability means more reproducible retention times, virtually no spiking, and longer column lifetimes.* This innovative Restek stabilization chemistry is currently applied to all fused silica and metal PLOT columns featured here.

*A particle trap is still recommended to protect valves and detectors.

PLOT Columns Available In:

Fused Silica

Rt Column Phases:


Rt-Alumina BOND/MAPD
Rt-Alumina BOND/ Na_2SO_4
Rt-Alumina BOND/KCl
Rt-Alumina BOND/CFC
Rt-MSieve 5A
Rt-Q-BOND
Rt-QS-BOND
Rt-S-BOND
Rt-U-BOND

Metal

MXT Column Phases:

MXT-Alumina BOND/MAPD
MXT-Alumina BOND/ Na_2SO_4
MXT-MSieve 5A
MXT-Q-BOND

Figure 1: Particles released from traditional PLOT columns can cause blockages.



Consistent Flow Restriction Factor (F) Guarantees Reproducible Flow

Thick layers of particles are difficult to deposit in a homogeneous layer, and in traditionally manufactured PLOT columns, this results in variable coating thicknesses. The positions where the layer is thicker act as restrictions and affect flow (Figure 2). Depending on the number and intensity of these restrictions, traditional PLOT columns often show greater variation in flow restriction than wall coated open tubular (WCOT) columns. In practice, conventional PLOT columns with the same dimensions can differ in flow by a factor of 4 to 6 when operated at the same nominal pressure. For applications where flow is important, such as with Deans switching, the nonreproducible flow behavior of most commercially available PLOT columns is a problem.

In order to measure flow restriction reproducibility, Restek introduced a new factor: the flow restriction factor (F). This factor is based on the retention time of an unretained marker compound as measured on both coated and uncoated tubing using the same back pressure setting (Equation 1). For quality control purposes, methane is used as the marker when evaluating porous polymer columns, and helium is used for testing molecular sieve 5A columns.

Flow restriction factor determination can be used to assess both the degree of column restriction and the reproducibility of the column coating process. Figure 3 shows typical results for PLOT columns manufactured using a conventional process. Because of the difference in flow restriction, individual columns have very different flow characteristics. In contrast, Figure 4 shows results for columns made using our Rt-QS-BOND (bonded porous polymer) PLOT column process. Clearly, our manufacturing process results in greater consistency in both column coating thickness and flow restriction, which results in more stable retention times and better performance in Deans and related flow switching techniques. Flow restriction factors are specified on the certificate of analysis (CofA) included with every Restek PLOT column, and the values are listed on the report.

Figure 2: Inconsistent coating thicknesses result in restrictions that cause significant variation in flow.





Figure 3: Traditional PLOT columns show significant flow variability, indicating inconsistent column coating thicknesses.




Equation 1: Flow restriction factor (F) is used to demonstrate coating consistency.

$$F = \frac{\text{tr}_1 \text{ of unretained component (uncoated tubing)}}{\text{tr}_2 \text{ of unretained component (coated column)}}$$

tr = retention time

Note: F values will always be <1 as the coated column always has more restriction than the uncoated column.

Figure 4: PLOT columns from Restek offer consistent flow restriction, giving more reproducible results column to column.



Restek's PLOT columns are exceptionally robust, featuring concentric stabilized coating layers. They allow for more consistent gas flows and are recommended for applications sensitive to variation in retention time or flow. These PLOT columns are a significant advance in technology and are ideal for efficient, reproducible analyses of permanent gases, solvents, and hydrocarbons.

PLOT Column Phase Cross-Reference: Similar Selectivity

Restek Rt and MXT Columns	Porous Layer	Supelco	Alltech	Agilent (J&W, Varian, Chrompack)	Quadrex
Alumina BOND/Na ₂ SO ₄	Aluminum oxide	Alumina-Sulfate	AT-Alumina	GS-Alumina, CP-Al ₂ O ₃ /Na ₂ SO ₄	_____
Alumina BOND/KCl*	Aluminum oxide	Alumina-Chloride	_____	GS-Alumina KCl, HP PLOT Al ₂ O ₃ , CP-Al ₂ O ₃ /KCl	_____
Alumina BOND/CFC*	Aluminum oxide			unique product	
Alumina BOND/MAPD	Aluminum oxide	_____	_____	Select Al ₂ O ₃ MAPD	_____
Msieve 5A	Molecular sieve 5A	Molsieve 5A	AT-Molesieve	HP PLOT Molesieve, CP-Molesieve 5A	PLT-5A
Q-BOND	100% Divinylbenzene	Supel-Q-PLOT	AT-Q	HP PLOT Q, CP-PoraPlot Q, PoraBond Q	_____
QS-BOND*	Intermediate polarity porous polymer	_____	_____	GS-Q	_____
S-BOND*	DVB vinylpyridine polymer	_____	_____	CP-PoraPlot S	_____
U-BOND*	DVB ethylene glycol-dimethylacrylate polymer	_____	_____	HP PLOT U, CP-PoraPlot U, CP-PoraBond U	_____

*Available only with fused silica tubing.

did you know?

Restek draws our own fused silica tubing and applies our own proprietary stationary phases. By fully managing our production streams, we are able to ensure unparalleled reliability and stability.

Rt-Alumina BOND Columns

Restek Rt-Alumina BOND columns are highly selective for C1–C5 hydrocarbons and separate all saturated and unsaturated hydrocarbon isomers above ambient temperatures. The reactivity of the aluminum oxide stationary phase is minimized by deactivation with inorganic salts, such as KCl or Na₂SO₄, to improve column response for polar unsaturates, like dienes, and the column's sensitivity (or response) ensures linear and quantitative chromatographic analysis for these compounds. Strong bonding minimizes particle generation and release, which allows valve switching with minimal risk to the injection or detection systems. And because they are stable up to at least 200 °C, Rt-Alumina BOND columns can be regenerated to restore full efficiency and selectivity by conditioning at their maximum temperature if water is adsorbed. High capacity and loadability give you exceptionally symmetrical peaks, making these columns ideal for volatile hydrocarbon separations at percent levels, as well as impurity analyses at ppm concentrations. Restek Rt-Alumina BOND PLOT columns are manufactured on fused silica tubing; select phases are also available on metal MXT tubing.

To ensure reproducible retention times and predictable flow behavior column to column, each Rt-Alumina BOND column is extensively tested. A hydrocarbon test mix confirms proper phase retention and selectivity. To calculate *k* (retention or capacity factor), which is a measure of phase retention, 1,3-butadiene is used, while selectivity is measured using retention indices for propadiene and methyl acetylene. The resolution of *trans*-2-butene and 1-butene is also verified and, to measure efficiency, plates per meter are checked using 1,3-butadiene. We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.

Rt-Alumina BOND/Na₂SO₄ Columns (fused silica PLOT)

Na₂SO₄ deactivation

- Acetylene and propadiene elute after butanes.
- Best separation for butene isomers (impurities in butene streams).
- Methyl acetylene elutes after 1,3-butadiene.
- Cyclopropane (impurity in propylene) elutes well before propylene.
- Temperature range: -60 to 200 °C.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	4 µm	30 m	-60 to 200 °C	ea.	19775
0.32 mm	5 µm	30 m	-60 to 200 °C	ea.	19757
	5 µm	50 m	-60 to 200 °C	ea.	19758
0.53 mm	10 µm	30 m	-60 to 200 °C	ea.	19755
	10 µm	50 m	-60 to 200 °C	ea.	19756

i tech tip

Traces of water in the carrier gas and samples will affect the retention and the selectivity of alumina. If exposed to water, the retention times will shorten. The column can be regenerated by conditioning for 15–30 minutes at 200 °C under normal carrier gas flow. Periodic conditioning ensures excellent run-to-run retention time reproducibility.

Unless noted, the maximum programmable temperature for an Rt-Alumina BOND column is 200 °C. Temperatures higher than the stated maximum temperature can cause irreversible changes to the porous layer adsorption properties.


similar phases

Alumina-Sulfate; AT-Alumina; CP-Al₂O₃/Na₂SO₄, GS-Alumina

Rt-Alumina BOND columns show unique retention characteristics for hydrocarbons.

also available! Metal MXT PLOT Columns

Refinery Gas on Rt-Alumina BOND (Na₂SO₄)



Peaks

1. Methane
2. Ethane
3. Ethylene
4. Propane
5. Propylene
6. Isobutane
7. *n*-Butane
8. Propadiene
9. Acetylene
10. *trans*-2-Butene
11. 1-Butene
12. Isobutylene
13. *cis*-2-Butene
14. *iso*-Pentane
15. *n*-Pentane
16. 1,3-Butadiene
17. *trans*-2-Pentene
18. 2-Methyl-2-butene
19. 1-Pentene
20. *cis*-2-Pentene
21. Hexanes

Column

Rt-Alumina BOND/Na₂SO₄, 50 m, 0.53 mm ID, 10 µm (cat.# 19756)
 Sample Injection
 Inj. Vol.: 10 µL split
 Liner: Taper (2 mm) (cat.# 20795)
 Inj. Temp.: 200 °C
 Split Vent
 Flow Rate: 80 mL/min
 Oven
 Oven Temp.: 45 °C (hold 1 min) to 200 °C at 10 °C/min (hold 3.5 min)
 Carrier Gas: H₂, constant pressure (8.0 psi, 55.2 kPa)
 Linear Velocity: 74 cm/sec @ 45 °C
 Detector: FID @ 200 °C

Rt-Alumina BOND/KCl Columns (fused silica PLOT)

KCl deactivation

- Restek's lowest polarity alumina column.
- Low moisture sensitivity reduces the need for frequent regeneration.
- Acetylene elutes before *n*-butane.
- Methyl acetylene (impurity in 1,3-butadiene) elutes before 1,3-butadiene.
- Temperature range: -60 to 200 °C.

similar phases

Alumina-Chloride; CP-Al₂O₃/KCl, GS-Alumina KCl, HP-PLOT Al₂O₃/KCl


Rt-Alumina BOND columns show unique retention characteristics for hydrocarbons.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	4 µm	30 m	-60 to 200 °C	ea.	19776
0.32 mm	5 µm	30 m	-60 to 200 °C	ea.	19761
	5 µm	50 m	-60 to 200 °C	ea.	19762
0.53 mm	10 µm	30 m	-60 to 200 °C	ea.	19759
	10 µm	50 m	-60 to 200 °C	ea.	19760

did you know?

All Restek PLOT columns come standard on a 7"-diameter, 11-pin cage. Metal MXT columns are also available coiled to 3.5" diameter for select phases by adding the suffix -273 to the part number. If you need more information, please call your local Restek representative.

Ethylene and C1-C5 Hydrocarbons by ASTM D6159-97 on Rt-Alumina BOND/KCl, Rtx-1



Column Rt-Alumina BOND/KCl *, 50 m, 0.53 mm ID, 10 µm (cat.# 19760)
Sample Injection Ethylene gas plus C1 through C5 hydrocarbons
Inj. Vol.: 1 µL split
Liner: 2 mm splitless (cat.# 20712)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 60 mL/min
Oven
Oven Temp.: 35 °C (hold 2 min) to 190 °C at 4 °C/min (hold 15 min)

Carrier Gas He, constant pressure (8.0 psi, 55.2 kPa)
Linear Velocity: 25.4 cm/sec @ 35 °C
Detector FID @ 200 °C
Make-up Gas
Type: N₂
Data Rate: 20 Hz
Instrument HP5890 GC
Notes * Rt-Alumina BOND/KCl, 50 m, 0.53 mm ID, 10.0 µm (cat.# 19760) in series with an Rtx-1, 30 m, 0.53 mm ID, 5.0 µm (cat.# 10179) connected using a universal Press-Tight connector (cat.# 20401).
 (conditions as per ASTM D6159-97)

Rt-Alumina BOND/CFC Columns (fused silica PLOT)

- Improved inertness for chlorofluorocarbon (CFC) compounds.
- Highly selective alumina-based column, separates most CFCs.
- High retention and capacity for CFCs.
- Temperature range: -60 to 200 °C.

The Alumina BOND/CFC adsorbent is ideal for retaining volatile halogenated compounds, especially CFCs (chlorinated fluorocarbons). It offers high selectivity, allowing a wide range of CFC isomers to be resolved at above ambient temperatures. The Rt-Alumina BOND/CFC column is thoroughly deactivated to reduce the reactivity of alumina. Even though there is still some residual reactivity for some mono- or di-substituted CFCs, the majority of these compounds can be accurately quantified from mainstream processes or in impurity analyses.

Restek
exclusive!




Especially when valve switching or backflushing is used, Restek recommends using particle traps to help prevent detector spikes and valve rotor scratches. Visit www.restek.com/plot for specialized PLOT column particle traps.

ID	df	Length	Temp. Limits	qty.	cat.#
0.53 mm	10 µm	30 m	-60 to 200 °C	ea.	19763

Impurity Analysis of 1,1,1,2-Tetrafluoroethane (CFC-134a) on Rt-Alumina BOND/CFC

Rt-Alumina BOND/CFC columns provide the best selectivity for impurity analysis of CFCs.



Column Rt-Alumina BOND/CFC, 30 m, 0.53 mm ID (cat.# 19763)
Sample 1,1,1,2-Tetrafluoroethane
Conc.: Neat
Injection
Inj. Vol.: 500 µL split
Oven
Oven Temp.: 80 °C (hold 6 min) to 140 °C at 10 °C/min (hold 2 min)
Carrier Gas He
Detector FID
Notes Gas sampling, purity analysis

GC_GN1155

Note that tailing peaks are common in CFC analyses due to overloading normally employed for this type of work.



For more chromatograms, search our extensive library at
www.restek.com/chromatograms



also available! Metal MXT PLOT Columns

See page 15 for more
information.

similar phases

Select Al₂O₃ MAPD

Rt-Alumina BOND/MAPD Columns (fused silica PLOT)

- Optimized deactivation produces maximum response when analyzing trace levels of acetylene, methyl acetylene, and propadiene.
- Stable response factors make this column ideal for process-type applications where recalibration must be minimized.
- High loadability reduces peak tailing and improves separations.
- Extended temperature range up to 250 °C for fast elution of high molecular weight (HMW) hydrocarbons and accelerated column regeneration following exposure to water.
- Temperature range: -60 to 250 °C.

Rt-Alumina BOND/MAPD PLOT columns are made specifically for the analysis of petrochemicals and downstream products such as ethylene, propylene, butylenes, and butadiene.

Also available in metal MXT-Alumina columns!

Restek's R&D chemists have optimized the deactivation technology applied to our Rt-Alumina BOND/MAPD column for improved analysis of trace concentrations of polar hydrocarbons like acetylene, methyl acetylene, and propadiene in hydrocarbon streams containing higher levels of C1-C5 hydrocarbons. Our alumina PLOT deactivation produces an incredibly inert column that offers superior reproducibility and stable response factors to maximize the number of analyses before recalibration is required. Its high sample capacity reduces peak tailing, thereby improving the separation of target compounds. In addition, a 250 °C maximum operating temperature lets you more quickly elute hydrocarbons up to dodecane and reduces regeneration time when the column is exposed to water from samples or carrier gases.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	4 µm	30 m	-60 to 250 °C	ea.	19781
	5 µm	30 m	-60 to 250 °C	ea.	19779
0.32 mm	5 µm	50 m	-60 to 250 °C	ea.	19780
0.53 mm	10 µm	50 m	-60 to 250 °C	ea.	19778

1,3-Butadiene on Rt-Alumina BOND/MAPD (Purity Analysis)

Peaks

1. Isobutane
2. *n*-Butane
3. Propadiene
4. *trans*-2-Butene
5. 1-Butene
6. Isobutene
7. *cis*-2-Butene
8. Isopentane
9. *n*-Pentane
10. 1,2-Butadiene
11. 1,3-Butadiene
12. Methyl acetylene

Column Rt-Alumina BOND/MAPD, 50 m, 0.53 mm ID, 10.0 µm (cat.# 19778)
Sample Crude 1,3-butadiene


Injection
Inj. Vol.: 10 µL split
Liner: 2.0 mm ID straight inlet liner (cat.# 20712)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 100 mL/min

Oven
Oven Temp.: 70 °C (hold 5 min) to 200 °C at 10 °C/min
(hold 0 min)

Carrier Gas He, constant pressure (20 psi, 137.9 kPa)
Detector FID @ 250 °C

Make-up Gas 30 mL/min
Flow Rate:

Make-up Gas N₂
Type:
Data Rate: 20 Hz
Instrument HP5890 GC



Rt-Alumina BOND/MAPD PLOT columns are made specifically for the analysis of petrochemicals and downstream products such as ethylene, propylene, butylenes, and butadiene.

Molecular Sieve 5A PLOT Columns

Restek's molecular sieve 5A PLOT columns are designed for efficient separation of argon/oxygen and other permanent gases, including carbon monoxide. Special coating and deactivation procedures ensure chromatographic efficiency and the integrity of the porous layer coating. Molecular sieves have very high retention, allowing separations of permanent gases at temperatures above ambient. Our deactivation technology also allows carbon monoxide to elute as a sharp peak. Additionally, our unique immobilization process guarantees that the uniform particles remain adhered to the tubing—even after continuous valve cycling.

Rt-Msieve 5A Columns (fused silica PLOT)

- Improve accuracy with sharp, symmetrical peaks for argon, oxygen, and carbon monoxide.
- Easily separate permanent gases at temperatures above ambient.
- Restek PLOT technology reduces particle release, improving flow reproducibility and reducing downtime for maintenance.
- Temperature range: -100 to 300 °C.

We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.

ID	df	Length	Temp. Limits	qty.	Similar to Part #	cat.#
0.25 mm	20 µm	15 m	-100 to 300 °C	ea.		19773
	30 µm	15 m	-100 to 300 °C	ea.	Agilent 19091P-MS7; Scion/Bruker/Varian CP7535	19720
0.32 mm	30 µm	30 m	-100 to 300 °C	ea.	Agilent 19091P-MS8; Scion/Bruker/Varian CP7536	19722
	50 µm	15 m	-100 to 300 °C	ea.	Agilent 19095P-MS9; Scion/Bruker/Varian CP7537	19721
0.53 mm	50 µm	30 m	-100 to 300 °C	ea.	Agilent 19095P-MS0; Scion/Bruker/Varian CP7538	19723

i tech tip


Because molecular sieve materials are very hydrophilic, they will adsorb water from the sample or carrier gas. Water contamination can have a detrimental effect on peak symmetry and can reduce the resolution of all compounds. If water contamination occurs, reactivate your Rt-Msieve 5A PLOT column by conditioning at 300 °C with dry carrier gas flow for 3 hours.

similar phases

AT-Molsieve; CP-Molsieve 5A, HP PLOT Molsieve; Molsieve 5A; PLT-5A

Rt-Msieve 5A PLOT columns are designed for efficient separation of Ar/O₂ and other permanent gases, including CH₄, C₂H₆, and CO.

Separation of Argon/Oxygen and Other Permanent Gases on Rt-Msieve 5A



Peaks	Conc. (µg/mL)
1. Hydrogen	40
2. Argon	30
3. Oxygen	50
4. Nitrogen	50
5. Methane	40
6. Carbon monoxide	50

Column Rt-Msieve 5A, 30 m, 0.53 mm ID, 50 µm (cat.# 19723)
Sample Permanent gases
Injection Sample valve
Sample Loop Vol.: 5 µL
Valve Name: 6-port Valco valve
Inj. Temp.: 200 °C
Valve Temp.: Ambient
Oven
Oven Temp.: 27 °C (hold 5 min) to 100 °C at 10 °C/min (hold 5 min)
Carrier Gas He, constant flow
Flow Rate: 5.0 mL/min
Detector Valco helium ionization detector @ 150 °C



also available!

Metal MXT PLOT Columns

See page 14 for more information.

Porous Polymer Columns

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications where moisture is of major concern.

similar phases

AT-Q; CP-PoraBOND Q,
CP-PoraPLOT Q, HP PLOT Q;
Supel-Q-PLOT

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications where moisture is of major concern.

Porous polymers are unique, highly retentive stationary phases with a wide application range that are able to elute both polar and nonpolar compounds. They are very hydrophobic, so water has no impact on retention times and even elutes as a good chromatographic peak. The Q-BOND is our most nonpolar and widely used porous polymer column; functional groups can be added to increase polarity (i.e., QS-, S-, and U-BOND). The process used to manufacture porous polymer PLOT columns causes the particles to adhere strongly to the walls of the tubing, so there is virtually no particle generation. You get reproducible performance from column to column, including selectivity and flow.

Rt-Q-BOND Columns (fused silica PLOT)

100% divinylbenzene

- Nonpolar PLOT column incorporating 100% divinylbenzene.
- Excellent for analysis of C1 to C3 hydrocarbons as well as isomers and alkanes up to C10.
- High retention for CO₂ simplifies gas analysis; CO₂ and methane separated from O₂/N₂/CO. (Note: O₂/N₂/CO not separated at ambient temperature.)
- Use for analysis of oxygenated compounds and solvents.
- Maximum temperature of 300 °C.

also available! Metal MXT PLOT Columns

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	8 µm	15 m	-60 to 280/300 °C	ea.	19764
	8 µm	30 m	-60 to 280/300 °C	ea.	19765
	10 µm	15 m	-60 to 280/300 °C	ea.	19743
0.32 mm	10 µm	30 m	-60 to 280/300 °C	ea.	19744
	20 µm	15 m	-60 to 280/300 °C	ea.	19741
0.53 mm	20 µm	30 m	-60 to 280/300 °C	ea.	19742

Restek porous
polymer PLOT
columns cover a
wide range of
polarities

least polar

Q-BOND


QS-BOND

S-BOND

U-BOND

most polar

Solvent Mixture on Rt-Q-BOND



Peaks

1. Methanol
2. Ethanol
3. Acetonitrile
4. Acetone
5. Dichloromethane
6. 1,1-Dichloroethene
7. Nitromethane
8. *trans*-1,2-Dichloroethylene
9. *cis*-1,2-Dichloroethylene
10. Tetrahydrofuran
11. Chloroform
12. Ethyl acetate
13. 1,2-Dichloroethane
14. 1,1,1-Trichloroethane
15. Benzene
16. 1,2-Dimethoxyethane
17. Trichloroethylene
18. 1,4-Dioxane
19. Pyridine
20. Dimethylformamide
21. Methylcyclohexane
22. Toluene
23. 2-Hexanone
24. Chlorobenzene
25. Ethylbenzene
26. *m*-Xylene
27. *p*-Xylene
28. *o*-Xylene

Column Rt-Q-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19742)
Sample Solvent mixture
Injection 1.0 µL split
Liner: Splitless taper (4 mm) (cat.# 20798)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 100 mL/min

Oven
Oven Temp.: 120 °C to 240 °C at 5 °C/min (hold 5.0 min)
Carrier Gas H₂, constant pressure (4.2 psi, 29.0 kPa)
Linear Velocity: 40 cm/sec @ 120 °C
Detector FID @ 240 °C

Rt-QS-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer

- Intermediate polarity porous polymer PLOT column incorporating low 4-vinylpyridine.
- Separates ethane, ethylene, and acetylene to baseline.
- Temperature range: -60 to 250 °C.


similar phases

GS-Q

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications in which moisture is of major concern.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	8 µm	15 m	-60 to 250 °C	ea.	19767
	8 µm	30 m	-60 to 250 °C	ea.	19768
0.32 mm	10 µm	15 m	-60 to 250 °C	ea.	19739
	10 µm	30 m	-60 to 250 °C	ea.	19740
0.53 mm	20 µm	15 m	-60 to 250 °C	ea.	19737
	20 µm	30 m	-60 to 250 °C	ea.	19738

Refinery Gas Mixture on Rt-QS-BOND



GC_PC01143

Peaks

1. Air
2. Methane
3. Carbon dioxide
4. Ethylene
5. Acetylene
6. Ethane
7. Propylene
8. Propane
9. Propadiene
10. Isobutane
11. Isobutylene

12. 1-Butene
13. 1,3-Butadiene
14. *n*-Butane
15. *cis*-2-Butene
16. *trans*-2-Butene
17. Isopentane
18. 1-Pentene
19. 2-Methyl-2-butene
20. *n*-Pentane
21. *cis*-2-Pentene
22. *n*-Hexane

Column

Rt-QS-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19738)
Refinery gas standard

Sample Injection
Inj. Vol.: 20 µL split
Liner: 2 mm (cat.# 20712)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 35 mL/min

Oven
Oven Temp.: 40 °C (hold 2 min) to 225 °C at 15 °C/min (hold 5 min)

Carrier Gas
He, constant pressure (11.5 psi, 79.3 kPa)

Linear Velocity: 68 cm/sec @ 40 °C

Detector
TCD @ 225 °C

Make-up Gas
Type: He
Data Rate: 20 Hz
Sensitivity Mode: He/H₂
Instrument HP5890 GC

Rt-S-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer

- Midpolarity porous polymer PLOT column, incorporating high 4-vinylpyridine.
- Use for the analysis of nonpolar and polar compounds.
- Temperature range: -60 to 250 °C.


similar phases

CP-PoraPLOT S

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications in which moisture is of major concern.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	8 µm	15 m	-60 to 250 °C	ea.	19769
	8 µm	30 m	-60 to 250 °C	ea.	19770
0.32 mm	10 µm	30 m	-60 to 250 °C	ea.	19748
	20 µm	15 m	-60 to 250 °C	ea.	19745
0.53 mm	20 µm	30 m	-60 to 250 °C	ea.	19746

Solvent Mixture on Rt-S-BOND



Column Rt-S-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19746)
Sample Solvent mixture
Injection
Inj. Vol.: 1.0 µL split
Liner: Taper (4 mm) (cat.# 20798)
Inj. Temp.: 200 °C
Split Vent Flow Rate: 100 mL/min
Oven
Oven Temp.: 120 °C to 220 °C at 5 °C/min (hold 5.0 min)
Carrier Gas H₂, constant pressure (4.2 psi, 29.0 kPa)
Linear Velocity: 40 cm/sec @ 120 °C
Detector FID @ 220 °C

Rt-U-BOND Columns (fused silica PLOT)

divinylbenzene ethylene glycol/dimethylacrylate

- Restek's highest polarity porous polymer column.
- Polar PLOT column, incorporating divinylbenzene ethylene glycol/dimethylacrylate.
- Highly inert for the analysis of polar and nonpolar compounds.
- Common applications include trace H₂S, COS, and mercaptans in hydrocarbon streams.
- Temperature range: -60 to 190 °C.


similar phases

CP-PoraBOND U, CP-Pora-PLOT U, HP-PLOT U

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications in which moisture is of major concern.

ID	df	Length	Temp. Limits	qty.	cat.#
0.25 mm	8 µm	15 m	-60 to 190 °C	ea.	19771
	8 µm	30 m	-60 to 190 °C	ea.	19772
0.32 mm	10 µm	15 m	-60 to 190 °C	ea.	19751
	10 µm	30 m	-60 to 190 °C	ea.	19752
0.53 mm	20 µm	15 m	-60 to 190 °C	ea.	19749
	20 µm	30 m	-60 to 190 °C	ea.	19750

Formaldehyde on Rt-U-BOND



Peaks

1. Air
2. Carbon dioxide
3. Formaldehyde
4. Water
5. Methanol

Column Rt-U-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19750)
Sample Formaldehyde (manual headspace)
Injection
Inj. Vol.: 10 µL split (split ratio 10:1)
Liner: 2 mm split Precision liner w/wool (cat.# 20823)
Inj. Temp.: 200 °C
Split Vent Flow Rate: 40 mL/min
Oven
Oven Temp.: 100 °C (hold 1 min) to 150 °C at 25 °C/min (hold 3 min)
Carrier Gas He, constant pressure (7.7 psi, 53.1 kPa)
Linear Velocity: 39 cm/sec @ 100 °C
Detector TCD @ 200 °C
Make-up Gas Type: He
Data Rate: 20 Hz
Sensitivity Mode: He/H₂
Instrument HP5890 GC

MXT-Alumina Columns



Molecular Sieve 5A PLOT Columns

Restek's molecular sieve 5A PLOT columns are designed for efficient separation of argon/oxygen and other permanent gases, including carbon monoxide. Special coating and deactivation procedures ensure chromatographic efficiency and the integrity of the porous layer coating. Molecular sieves have very high retention, allowing separations of permanent gases at temperatures above ambient. Our deactivation technology also allows carbon monoxide to elute as a sharp peak. Additionally, our unique immobilization process guarantees that the uniform particles remain adhered to the tubing—even after continuous valve cycling.

Because molecular sieve materials are very hydrophilic, they will adsorb water from the sample or carrier gas. Water contamination can have a detrimental effect on peak symmetry and can reduce the resolution of all compounds. If water contamination occurs, reactivate your Restek Msieve 5A PLOT column by conditioning at 300 °C with dry carrier gas flow for 3 hours.

MXT-Msieve 5A Columns (Siltek-treated stainless-steel PLOT)

Advantages of metal MXT PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Rugged material withstands rough handling and shock.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5"-coil diameter or 7"-diameter, 11-pin cage.*

*We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.

ID	df	Length	Temp. Limits	Column Config	qty.	cat.#
0.53 mm	50 µm	30 m	-100 to 300 °C	7" 11-pin cage	ea.	79723
	50 µm	30 m	-100 to 300 °C	3.5" Coil	ea.	79723-273

MXT-Alumina BOND/Na₂SO₄ Columns (Siltek-treated stainless-steel PLOT)

Advantages of metal MXT PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Rugged material withstands rough handling and shock.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5"-coil diameter or 7"-diameter, 11-pin cage.*

*We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.

ID	df	Length	Temp. Limits	Column Config	qty.	cat.#
0.53 mm	10 µm	30 m	-60 to 200 °C	7" 11-pin cage	ea.	79714
	10 µm	30 m	-60 to 200 °C	3.5" Coil	ea.	79714-273

MXT-Alumina BOND/MAPD Columns (Siltek-treated stainless-steel PLOT)

- Optimized deactivation produces maximum response when analyzing trace levels of acetylene, methyl acetylene, and propadiene.
- Stable response factors make this column ideal for process-type applications where recalibration must be minimized.
- High loadability reduces peak tailing and improves separations.
- Extended temperature range up to 250 °C for fast elution of high molecular weight (HMW) hydrocarbons and accelerated column regeneration following exposure to water.
- Temperature range: -60 to 250 °C.

MXT-Alumina BOND/MAPD PLOT columns are made specifically for the analysis of petrochemicals and downstream products, such as ethylene, propylene, butylenes, and butadiene!

Also available in fused silica Rt-Alumina columns!

Restek's R&D chemists have optimized the deactivation technology applied to our MXT-Alumina BOND/MAPD column for improved analysis of trace concentrations of polar hydrocarbons like acetylene, methyl acetylene, and propadiene in hydrocarbon streams containing higher levels of C1-C5 hydrocarbons. Our new Alumina PLOT deactivation produces an incredibly inert column that offers superior reproducibility and stable response factors to maximize the number of analyses before recalibration is required. Its high sample capacity reduces peak tailing, thereby improving the separation of target compounds. In addition, a 250 °C maximum operating temperature lets you more quickly elute hydrocarbons up to dodecane and reduces regeneration time when the column is exposed to water from samples or carrier gases. We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.



ID	df	Length	Temp. Limits	Column Config	qty.	cat.#
0.53 mm	10 µm	30 m	-60 to 250 °C	7" 11-pin cage	ea.	79728
	10 µm	30 m	-60 to 250 °C	3.5" Coil	ea.	79728-273

MXT-Q-BOND Columns (Siltek-treated stainless-steel PLOT)

Advantages of metal MXT PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Rugged material withstands rough handling and shock.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5"-coil diameter or 7"-diameter, 11-pin cage.*

*We do not recommend recoiling PLOT columns as this may cause particles to dislodge from the side of the tubing.

ID	df	Length	Temp. Limits	Column Config	qty.	cat.#
0.25 mm	8 µm	15 m	-60 to 300 °C	7" 11-pin cage	ea.	79718
	8 µm	15 m	-60 to 300 °C	3.5" Coil	ea.	79718-273
0.53 mm	20 µm	30 m	-60 to 280/300 °C	7" 11-pin cage	ea.	79716
	20 µm	30 m	-60 to 280/300 °C	3.5" Coil	ea.	79716-273



20394

MXT Low-Dead-Volume Connector Kits for Metal Columns

- Connect a guard column/transfer line to an MXT stainless-steel column.
- Low thermal mass tracks rapid oven temperature programming.
- Stainless-steel ferrules and nuts.
- Available in “Y” and union configurations.
- Siltek treatment ensures ultimate inertness.

Each kit contains the MXT union, two 1/32-inch ferrules and nuts.

Description	Includes	Fits Column ID	Union Bore	qty.	cat. #
MXT Low-Dead-Volume Connector Kit	MXT union, 1/32-inch ferrules (2) and nuts (2)	For 0.28 mm ID MXT Columns	0.25 mm	kit	20397
	MXT union, 1/32-inch ferrules (2) and nuts (2)	For 0.18/0.25/0.32 mm ID MXT Columns	0.25 mm	kit	20536
	MXT union, 1/32-inch ferrules (2) and nuts (2)	For 0.53 mm ID MXT Columns	0.25 mm	kit	20394



19754

PLOT Column Particle Trap

- Includes two Press-Tight connectors, a 2.5 m column, and high-temperature string.
- Protects detector and valves; connects between column and detector or valve.
- Eliminates detector spikes and scratches in valve rotors.

The technology used to adhere particles in PLOT columns is excellent; however, it is still possible for particles to dislodge when extreme pressure shocks and gas flow changes occur. This sometimes happens when valve switching or backflushing is used. In those cases, using particle traps is recommended.

Description	Includes	Size	qty.	cat. #
PLOT Column Particle Trap	Press-Tight connectors (2); 2.5 m, 0.25 mm ID column; high-temperature string	0.25 mm ID	ea.	19774
	Press-Tight connectors (2); 2.5 m, 0.32 mm ID column; high-temperature string	0.32 mm ID	ea.	19753
	Press-Tight connectors (2); 2.5 m, 0.53 mm ID column; high-temperature string	0.53 mm ID	ea.	19754

Learn more at www.restek.com/petro

RESTEK
Pure Chromatography

Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# PCSS1163K-UNV