

Rapid quantification of PFAS compounds in water utilizing a weak anion exchange coated blade spray mass spectrometry (CBS-MS)

Thomas Kane¹; Ryan Micklitsch¹; Shane Stevens¹; Diego Lopez¹; Mathew Lininger¹; Gary Stidsen¹; Tracey Peters¹

¹Restek Corporation, Bellefonte, PA

PFAS and EPA 1633

PFAS compounds are synthetic organofluorine chemicals that have been produced since the 1940's and have been used in a variety of consumer products and industries.

EPA 1633 analyzes **40 PFAS compounds** in aqueous, solid, and tissue samples. This method utilizes a weak anion exchange (WAX) SPE cartridge for clean-up, with LC-MS/MS analysis on a C18 column.

Statement of Problem

EPA Method 1633 analyzes **large volume** samples of **complex matrices** with significant solids content.

- Sample preparation is difficult, time consuming, and has high potential for contamination.
- LC-MS analysis requires **long analysis times** and is also rife with potential sources of contamination.

The Solution = Coated Blade Spray

Coated Blade Spray (CBS) is an **easy to use** open-bed SPME device that can be directly interfaced with mass spectrometry (MS) instrumentation as a substrate electrospray ionization device.

- No cartridge, no clogging!**
- Less than **15 second** analysis time.
- Green:** consumes a fraction of the solvent needed for a traditional LC-MS analysis.
- CBS is an excellent screening analysis** to reduce laboratory workload.





Figure 1. CBS-MS Total Solution employing an Opentrons OT-2 and a Thermo TSQ Altis.

Coated Blade Spray is a simple process!

Capture the target analytes on a blade, then **release** them into the mass spectrometer.




Experimental Conditions

40 PFAS compounds from the EPA 1633 list were extracted from **substitute wastewater** (ASTM D5905-98) with a **WAX-coated blade** and analyzed with an Opentrons OT-2 liquid handler attached to a Thermo TSQ Altis.

Capture		Release	
Sample Volume:	500 µL	Spray Solvent:	9:1 MeOH:CHCl ₃ , 1% NH ₄ OH
Extraction Time:	60 minutes	Elution Volume:	10 µL
Rinse:	Water	Elution Time:	3 s
Rinse Time:	10 s	Electrospray voltage:	-2800 V
Shaker Speed:	1200 rpm		

> 70 Transitions in a five second spray

A **five second spray** is enough to achieve **quantitative** results, even for large panels. This work monitored over 70 MRM transitions utilizing a 2 ms dwell time to achieve **20 scans** per compound.



Results and Discussion

The table below shows signal-to-noise (S:N) values at six different concentrations. "Cal 1" represents samples at one quarter the minimum concentration level described in EPA 1633. **34 compounds at this level yield S:N values greater than three.**

Compound	Cal 1	Cal 2	Cal 3	Cal 4	Cal 5	Cal 6
PFEESA	66	112	246	411	596	1,277
6:2FTS	52	91	211	325	451	1,532
8:2FTS	50	91	210	318	460	2,063
9CI-PF3ONS	49	94	206	321	463	1,157
11CI-PF3Ouds	39	69	155	266	360	1,069
4:2FTS	36	57	128	201	281	797
PFHxS	29	48	106	173	243	509
PFPeS	25	43	94	161	216	444
ADONA	21	37	84	134	182	556
PFOSHA	21	34	75	119	170	502
PFHpA	20	28	66	119	175	939
PFHpS	18	33	74	111	154	348
PFnPA	16	31	71	138	189	1,177
3,6-OPFOA	16	23	54	81	112	343
PFBS	14	24	51	84	118	236
PF4OPeA	13	21	46	72	105	301
N-EtFOSA	12	20	48	76	107	801
PFNA	12	22	50	73	105	403
L-PFDS	11	18	43	63	87	219
HFPO-DA	10	16	35	59	90	288
PFDA	9	16	35	51	74	315
PFHpA	7	15	34	59	82	607
PFTrDA	7	10	24	36	48	199
PFHpA	7	12	26	42	58	185
PFPeA	7	11	22	35	51	141
PFUnA	6	11	23	34	47	214
PFTeDA	6	9	21	28	39	176
PFHxA	6	9	20	32	45	138
PFDoA	5	8	16	27	35	148
L-PFNS	5	8	15	24	32	93
N-MeFOSA	5	9	18	28	42	303
PFOA	5	8	16	26	34	120
PFOS	4	7	14	21	30	72
FOSA	3	5	12	19	30	220
PFBA	3	4	7	11	15	38
L-PFDoS	2	4	7	15	17	44
N-MeFOSA	2	2	5	8	12	59
N-EtFOSA	1	1	3	4	5	23

Figure 2. Table of S:N values at six concentrations (left). Calibration curves for four prominent compounds (right).

Quantitative results were promising:

- 19 compounds with $R^2 > 0.99$
- 34 compounds with single or double digit ppt MDL blank

During tuning two compounds underwent in-source fragmentation, resulting in MRMs indistinguishable from N-MeFOSA and N-EtFOSA. Therefore, results for these pairs of compounds are summed together.

CBS-MS quantitatively identified **37 PFAS compounds** extracted from substitute wastewater at **1/1000th** the typical sample volume of EPA Method 1633.

Acknowledgements

Janusz Pawliszyn is the inventor of CBS and a collaborator of Restek.

- Gómez-Ríos, G.A.; Pawliszyn, J., *Angew. Chem.*, **2014**, 53, 14503-14507
- Gómez-Ríos, G.A., *The Analytical Scientist*, **2017**, Blade Stunner

WAX sorbent was purchased from Suzhou Nanomicro Technology Co. The authors declare no competing financial interest.