

Simultaneous Analysis of Ultrashort-Chain to Long-Chain (C1 to C10) and Alternative PFAS in Human Plasma and Serum

Shun-Hsin Liang, Justin Steimling; Restek Corporation

Introduction

Ultrashort-chain (USC) per- and polyfluoroalkyl substances (PFAS) are small and very polar compounds with carbon chain lengths shorter than C4 (**Figure 1**). Their ubiquitous and high levels of occurrence in environmental aquatic systems have raised significant concern in conjunction with long-chain PFAS contamination. Measuring USC PFAS in blood can not only monitor human exposure but also serves as a valuable tool for studying the potential risks associated with USC PFAS exposure. The high polarity of USC PFAS poses a challenge to current analytical practices based on the reverse-phased liquid chromatography, primarily due to insufficient chromatographic retention. In this study, a simple and reliable workflow was developed for the simultaneous analysis of C1 to C10 perfluoroalkyl carboxylic and sulfonic acids, along with four alternative PFAS, in human plasma and serum. The samples underwent a single-step protein precipitation procedure and were analyzed with a user-friendly LC-MS/MS method, implementing a polar-embedded LC column. Method accuracy and precision were evaluated with fortified fetal bovine serum. The method’s validity was confirmed by accurately measuring targeted PFAS with known concentrations in NIST standard reference human plasma (1950) and serum (1957).

Figure 1: Structures of C1 to C3 PFAS

