

Simultaneous Determination of *Alternaria* Toxins, Ergot Alkaloid Epimers, and Other Major Mycotoxins in Various Food Matrices by LC-MS/MS

Diego A. Lopez, Shun-Hsin Liang, and Justin Steimling; Restek Corporation, Bellefonte, PA, USA.

Abstract & Introduction

Various food commodities are vulnerable to different types of fungal pathogens and could be contaminated with differential classes of mycotoxins as a result. It is ideally to implement a generic method for simultaneous determination of multi-mycotoxins in different food matrices or agricultural products. In this study, a simplified sample preparation procedure and a reliable LC-MS/MS analytical method was developed for comprehensive measurement of 37 regulated and emerging mycotoxins including 5 *Alternaria* toxins, 6 major ergot alkaloids and their corresponding epimers. Four different food matrices (baby wheat cereal, peanut, tomato puree, and blended flour) were chosen for method validation to demonstrate the applicability of this analytical method to a wide range of food types. Sample extraction was performed using a formic acid-acidified 80:20 acetonitrile:water solution followed by extract dry-down and reconstitution in a 50:50 water:methanol solution for injection analysis on a Biphenyl LC column. Chromatographic analysis was performed using MS-friendly acidic mobile phases and completed with a short 11-minute cycling time for proper separation of ergot alkaloid epimers. Method accuracy and precision was evaluated by fortification of food samples at 3 different levels. Accurate quantification was achieved using matrix-matched calibration standards at the range of 0.4 to 400 µg/kg. The recoveries of all mycotoxins (except citrinin) in fortified samples were from 70% to 120%, and the relative standard deviation was less than 20%. The established workflow was simple and fast for multi-mycotoxin determination with a unique benefit of simultaneous analysis of *Alternaria* toxins and ergot alkaloids. Furthermore, a novel inert Biphenyl LC column demonstrated the high degree of Non-Specific Binding (NSB) that occurs between the column's stainless-steel hardware and certain mycotoxins. The implementation of the inert column offers a robust and improved chromatographic performance as it mitigates the NSB for highly adsorptive analytes (e.g. Fumonisin, Aflatoxins, and Tenuazonic acid) leading to better sensitivity and peak shapes without the need of mobile phase additives or sample passivation.

Methods

Table 1: Analytical Conditions (Waters Xevo TQ-S with Acquity UPLC)

Analytical Column	Raptor Biphenyl 2.7µm 100 mm x 2.1 mm id or Raptor Inert Biphenyl 2.7µm 100 mm x 2.1 mm id	
Guard Column	Raptor Biphenyl EXP Guard Column Cartridge 2.7µm, 5 mm x 2.1 mm id	
Mobile Phase A	0.05% formic acid in water	
Mobile Phase B	0.05% formic acid in methanol	
Gradient	Time (min)	%B
	0.00	25
	5.00	50
	9.00	100
	9.01	25
	11.00	25
Flow Rate	0.4 mL/min	
Injection Volume	5 µL	
Column Temp.	60°C	
Ion Mode	Scheduled MRM in positive ESI	

Food Products

Baby wheat cereal, raw peanut, tomato puree, and flours were purchased from local grocery stores. Baby wheat cereal and tomato puree were used as their original forms. Raw peanut was grinded and stored in the refrigerator. A blended flour was prepared by mixing white rice flour (75%), brown rice flour (5%), millet flour (5%), oat flour (5%), all-purpose wheat flour (5%), and all-purpose gluten free flour (5%) with a handheld blender.

PATENTS & TRADEMARKS

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

Sample and Matrix-Matched Standards Preparation

Two grams of the sample were weighed into a 50-mL polypropylene centrifuge tube and fortified at 5, 50, and 200 µg/kg with stock standard solution. After sitting at room temperature for 10 minutes, 16 mL of extraction solution (80:20 acetonitrile:water) containing 0.5% formic acid (no formic acid for tomato puree) were added and the tube was stirred to gain homogenous suspension. The extraction was carried out by shaking horizontally on a digital pulse mixer (Glas-Col LLC, Terre Haute, IN) at 800 rpm for 20 minutes. After centrifuging for 5 minutes at 4000 rpm, 1 mL of extract was evaporated to dryness at 45°C under a gentle stream of nitrogen. The dried extract was reconstituted with 1 mL of 50:50 water:methanol solution and a 0.4 mL aliquot was transferred to and filtered using a Thomson SINGLE StEP filter vial with a 0.2 µm PTFE filter. To prepare matrix-matched calibration standards, the non-fortified matrices were extracted and dried down as described for the sample preparation procedure followed by reconstitution in 50:50 water/methanol solution containing 0.05 – 50 ng/mL of analytes which equals to 0.4 – 400 µg/kg of sample concentration.

Results & Discussion

- (1) **Chromatographic Performance:** A fast chromatographic method using the Raptor Biphenyl column was established (see **Table 1**) for simultaneous analysis of 38 mycotoxins with a 11-minute total cycling time (**Figure 1**). Analytes were detected with ESI+ and the MRMms were shown in **Table 2**. All epimer pairs of ergot alkaloids were chromatographically separated for definitive and accurate quantification. It was noted that whenever a new Biphenyl column was used, it would need to be rinsed and maintained under the mobile phase overnight to gain an acceptable and quantifiable peak shape for tenuazonic acid.
- (2) **Linearity:** It was shown that a consistent and most suitable linearity of all analytes could be obtained with a quadratic regression ($1/x$ weighted). The lowest concentrated standards were varied due to the differential MS ionization of analytes and specific matrix effect of different food matrices. Nevertheless, most analytes were quantifiable at the full range of 0.4 – 400 $\mu\text{g/kg}$ and all compounds showed proper linearity with $r^2 > 0.997$ and deviations $< 30\%$ (**Table 3**).
- (3) **Accuracy & Precision:** For each food sample, 3 batches of analyses were performed on different days with a total of 9 repetition of each fortified level. The average recovery and relative standard deviation (RSD) were shown in **Table 4**. Except citrinin in solid samples, all analytes had the recovery of 72 – 112% for 3 fortification levels among 4 different types of food matrices. The satisfactory method precision was demonstrated with the %RSD of within 0.5 – 12%. For solid samples, the use of formic acid-containing extraction solution was necessary to obtain adequate recovery for fumonisin Bs but resulted in low recovery (24 – 36%) of citrinin. For food with high water content such as tomato puree, acceptable recovery of both fumonisin Bs (90 – 94%) and citrinin (72 – 77%) were achievable without the addition of formic acid. Due to specific matrix interference, nivalenol could not be measured in baby wheat cereal. The negative impact of matrix interference could also be observed for deoxynivalenol, fusarenon X, and patulin for tomato puree analysis in which the 5 $\mu\text{g/kg}$ fortification sample was not quantifiable.

Figure 1. Chromatogram of Fortified Blended Flour at 50 µg/kg




Table 2. MS Transition and Retention Time

Compounds	Retention time		Precursor ion (m/z)	Product 1 (m/z)	Product 2 (m/z)
	(min)	(min)			
Atorvastatin B1	8.20	31.22	241	284.3	284.3
Atorvastatin B2	7.81	31.51	104(H ⁺)	287.0	259.0
Atorvastatin G1	7.62	32.91	104(H ⁺)	199.7	243.0
Atorvastatin G2	8.15	31.12	241	289.0	289.0
Sichuan A1	8.31	40.41	104(H ⁺)	239.0	358.0
Atorvastatin/deslorvastatin	3.21	339.2	104(H ⁺)	213.1	231.1
10-Acetyl-deslorvastatin	3.14	339.2	104(H ⁺)	137.1	232.2
Deslorvastatin	3.10	297.2	104(H ⁺)	231.0	249.0
Deslorvastatin/pirfenidone	5.73	384.2	384.2	247.1	307.2
Fluoxetine B1	5.63	72.25	104(H ⁺)	152.3	314.2
Fluoxetine B2	5.63	70.64	104(H ⁺)	136.8	318.1
Fluoxetine B3	6.32	70.64	104(H ⁺)	136.2	318.1
Fluoxetine X1	1.98	355.1	104(H ⁺)	237.1	247.1
HT-2	6.20	447.2	104(H ⁺)	345.1	285.1
Nivalenol	0.92	295.1	[M-H] ⁺ O	137.1	91.0
T2	7.14	489.2	104(H ⁺)	387.1	245.1
2-Zealandin	6.96	305.1	[M-H] ⁺ O	285.1	175.0
Zealandin	7.56	319.2	104(H ⁺)	283.1	187.0
Citrinin	5.43	251.2	104(H ⁺)	233.1	205.0
Patulin	1.03	155.0	104(H ⁺)	99.1	81.0
Patulin B1	3.20	259.0	104(H ⁺)	157.1	187.0
Alerianol					
Alsterlundether	7.69	273.0	173.0	199.1	108.0
Tertenn	4.70	293.2	104(H ⁺)	257.1	275.2
Atenolol	6.70	415.2	104(H ⁺)	312.2	306.1
Atenolol B1	6.70	398.2	104(H ⁺)	159.2	155.0
Ergocristine	6.03	562.0	562.0	268.2	223.2
Ergocristine	7.07	562.0	562.0	268.2	233.2
Ergocristine	6.13	610.0	610.0	232.2	592.2
Ergocristine	7.53	576.0	576.0	232.2	592.2
Ergocryptine	6.32	576.0	576.0	268.2	232.2
Ergocryptine	6.13	574.0	574.0	232.2	592.2
Ergometrine	1.27	326.2	104(H ⁺)	232.2	208.1
Ergometrine	1.83	326.2	104(H ⁺)	232.2	208.1
Ergometrine	5.47	548.0	548.0	232.2	208.1
Ergometrine	5.67	548.0	548.0	232.2	208.1
Ergometrine	5.90	548.0	548.0	232.2	208.2
Ergometrine	5.67	548.0	548.0	232.2	208.2

Table 4. Recovery & Precision


Concentration, µg/kg	Baby White Cereal				Peanut				Average Recovery (RSD, %)				Tonnato Pure				Blended Flour			
	5	50	100	200	5	50	100	200	5	50	100	200	5	50	100	200	5	50	100	200
Aflatoxin B1	105 (4.8)	100 (0.0)	79.8 (2.6)	82 (1.6)	98 (2.6)	99 (0.2)	89 (0.5)	87 (2.7)	92 (7.8)	97 (6.2)	103 (3.0)	101 (2.8)	101 (2.8)	95.5 (1.3)	89 (0.15)		101 (2.8)	95.5 (1.3)	89 (0.15)	
Aflatoxin B2	110 (4.1)	109 (2.8)	106 (2.3)	102 (5.8)	105 (8.1)	97 (3.4)	91 (2.9)	91 (2.9)	91 (2.4)	93 (3.9)	94 (7.0)	100 (2.3)	100 (2.3)	101 (0.9)	88 (7.1)		100 (2.3)	101 (0.9)	88 (7.1)	
Aflatoxin G1	105 (6.1)	107 (1.7)	102 (2.1)	98 (2.4)	97 (3.8)	91 (2.4)	91 (2.4)	91 (2.4)	91 (2.1)	92 (2.6)	93 (3.5)	99 (3.1)	92 (2.6)	96 (1.4)	93 (2.5)		99 (3.1)	100 (1.6)	93 (2.5)	
Aflatoxin G2	108 (3.0)	109 (1.3)	104 (2.2)	104 (1.2)	100 (1.8)	95 (1.2)	91 (2.1)	91 (2.1)	86 (8.8)	94 (2.5)	98 (2.3)	96 (7.1)	101 (2.6)	100 (1.4)	94 (2.0)		96 (7.1)	100 (1.4)	94 (2.0)	
Ochratoxin A	109 (1.8)	105 (2.1)	94.5 (1.0)	101 (1.1)	97 (2.0)	91 (1.1)	91 (1.1)	91 (1.1)	92 (2.1)	92 (1.8)	92 (1.3)	92 (1.3)	92 (1.3)	92 (1.3)	92 (1.3)		92 (1.3)	92 (1.3)	92 (1.3)	
3- < 15-Acetyldeoxynivalenol	104 (3.6)	108 (1.8)	104 (3.3)	101 (1.5)	95 (5.8)	91 (0.4)	91 (0.4)	91 (0.4)	91 (9.4)	98 (12.7)	95 (0.18)	98 (4.2)	101 (2.9)	101 (2.9)	100 (0.9)		98 (4.2)	101 (2.9)	100 (0.9)	
Deoxynivalenol	112 (4.0)	102 (1.6)	95 (7.1)	98 (1.5)	93 (7.48)	88 (2.3)				90 (3.6)	94 (5.2)	102 (3.5)	97 (5.2)	97 (5.2)	96 (1.0)		102 (3.5)	97 (5.2)	96 (1.0)	
Cytosineguanylinol	105 (14.0)	107 (1.5)	103 (1.2)	93 (2.43)	95 (4.39)	93 (8.50)				90 (9.38)	94 (5.47)	94 (0.19)	98 (1.63)	101 (3.1)	98 (7.18)		98 (1.63)	101 (3.1)	98 (7.18)	
Fumonisin B1	94 (3.64)	94 (0.28)	92 (3.26)	87 (2.31)	88 (2.45)	87 (8.66)				91 (8.36)	91 (5.19)	91 (9.70)	100 (2.3)	99 (6.17)	96 (1.21)		100 (2.3)	99 (6.17)	96 (1.21)	
Fumonisin B2	93 (3.41)	93 (0.41)	93 (1.29)	94 (7.32)	94 (7.32)	88 (8.19)				89 (8.41)	92 (3.23)	92 (8.08)	104 (2.7)	99 (6.14)	94 (1.21)		104 (2.7)	99 (6.14)	94 (1.21)	
Fumonisin B3	91 (8.49)	94 (4.49)	91 (2.48)	96 (0.27)	90 (1.28)	91 (8.49)				91 (8.49)	90 (1.28)	90 (1.28)	104 (2.7)	99 (6.14)	95 (1.21)		104 (2.7)	99 (6.14)	95 (1.21)	
Fusarenen X	99 (9.39)	100 (2.9)	91 (2.48)	86 (9.70)	93 (1.10)	88 (3.10)				92 (0.68)	93 (9.39)	93 (9.39)	100 (2.3)	99 (6.17)	98 (1.6)		100 (2.3)	99 (6.17)	98 (1.6)	
HT-2	110 (2.4)	111 (1.0)	108 (1.1)	100 (1.7)	100 (2.0)	94 (3.30)				96 (8.31)	95 (1.21)	99 (0.14)	101 (1.6)	101 (2.3)	98 (1.3)		101 (1.6)	101 (2.3)	98 (1.3)	
Nivalenol					98 (3.62)	89 (0.16)				92 (5.45)	93 (7.15)						95 (5.47)	92 (2.3)	92 (2.3)	
T-2	111 (2.1)	110 (1.8)	108 (2.4)	99 (1.27)	91 (1.7)	95 (9.21)				92 (0.63)	94 (7.13)	98 (6.13)	102 (1.3)	101 (1.3)	96 (1.6)		102 (1.3)	101 (1.3)	96 (1.6)	
Zearalenol	100 (4.9)	102 (2.2)	98 (1.54)	89 (8.16)	97 (4.73)	94 (7.93)				94 (8.95)	96 (6.17)	96 (6.17)	100 (2.3)	99 (6.17)	94 (1.21)		100 (2.3)	99 (6.17)	94 (1.21)	
Zearalenone	110 (6.7)	105 (1.0)	105 (3.7)	98 (3.73)	97 (4.78)	93 (3.15)				95 (0.45)	93 (6.27)	95 (7.20)	101 (2.3)	101 (2.3)	92 (1.4)		101 (2.3)	101 (2.3)	92 (1.4)	
Citralin	26 (1.2)	26 (6.1)	301 (3.8)	214 (1.87)	25 (1.19)	25 (8.35)				71 (9.47)	76 (4.16)	77 (1.77)	32 (3.95)	32 (2.63)	35 (8.45)					
Patulin	166 (4.6)	95 (5.6)	89 (2.51)	88 (8.120)	83 (6.09)	80 (7.72)				98 (9.36)	101 (4.5)	93 (6.44)	86 (1.3)	92 (2.29)						
Alternariol	108 (4.1)	108 (1.6)	104 (1.0)	94 (2.14)	95 (4.24)	96 (2.72)				89 (3.46)	91 (8.25)	91 (4.13)	98 (4.23)	101 (2.5)	96 (3.32)					
Alternariol monomethyl ether	108 (4.1)	109 (2.2)	99 (3.12)	93 (5.13)	93 (5.87)	89 (8.24)				91 (3.66)	88 (7.51)	93 (9.39)	104 (2.9)	101 (1.7)	93 (7.19)					
Patulinol	101 (2.1)	109 (1.2)	105 (2.1)	95 (1.2)	95 (5.12)	92 (2.72)				88 (4.34)	92 (4.16)	92 (8.18)	100 (2.3)	99 (6.17)	98 (1.6)					
Tenutoxin	111 (3.6)	109 (2.5)	103 (1.4)	104 (1.29)	101 (1.1)	95 (3.14)				92 (5.62)	94 (2.72)	95 (8.14)	104 (2.4)	105 (2.1)	98 (2.19)					
Tricenoic acid		85 (8.17)	87 (4.63)	92 (5.47)	91 (0.21)	88 (5.24)				89 (3.43)	88 (5.20)						92 (5.88)	90 (0.95)		
Ergocornine	109 (1.5)	109 (1.4)	102 (1.3)	93 (8.15)	93 (2.44)	91 (2.33)				91 (5.30)	93 (1.19)	92 (9.06)	102 (2.5)	101 (1.9)	97 (6.17)					
Ergocornine	109 (1.5)	109 (1.2)	101 (1.9)	105 (3.0)	104 (2.4)	95 (5.81)				89 (9.38)	93 (2.12)	92 (5.31)	101 (2.5)	102 (2.6)	95 (7.24)					
Ergosinine	108 (2.3)	108 (2.9)	103 (4.4)	92 (1.88)	91 (7.51)	92 (0.22)				93 (1.29)	94 (2.00)	94 (3.08)	105 (1.7)	99 (8.20)	96 (7.18)					
Ergosinine	106 (5.3)	105 (1.4)	104 (1.0)	100 (6.63)	94 (4.41)	90 (6.16)				91 (6.44)	92 (4.41)	92 (6.16)	100 (2.3)	99 (6.17)	98 (1.6)					
Ergosyringine	107 (2.0)	109 (1.9)	104 (3.4)	90 (2.00)	94 (7.61)	92 (1.77)				93 (1.30)	93 (5.22)	93 (2.07)	95 (2.79)	99 (9.12)	97 (4.14)					
Ergosyringine	106 (1.7)	108 (2.0)	101 (1.1)	103 (5.3)	105 (4.40)	101 (4.2)				91 (1.43)	95 (1.15)	98 (1.16)	101 (2.0)	101 (1.8)	95 (4.19)					
Ergotrine	92 (8.73)	90 (0.42)	88 (3.36)	101 (2.3)	96 (2.16)	86 (7.19)				90 (7.36)	88 (9.61)	87 (6.15)	100 (1.8)	99 (7.32)	95 (3.13)					
Ergotristine	101 (1.2)	99 (1.19)	94 (3.07)	93 (2.43)	95 (5.17)	89 (3.12)				90 (2.06)	90 (1.44)	89 (1.13)	100 (1.5)	98 (5.19)	91 (1.19)					
Ergine	106 (2.6)	106 (5.6)	103 (3.2)	90 (8.09)	93 (8.22)	89 (2.24)				92 (1.32)	90 (4.31)	90 (7.19)	99 (9.27)	99 (1.80)	98 (1.1)					
Ergine	111 (1.8)	109 (0.9)	100 (1.2)	100 (1.2)	97 (1.63)	97 (1.63)				92 (1.63)	92 (1.63)	92 (1.63)	92 (1.63)	92 (1.63)	92 (1.63)					
Ergotamine	109 (1.9)	108 (1.7)	102 (2.8)	91 (0.28)	92 (6.18)	89 (8.16)				91 (1.21)	90 (6.17)	90 (7.13)	101 (1.9)	100 (1.3)	96 (4.22)					
Ergotamine	109 (1.9)	109 (0.7)	101 (0.8)	88 (2.20)	101 (1.15)	96 (3.13)				93 (6.15)	94 (7.17)	94 (5.06)	101 (2.3)	99 (7.13)	97 (1.1)					

Table 3. Calibration Ranges

Compound	Baby Wheat Cereal			Peanut			Tomato Puree			Blended Food		
	Unlabeled (μg/L)	<i>t</i> ¹	<i>p</i> ²	Unlabeled (μg/L)	<i>t</i> ¹	<i>p</i> ²	Unlabeled (μg/L)	<i>t</i> ¹	<i>p</i> ²	Unlabeled (μg/L)	<i>t</i> ¹	<i>p</i> ²
Aflatoxin B1	0.4	400	0.9996	0.4	400	0.9998	0.4	400	0.9995	0.4	400	0.9998
Aflatoxin B2	0.4	400	0.9997	0.4	400	0.9998	0.4	400	0.9996	0.4	400	1.000
Aflatoxin G1	0.4	400	0.9999	0.4	400	0.9997	0.4	400	0.9979	0.4	400	1.000
Aflatoxin G2	0.4	400	0.9997	0.4	400	0.9998	0.4	400	0.9996	0.4	400	0.9998
Ochratoxin	0.4	400	0.9998	0.4	400	0.9993	0.4	400	0.9993	0.4	400	1.000
3-Acetyldeoxymiscelol	0.4	400	0.9994	2.0	400	0.9997	0.4	400	0.9982	2.0	400	0.9998
Deoxynivalenol	0.4	400	0.9998	0.4	400	0.9994	0.4	400	0.9991	2.0	400	1.000
Diacetoxysynepicol	0.8	400	0.9998	0.4	400	0.9995	0.8	400	0.9993	0.4	400	0.9998
Fumonisin B1	0.4	400	0.9999	0.4	400	0.9994	0.4	400	0.9999	0.4	400	0.9999
Fumonisin B2	0.4	400	0.9997	0.4	400	0.9997	0.4	400	0.9998	0.4	400	0.9999
Fumonisin B3	0.4	400	0.9999	0.4	400	0.9997	0.4	400	1.000	0.4	400	0.9999
Fusarenon X	0.4	400	0.9971	2.0	400	0.9971	8.0	400	0.9974	2.0	400	0.9995
HT-2	0.4	400	0.9999	0.4	400	0.9997	0.4	400	0.9997	0.4	400	0.9999
Neovialenol	0.4	400	0.9999	8.0	400	0.9990	0.4	400	0.9997	0.4	400	0.9997
T-2	0.4	400	0.9998	0.4	400	0.9998	0.4	400	0.9992	0.4	400	1.000
α-Zearalenol	0.4	400	0.9985	2.0	400	0.9992	2.0	400	0.9973	2.0	400	0.9994
Zearalenone	0.4	400	0.9998	0.4	400	0.9996	0.8	400	0.9999	2.0	400	0.9998
Citralin	0.4	400	0.9996	0.4	400	0.9986	0.4	400	0.9984	0.4	400	0.9999
Patulin	0.4	400	0.9991	0.4	400	0.9995	8.0	400	0.9997	0.4	400	0.9993
Alternariol	0.4	400	0.9998	0.4	400	0.9990	0.4	400	0.9996	0.4	400	0.9997
Alternariol monomethyl ether	0.4	400	0.9996	0.4	400	0.9996	0.4	400	0.9993	0.4	400	0.9999
Altrenurone	0.4	400	0.9999	0.4	400	0.9997	2.0	400	0.9999	0.4	400	1.000
Tenoxatin	0.2	400	0.9998	0.4	400	0.9998	0.8	400	0.9998	0.4	400	0.9999
Trichothecene	0.4	400	0.9999	0.4	400	0.9997	8.0	400	0.9987	0.4	400	0.9992
Ergosterione	0.4	400	0.9999	0.4	400	0.9998	0.4	400	0.9998	0.4	400	0.9999
Ergocornine	0.4	400	0.9998	0.4	400	0.9997	0.4	400	0.9996	0.4	400	0.9999
Ergosinine	0.4	400	0.9998	0.4	400	0.9997	0.4	400	0.9997	0.4	400	1.000
Ergocristine	0.4	400	0.9999	0.4	400	0.9997	0.4	400	0.9999	0.4	400	0.9999
Ergocryptine	0.4	400	0.9999	0.4	400	0.9998	0.4	400	0.9998	0.4	400	0.9999
Ergocryptine	0.4	400	0.9999	0.4	400	0.9998	0.4	400	0.9988	0.4	400	1.000
Ergonovine	0.4	400	0.9998	0.4	400	0.9999	0.4	400	0.9973	0.4	400	0.9999


Further Findings

Non-Specific Binding (NSB) was leading to unstable signals, peak tailing, and low sensitivity as the injections progressed for most of the mycotoxins; including Fumonisin (Figure 2), Tenuazonic acid, Citrinin, Ochratoxin A, and Diacetoxycispermol. NSB was attributed to the presence of active metal surfaces within the column's hardware and instrument.




NSB is aggravated under low pH and low ionic strength mobile phases, leading to irreproducible chromatography. Analytes with acidic functional groups or chelating moieties are particularly prone to this surface metal interaction.

Figure 2. Chromatograms of Fumonisin B1 at 10 ng/mL



In order to mitigate NSB, an inert column hardware was implemented; which allows for the efficient and reproducible analysis of metal sensitive compounds. The inert hardware consists of a premium inert coating applied to the stainless-steel surface of the column that guarantees a more consistent chromatography. **Figure 3** shows the overall benefit of switching from a conventional to an inert column without any change in the method itself.

Figure 3. Effect of Inert Column on Selected Mycotoxins at 10 ng/mL



Conclusions

A workflow was established in this study to provide a unique solution for simultaneous determination of *Alternaria* toxins, ergot alkaloid epimers, and other major mycotoxins produced by fungal genus of *Aspergillus*, *Fusarium*, and *Penicillium*. The reported method was rugged, accurate, and precise using a combination of convenient sample preparation procedure and a fast 11-minute chromatographic analysis. Most importantly, this solution could be applied to multi-mycotoxin quantification in a wide variety of food products. Furthermore, the application of inert column hardware aided in the consistent analysis of several mycotoxins that tend to interact with metal surfaces and therefore increasing the sensitivity of such.