

New Alumina Column Shows Promise for Analysing Chlorofluorocarbons

Jaap de Zeeuw* and Tom Vezza

Restek Corporation, 110 Benner Circle, Bellefonte, PA, US

* Restek Corporation, Weerhaan 9, Middelburg, The Netherlands

A unique deactivation for alumina PLOT columns has recently been developed to control the high activity of alumina surfaces. Rt®-Alumina BOND/CFC columns utilize this deactivation and provide high retention and high selectivity for volatile halogenated hydrocarbons, without the activity that is usually observed with alumina adsorbents. Chlorofluorocarbon analysis, in particular, can benefit from this new technology.

Chlorofluorocarbons (CFCs) are a group of halogenated hydrocarbons that were heavily used as refrigerants and propellants, until their role in the depletion of the earth's ozone layer became known. In addition to being active ozone depleting substances, CFCs are also greenhouse gases that contribute to global warming. Now, the manufacture of CFCs is prohibited in most countries and the use of recycled material is restricted. Because of this, CFC analysis has become a priority, particularly in the environmental field and for industry as companies develop alternative compounds.

CFCs are highly volatile and amenable to gas chromatography, but an extremely retentive GC column is needed for successful analysis. Thick film nonpolar capillaries can be used, but these are low in efficiency and generally do not perform optimally. Separations of CFCs have also been done with packed columns, as they can be made with adsorbents or liquid phases with a low phase ratio which allows the resolution of volatile compounds with very specific selectivity. Separations on packed columns are mainly possible due to selectivity combined with high *k* values; efficiency is still relatively poor as the number of theoretical plates is very low.

The introduction of chemically bonded phases and porous layer open tubular (PLOT) capillary columns allowed selectivity to be combined with a high number of theoretical plates, which provided more opportunity for resolution. Alumina, in particular, offers high retention and good selectivity for a large range of CFCs, but historically these columns have been limited by high activity. Restek recently developed the Rt®-Alumina BOND/CFC column, a new alumina PLOT column which uses a unique deactivation that is quite promising for CFC analysis.

Alumina: Limited by Reactivity

In order to benefit from the retention and selectivity of alumina columns for volatiles analysis, the reactivity of the alumina must be reduced. This is usually accomplished through salt deactivations, often using sodium sulfate or potassium chloride. Even with these conventional deactivations, reactivity is often a problem. For example, using 1,2-butadiene as a test probe, column activity is revealed by the formation of reaction products (Figures 1 and 2). In contrast, the new CFC deactivation developed by Restek results in a more fully passivated surface (Figure 3). The deactivation is also quite robust, as the column can be held at the maximum temperature for alumina (200°C) for 16 hours and no change in the retention or response of volatile halogenated hydrocarbons is seen (Figure 4).




Figure 1: Butadiene on conventional alumina PLOT column with sodium sulfate deactivation. Butadiene, used here as an indicator of inertness, decomposes rapidly on a conventional alumina column because the sodium sulfate deactivation does not effectively prevent reactivity. Column: 50 m x 0.32 mm; Sample: 1,2-butadiene, 40 µL; Oven: 150°C; Carrier gas: H₂, 156 kPa, 23 psi; Injection: split (1:30).




Figure 2: Butadiene on conventional alumina PLOT column with potassium chloride deactivation. Inertness is somewhat improved, but the column is still too reactive as evidenced by the formation of reaction products. See Figure 1 for conditions.




Figure 3: Butadiene on Rt®-Alumina BOND/CFC columns. These new alumina PLOT columns are significantly more inert than other alumina columns, as shown by the reduced reactivity with butadiene. See Figure 1 for conditions.




Figure 4: Impact of high temperature activation of alumina. Highly inert Rt®-Alumina BOND/CFC columns produce stable responses and retention times, even after 16 hours at maximum temperature. Column: Rt®-Alumina BOND/CFC, 30 m x 0.53 mm; Sample: halogenated hydrocarbons; Oven: 135°C; Carrier gas: H₂, 34 kPa, 5 psi; Injection: split (20:1); Peaks: 1. Dichlorofluoromethane, 2. Chloromethane, 3. Vinyl chloride, 4. Trichlorofluoromethane, 5. Chloroethane.

A Promising New Alumina for CFCs

The Rt®-Alumina BOND/CFC column provides excellent retention and selectivity for CFC analysis. Figures 5 and 6 show practical applications of this new technology for CFC impurities analysis. Different components are fully resolved and show no signs of decomposition. Although this new column is highly inert, residual reactivity for some mono- or di-substituted halogenated hydrocarbons remains, however the majority of these compounds can be accurately quantified from mainstream processes or in impurity analyses.

Summary

Overall, early results from our lab, as well as feedback from users in the field, show the new Rt®-Alumina BOND/CFC column to be quite promising for CFC separations. The unique deactivation effectively reduces column reactivity and allows the benefits of alumina retention and selectivity to be realised for laboratories testing these priority compounds.

Restek Corporation

Restek is a leading developer and manufacturer of chromatography columns, standards, and accessories. We provide analysts around the world with the innovative tools they need to monitor the quality of air, water, soil, food, pharmaceutical, chemical, and petroleum products.

www.restek.com

RESTEK
A Company of Owners




Figure 5: Separation of chlorofluorocarbons on an Rt®-Alumina BOND/CFC column. Column: Rt®-Alumina BOND/CFC, 30 m x 0.53 mm; Sample: CFC-134a; Oven: 80°C (hold 6 min.) to 140°C at 10°C/min. (hold 2 min.); Injection: gas sampling, 500 µL; Carrier gas: helium; Peaks: 1. Chloropentafluoroethane (CFC-115), 2. Dichlorodifluoromethane (CFC-12), 3. Chlorodifluoromethane (CFC-22), 4. 1,1,1,2-Tetrafluoroethane (CFC-134a). Chromatogram courtesy of André Hähnel, Westfalen AG.




Figure 6: Impurity analysis of CFC-134a on an Rt®-Alumina BOND/CFC column. Symmetric peaks with no sign of decomposition are obtained due to high column inertness. See Figure 5 for conditions. Peaks: 1. Difluoromethane (CFC-32), 2. Pentafluoroethane (CFC-125), 3. 1,1,1,2-Tetrafluoroethane (CFC-134a). Chromatogram courtesy of André Hähnel, Westfalen AG.