

Highly Efficient LC-MS/MS Analysis of Organophosphorous Pesticides Utilizing ARC-18 Column Selectivity with Inert Column Technology

Tina Brandscher¹, Hansjörg Majer¹, Paul Connolly² and Diego Lopez²; ¹ Restek GmbH, Bad Homburg, DE; ² Restek, Bellefonte, PA, USA

Introduction

Pesticides are the natural or synthetic substances created to eliminate all types of pests. They are used ubiquitously to help increase crop yields; however, they can pose risks for public health and pollinators (honeybees) and are therefore strictly monitored.

However, monitoring this diverse, ever-growing group of chemicals is not always easy. Pesticides can often be modified by just a few atoms, leading to significant changes in their properties. To accurately quantify the resulting isomers and isobars, chromatographic separation is essential, which requires the careful selection of stationary phases, such as the Raptor ARC-18 demonstrated here.

Moreover, subgroups like e.g. organophosphorous pesticides are particularly difficult, as they tend to form nonspecific adsorption (NSA) or nonspecific binding (NSB) to metal surfaces inside analytical columns resulting in poor peak shape, low response and sensitivity.

In this work, Restek looks to establish the benefits of coated column technologies by comparing methods developed on Inert and standard hardware to increase sensitivity and peak shape on the one hand and decrease time and effort for matrix or chemical based passivation techniques on the other for a wide panel of organophosphorus pesticides.

Materials and Method

Table I: Method developed by Restek to analyze Pesticides with inert hardware for better peak shape and response.

| Columns | Raptor (Inert) ARC-18 | | |
|------------------|---|----|-----|
| Dimensions: | 100 mm x 2.1 mm ID | | |
| Particle Size: | 2.7 µm | | |
| Temp.: | 50 °C | | |
| Standard/Sample: | LC multiresidue pesticide standard #1 (cat. # 31972) | | |
| Diluent: | Water, 0.1% Formic Acid | | |
| Conc.: | 1 ng/ mL | | |
| Inj. Vol.: | 5 µL | | |
| Detector: | MS/MS | | |
| Ion Mode: | ESI ⁺ | | |
| | MRM | | |
| Mobile Phase A: | Water, 2 mM NH ₄ CH ₃ COO 0.1% formic acid | | |
| Mobile Phase B: | MeOH, 2 mM NH ₄ CH ₃ COO, 0.1% formic acid | | |
| Time (min) | Flow (mL/min) | %A | %B |
| 0.00 | 0.4 | 95 | 5 |
| 2.00 | 0.4 | 40 | 60 |
| 4.00 | 0.4 | 25 | 75 |
| 6.00 | 0.4 | 0 | 100 |
| 7.50 | 0.4 | 0 | 100 |
| 7.51 | 0.4 | 95 | 5 |
| 9.00 | 0.4 | 95 | 5 |

Results and Discussion

In the data table below, we analyse the area and height ratio between the analytes on each of the two hardware types. Values greater than 1 describing a positive improvement in response, and less than 1 describing a loss of sensitivity.

Across the full panel of 15 organophosphorous pesticides analyzed, **almost all analytes** showed an **improvement** in sensitivity for both peak area and peak height. Of the panel, **Carbaryl** was the **sole** compound which showed a **slight reduction** in both metrics, while Mevinphos Isomer 2 showed comparable results. All other components showed ratios up to **2.0**, meaning that **signal intensity is doubled** using inert columns.

The compounds Methamidophos, Acephate, Omethoate, Monocrotophos, Vamidothion and Temephos are typically more problematic to analyse in matrix samples, with a **tendency for NSA/NSB** interactions leading to poor sample recoveries and sensitivity issues.

The use of inert hardware here resulted in an **average of 1.6-fold improvement** in peak area, and **1.7-fold improvement** for peak height across the six compounds, along with significant improvement on signal/noise ratio and a reduction of tailing. This demonstrates the improvement in NSA/NSB interactions when switching to Inert columns.

Table II: Analysis of 15 organophosphorous pesticides with stainless steel (SS) and Inert Raptor ARC-18 column. Comparison of area and height ratios. Values > 1 = better sensitivity, Values < 1 = worse selectivity (only Carbaryl, marked in orange).

| Compound | Peak Area | | | Peak Height | | |
|-----------------------|-----------------|---------|------------------------|-----------------|--------|-------------------------|
| | Stainless Steel | Inert | Areas Ratio (Inert/SS) | Stainless Steel | Inert | Height Ratio (Inert/SS) |
| Methamidophos | 254969 | 428941 | 1.68 | 52553 | 105189 | 2.00 |
| Acephate | 168776 | 300642 | 1.78 | 58418 | 104729 | 1.79 |
| Omethoate | 579502 | 892008 | 1.54 | 216157 | 337690 | 1.56 |
| Monocrotophos | 140095 | 215810 | 1.54 | 51402 | 78425 | 1.53 |
| Dicrotophos | 340978 | 404916 | 1.19 | 135380 | 159292 | 1.18 |
| Dimethoate | 461156 | 807805 | 1.75 | 188746 | 342939 | 1.82 |
| Trichlorfon | 84233 | 173942 | 2.07 | 34793 | 63266 | 1.82 |
| Vamidothion | 913264 | 1333829 | 1.46 | 354311 | 547308 | 1.54 |
| Mevinphos Isomer 1 | 213632 | 311274 | 1.46 | 82105 | 129961 | 1.58 |
| Mevinphos Isomer 2 | 56093 | 74030 | 1.32 | 29070 | 29802 | 1.03 |
| Carbaryl | 43590 | 39671 | 0.91 | 14563 | 11924 | 0.82 |
| Isocarbophos | 21587 | 33294 | 1.54 | 9062 | 11941 | 1.32 |
| Dimethomorph Isomer 1 | 462425 | 511766 | 1.11 | 166990 | 172977 | 1.04 |
| Dimethomorph Isomer 2 | 896109 | 877031 | 0.98 | 311657 | 328826 | 1.06 |
| Temephos | 98793 | 164310 | 1.66 | 35383 | 64751 | 1.83 |




Figure 1: Comparison of different organophosphorous pesticides analyzed on Raptor Inert ARC-18 and traditional stainless steel Raptor ARC-18. Right: Chemical Structures of pesticides

Conclusion

Choosing the Raptor ARC-18 stationary phase resulted in a balanced retention profile over the whole gradient. Isobars like the Mevinphos or Dimetomorph isomers could be separated with high target resolution.

In combination with the inert column hardware technology enabled low detection limits, improved peak shapes and high response rates. Laborious conditioning and complicated passivation are not necessary with these coated inert columns, saving time and enables to measure more samples/time.

VISIT US AT BOOTH NO. 26

For more information contact Jan Pschierer; Jan.Pschierer@restek.com - LC Specialist • Restek GmbH, Schaberweg 23, 61348 Bad Homburg v. d. H.