

Automated In-Line Pigment Removal for Multi-Residue Pesticide Analysis in Spinach by LC-MS/MS

Sharon Lupo, Randy Romesberg, Xiaoning Lu

Restek Corporation

Introduction

- The presence of pigments in sample extracts can compromise the LC-MS/MS analysis of pesticides in green vegetables through ionization suppression or enhancement.
- A faster and simpler workflow has been developed which effectively removes pigments from sample extracts without the loss of pesticides inherent to dispersive SPE (dSPE) which contain graphitized carbon black (GCB).

Methods

Table 1: LC instrument conditions for ILSP workflow



Instrument:	Shimadzu Nexera UHPLC and 8045 MS/MS
ILSP-PR Cartridge:	5 x 2.1 mm
LC Column:	Raptor ARC-18 2.7 μ m, 100 x 2.1 mm (Cat# 9314A12)
Pump A:	Water, 2mM NH ₄ HCO ₂ + 0.2% HCOOH
Pump B/Aux.:	Methanol, 2mM NH ₄ HCO ₂ + 0.2% HCOOH
Flow A/B, Aux.:	0.5 mL/min, 1.0 mL/min
Time Program:	
Oven Temp.:	50 °C
Injection Volume:	4 μ L of spinach extracts fortified with LC Multiresidue Pesticide Standard #4 (Cat# 31975)

Figure 1: ILSP valve configuration for: a) load, and b) wash



Results & Discussion

- Matrix constituents were removed using an UHPLC system equipped with an auxiliary pump, 6-port switching valve, and dual-directional ILSP cartridge (Table 1 and Figure 1).
- The ILSP procedure was compared to an AOAC QuEChERS extraction (Figure 2) and evaluated for recovery (Figure 3) and matrix effects (Figures 4-5) for 63 pesticides in spinach.

Figure 2: Comparison of the number of steps and the time required to prepare 14 samples: ILSP vs. QuEChERS




Figure 3: Analytes with recoveries of 70-120% and RSD <20% for fortified spinach samples extracted by ILSP vs. QuEChERS




Figure 4: Impact of cleanup procedure (dSPE vs. ILSP) on calculated matrix effects for QuEChERS extracted spinach samples: (a) overall distribution, (b) planar pesticides





Figure 5: XIC of pigments remaining in QuEChERS extracted spinach samples cleaned up with either (a) dSPE or (b) ILSP



Conclusion

- The ILSP workflow has fewer steps, takes 1/5 the time, and is more amenable to automation than QuEChERS.
- ILSP removes more pigments than dSPE with GCB, reducing MS contamination without loss of planar pesticides.
- Of the 63 pesticides tested, 86% met the acceptance criteria for samples fortified at 0.005 mg/kg and 100% met the criteria at 0.05 mg/kg.

PATENTS & TRADEMARKS
Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.