

Blood Alcohol and Inhalants of Abuse Analysis Using Gas Chromatography Optimized with Computer Software.

Ramkumar Dhandapani, Corby Hilliard, Jamie York, Jared Burkhart.
Restek Corporation, 110 Benner Circle, Bellefonte, PA 16823

Introduction

Screening for volatile inhalants of abuse, as well as analyzing blood alcohol content, is commonly performed in forensic toxicology laboratories using headspace gas chromatography with flame ionization detection (HS-GC-FID). The analyses are generally performed using dual columns with specialized stationary phases that can actually separate these volatile compounds. While separation profiles of standard blood alcohol screening compounds are usually well established by manufacturers of these specialized specific columns, elution profiles of inhalants may not be as readily available. In addition, providing example chromatograms with static run conditions may not be sufficient, who want to experiment with factors like column dimensions, carrier gases, and flow rates.

These issues can be addressed using computer modeling software to predict retention times of compounds of interest on various stationary phases. In addition to the ability of the web-based software to help select a column and provide an optimized separation of compounds of interest on a specific stationary phase, the software can also be used to make changes to analytical conditions and observe the effect on retention, making it a valuable tool for method development and optimization.

Figure 1: Comparison of instrument confirmation run retention times with computer predicted retention times for BAC and for BAC. There is excellent agreement, with all predicted retention times being ± 0.5 difference from actual.

BAC-1				BAC-2			
Compound	Instrument Retention Time (min)	Model Predicted Retention Time (min)	Difference (min)	Compound	Instrument Retention Time (min)	Model Predicted Retention Time (min)	Difference (min)
1,1,1-trichloroethane	1.22	1.24	-0.02	1,1,1-trichloroethane	1.22	1.24	-0.02
1,2-dichloroethane	1.22	1.24	-0.02	1,2-dichloroethane	1.22	1.24	-0.02
1,3-dichloroethane	1.22	1.24	-0.02	1,3-dichloroethane	1.22	1.24	-0.02
1,4-dichloroethane	1.22	1.24	-0.02	1,4-dichloroethane	1.22	1.24	-0.02
1,5-dichloroethane	1.22	1.24	-0.02	1,5-dichloroethane	1.22	1.24	-0.02
1,6-dichloroethane	1.22	1.24	-0.02	1,6-dichloroethane	1.22	1.24	-0.02
1,7-dichloroethane	1.22	1.24	-0.02	1,7-dichloroethane	1.22	1.24	-0.02
1,8-dichloroethane	1.22	1.24	-0.02	1,8-dichloroethane	1.22	1.24	-0.02
1,9-dichloroethane	1.22	1.24	-0.02	1,9-dichloroethane	1.22	1.24	-0.02
1,10-dichloroethane	1.22	1.24	-0.02	1,10-dichloroethane	1.22	1.24	-0.02
1,11-dichloroethane	1.22	1.24	-0.02	1,11-dichloroethane	1.22	1.24	-0.02
1,12-dichloroethane	1.22	1.24	-0.02	1,12-dichloroethane	1.22	1.24	-0.02
1,13-dichloroethane	1.22	1.24	-0.02	1,13-dichloroethane	1.22	1.24	-0.02
1,14-dichloroethane	1.22	1.24	-0.02	1,14-dichloroethane	1.22	1.24	-0.02
1,15-dichloroethane	1.22	1.24	-0.02	1,15-dichloroethane	1.22	1.24	-0.02
1,16-dichloroethane	1.22	1.24	-0.02	1,16-dichloroethane	1.22	1.24	-0.02
1,17-dichloroethane	1.22	1.24	-0.02	1,17-dichloroethane	1.22	1.24	-0.02
1,18-dichloroethane	1.22	1.24	-0.02	1,18-dichloroethane	1.22	1.24	-0.02
1,19-dichloroethane	1.22	1.24	-0.02	1,19-dichloroethane	1.22	1.24	-0.02
1,20-dichloroethane	1.22	1.24	-0.02	1,20-dichloroethane	1.22	1.24	-0.02
1,21-dichloroethane	1.22	1.24	-0.02	1,21-dichloroethane	1.22	1.24	-0.02
1,22-dichloroethane	1.22	1.24	-0.02	1,22-dichloroethane	1.22	1.24	-0.02
1,23-dichloroethane	1.22	1.24	-0.02	1,23-dichloroethane	1.22	1.24	-0.02
1,24-dichloroethane	1.22	1.24	-0.02	1,24-dichloroethane	1.22	1.24	-0.02
1,25-dichloroethane	1.22	1.24	-0.02	1,25-dichloroethane	1.22	1.24	-0.02
1,26-dichloroethane	1.22	1.24	-0.02	1,26-dichloroethane	1.22	1.24	-0.02
1,27-dichloroethane	1.22	1.24	-0.02	1,27-dichloroethane	1.22	1.24	-0.02
1,28-dichloroethane	1.22	1.24	-0.02	1,28-dichloroethane	1.22	1.24	-0.02
1,29-dichloroethane	1.22	1.24	-0.02	1,29-dichloroethane	1.22	1.24	-0.02
1,30-dichloroethane	1.22	1.24	-0.02	1,30-dichloroethane	1.22	1.24	-0.02
1,31-dichloroethane	1.22	1.24	-0.02	1,31-dichloroethane	1.22	1.24	-0.02
1,32-dichloroethane	1.22	1.24	-0.02	1,32-dichloroethane	1.22	1.24	-0.02
1,33-dichloroethane	1.22	1.24	-0.02	1,33-dichloroethane	1.22	1.24	-0.02
1,34-dichloroethane	1.22	1.24	-0.02	1,34-dichloroethane	1.22	1.24	-0.02
1,35-dichloroethane	1.22	1.24	-0.02	1,35-dichloroethane	1.22	1.24	-0.02
1,36-dichloroethane	1.22	1.24	-0.02	1,36-dichloroethane	1.22	1.24	-0.02
1,37-dichloroethane	1.22	1.24	-0.02	1,37-dichloroethane	1.22	1.24	-0.02
1,38-dichloroethane	1.22	1.24	-0.02	1,38-dichloroethane	1.22	1.24	-0.02
1,39-dichloroethane	1.22	1.24	-0.02	1,39-dichloroethane	1.22	1.24	-0.02
1,40-dichloroethane	1.22	1.24	-0.02	1,40-dichloroethane	1.22	1.24	-0.02
1,41-dichloroethane	1.22	1.24	-0.02	1,41-dichloroethane	1.22	1.24	-0.02
1,42-dichloroethane	1.22	1.24	-0.02	1,42-dichloroethane	1.22	1.24	-0.02
1,43-dichloroethane	1.22	1.24	-0.02	1,43-dichloroethane	1.22	1.24	-0.02
1,44-dichloroethane	1.22	1.24	-0.02	1,44-dichloroethane	1.22	1.24	-0.02
1,45-dichloroethane	1.22	1.24	-0.02	1,45-dichloroethane	1.22	1.24	-0.02
1,46-dichloroethane	1.22	1.24	-0.02	1,46-dichloroethane	1.22	1.24	-0.02
1,47-dichloroethane	1.22	1.24	-0.02	1,47-dichloroethane	1.22	1.24	-0.02
1,48-dichloroethane	1.22	1.24	-0.02	1,48-dichloroethane	1.22	1.24	-0.02
1,49-dichloroethane	1.22	1.24	-0.02	1,49-dichloroethane	1.22	1.24	-0.02
1,50-dichloroethane	1.22	1.24	-0.02	1,50-dichloroethane	1.22	1.24	-0.02
1,51-dichloroethane	1.22	1.24	-0.02	1,51-dichloroethane	1.22	1.24	-0.02
1,52-dichloroethane	1.22	1.24	-0.02	1,52-dichloroethane	1.22	1.24	-0.02
1,53-dichloroethane	1.22	1.24	-0.02	1,53-dichloroethane	1.22	1.24	-0.02
1,54-dichloroethane	1.22	1.24	-0.02	1,54-dichloroethane	1.22	1.24	-0.02
1,55-dichloroethane	1.22	1.24	-0.02	1,55-dichloroethane	1.22	1.24	-0.02
1,56-dichloroethane	1.22	1.24	-0.02	1,56-dichloroethane	1.22	1.24	-0.02
1,57-dichloroethane	1.22	1.24	-0.02	1,57-dichloroethane	1.22	1.24	-0.02
1,58-dichloroethane	1.22	1.24	-0.02	1,58-dichloroethane	1.22	1.24	-0.02
1,59-dichloroethane	1.22	1.24	-0.02	1,59-dichloroethane	1.22	1.24	-0.02
1,60-dichloroethane	1.22	1.24	-0.02	1,60-dichloroethane	1.22	1.24	-0.02
1,61-dichloroethane	1.22	1.24	-0.02	1,61-dichloroethane	1.22	1.24	-0.02
1,62-dichloroethane	1.22	1.24	-0.02	1,62-dichloroethane	1.22	1.24	-0.02
1,63-dichloroethane	1.22	1.24	-0.02	1,63-dichloroethane	1.22	1.24	-0.02
1,64-dichloroethane	1.22	1.24	-0.02	1,64-dichloroethane	1.22	1.24	-0.02
1,65-dichloroethane	1.22	1.24	-0.02	1,65-dichloroethane	1.22	1.24	-0.02
1,66-dichloroethane	1.22	1.24	-0.02	1,66-dichloroethane	1.22	1.24	-0.02
1,67-dichloroethane	1.22	1.24	-0.02	1,67-dichloroethane	1.22	1.24	-0.02
1,68-dichloroethane	1.22	1.24	-0.02	1,68-dichloroethane	1.22	1.24	-0.02
1,69-dichloroethane	1.22	1.24	-0.02	1,69-dichloroethane	1.22	1.24	-0.02
1,70-dichloroethane	1.22	1.24	-0.02	1,70-dichloroethane	1.22	1.24	-0.02
1,71-dichloroethane	1.22	1.24	-0.02	1,71-dichloroethane	1.22	1.24	-0.02
1,72-dichloroethane	1.22	1.24	-0.02	1,72-dichloroethane	1.22	1.24	-0.02
1,73-dichloroethane	1.22	1.24	-0.02	1,73-dichloroethane	1.22	1.24	-0.02
1,74-dichloroethane	1.22	1.24	-0.02	1,74-dichloroethane	1.22	1.24	-0.02
1,75-dichloroethane	1.22	1.24	-0.02	1,75-dichloroethane	1.22	1.24	-0.02
1,76-dichloroethane	1.22	1.24	-0.02	1,76-dichloroethane	1.22	1.24	-0.02
1,77-dichloroethane	1.22	1.24	-0.02	1,77-dichloroethane	1.22	1.24	-0.02
1,78-dichloroethane	1.22	1.24	-0.02	1,78-dichloroethane	1.22	1.24	-0.02
1,79-dichloroethane	1.22	1.24	-0.02	1,79-dichloroethane	1.22	1.24	-0.02
1,80-dichloroethane	1.22	1.24	-0.02	1,80-dichloroethane	1.22	1.24	-0.02
1,81-dichloroethane	1.22	1.24	-0.02	1,81-dichloroethane	1.22	1.24	-0.02
1,82-dichloroethane	1.22	1.24	-0.02	1,82-dichloroethane	1.22	1.24	-0.02
1,83-dichloroethane	1.22	1.24	-0.02	1,83-dichloroethane	1.22	1.24	-0.02
1,84-dichloroethane	1.22	1.24	-0.02	1,84-dichloroethane	1.22	1.24	-0.02
1,85-dichloroethane	1.22	1.24	-0.02	1,85-dichloroethane	1.22	1.24	-0.02
1,86-dichloroethane	1.22	1.24	-0.02	1,86-dichloroethane	1.22	1.24	-0.02
1,87-dichloroethane	1.22	1.24	-0.02	1,87-dichloroethane	1.22	1.24	-0.02
1,88-dichloroethane	1.22	1.24	-0.02	1,88-dichloroethane	1.22	1.24	-0.02
1,89-dichloroethane	1.22	1.24	-0.02	1,89-dichloroethane	1.22	1.24	-0.02
1,90-dichloroethane	1.22	1.24	-0.02	1,90-dichloroethane	1.22	1.24	-0.02
1,91-dichloroethane	1.22	1.24	-0.02	1,91-dichloroethane	1.22	1.24	-0.02
1,92-dichloroethane	1.22	1.24	-0.02	1,92-dichloroethane	1.22	1.24	-0.02
1,93-dichloroethane	1.22	1.24	-0.02	1,93-dichloroethane	1.22	1.24	-0.02
1,94-dichloroethane	1.22	1.24	-0.02	1,94-dichloroethane	1.22	1.24	-0.02
1,95-dichloroethane	1.22	1.24	-0.02	1,95-dichloroethane	1.22	1.24	-0.02
1,96-dichloroethane	1.22	1.24	-0.02	1,96-dichloroethane	1.22	1.24	-0.02
1,97-dichloroethane	1.22	1.24	-0.02	1,97-dichloroethane	1.22	1.24	-0.02
1,98-dichloroethane	1.22	1.24	-0.02	1,98-dichloroethane	1.22	1.24	-0.02
1,99-dichloroethane	1.22	1.24	-0.02	1,99-dichloroethane	1.22	1.24	-0.02
1,100-dichloroethane	1.22	1.24	-0.02	1,100-dichloroethane	1.22	1.24	-0.02

Objectives

- Develop retention time models for inhalants of abuse and blood alcohol analysis of interest on four unique GC stationary phases, using web-based modeling software.
- Verify accuracy of the models against actual analyses.
- Demonstrate the utility of the modeling software for optimizing separations, speeding up methods, and transferring to other carrier gases and column dimensions.

Figure 2: The ProZDGC Chromatogram modeler can be accessed at www.restek.com/prozdgcm. Modeling can be done by entering in compounds of interest or searching by stationary phase and selecting compounds of interest that are modeled on that phase.



Figure 3: By default, the modeler will develop a speed optimized separation on the column you select. Conditions will be presented that the analyst can optimize the chromatogram. In addition, resolution for each compound is calculated. Below are some examples of various compound classes that are modeled on the BAC-1 and BAC-2 GC stationary phases. "a" Select solvent, "b" "Propanol" and "c" BAC resolution compounds with internal standards.

