

# Oven Temperature Sensor Assembly

cat.# 23039, 23040

Always follow the instrument manufacturer's instructions when installing or removing a GC oven temperature sensor. Before installing, carefully remove the protective cap from the sensor. Make sure that the tip of the sensor is not covered by the protective sleeve once it is installed.

**Questions about this or any other Restek product? Contact us or your local Restek representative ([www.restek.com/contact-us](http://www.restek.com/contact-us)).**

Restek patents and trademarks are the property of Restek Corporation.  
(See [www.restek.com/Patents-Trademarks](http://www.restek.com/Patents-Trademarks) for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2018 Restek Corporation.  
All rights reserved.  
Printed in the U.S.A.  
#203-03-043  
Rev. date: 04/18

[www.restek.com](http://www.restek.com)



**RESTEK**  
Pure Chromatography