

Personal Sample Inlet

cat.# 22123

Includes:

- 3' x 1/16" OD PTFE tubing
- Adjustable alligator clip
- PTFE reducing ferrule
- 1/4" stainless steel nut




PTFE reducing ferrule allows direct connection of 1/16" tubing and 1/4" flow controller.




Convenient clip can be moved along entire tubing length, for proper attachment in breathing zone.



open



closed



Contact Technical Service at 1-800-356-1688, 1-814-353-1300, ext. 4, or support@restek.com (or contact your Restek representative) if you have any questions about this product or any other Restek product.

RESTEK



© 2009 Restek Corporation.
All rights reserved.
Printed in the U.S.A.
#500-13-002
Rev. date: 03/12

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823 • www.restek.com
phone: 1-814-353-1300 • 1-800-356-1688 • fax: 1-814-353-1309 • e-mail: support@restek.com
Restek France • phone: +33 (0)1 60 78 32 10 • e-mail: restek@restekfrance.fr
Restek GmbH • phone: +49 (0)6172 2797 0 • e-mail: info@restekgmbh.de
Restek Ireland • phone: +44 (0)2890 814576 • e-mail: restekireland@aol.com
Restek Japan • phone: +81 (3)6459 0025 • e-mail: restekjapan@restek.com
Thames Restek U.K. LTD • phone: +44 (0)1494 563377 • e-mail: sales@thamesrestek.co.uk