

Mini-Can Stand

cat.# 22124

- Small, space-saving footprint—12" diameter x 6.5" high.



- Securely holds 2-³/₄" diameter cans.



- Collapsible, for easy storage and transport.



- 2 out of 3 stainless steel legs are movable, for reliable stability on uneven surfaces.

**Call Technical Service at 800-356-1688 or 814-353-1300, ext. 4 (or your Restek representative)
if you have any questions about this product or any other Restek product.**



© 2009 Restek Corporation



#500-13 [003]
Rev. date: 7/09

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823
phone: 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: +33 (0)1 60 78 32 10 • e-mail: restek@restekfrance.fr

Restek GmbH • phone: +49 (0)6172 2797 0 • e-mail: info@restekgmbh.de

Restek Ireland • phone: +44 (0)2890 814576 • e-mail: restekeurope@aol.com

Restek Japan • phone: +81 (3)6459 0025 • e-mail: ryosei.kanaguchi@restek.com

Thames Restek U.K. LTD • phone: +44 (0)1494 563377 • e-mail: sales@thamesrestek.co.uk