

Restek Soil Gas Samplers: Fewer Connections, Fewer Leaks, Fewer Errors



Restek soil gas samplers are designed for simplicity, reliability, and flexibility, so you can maximize data accuracy while minimizing time, expense, and user error. The rugged, unibody design has fewer connections, so it's quick and easy to set up and has fewer places where leaks could occur compared to typical samplers. The all-stainless-steel construction makes cleaning simple and means no contamination from degrading O-rings. Samplers ship leak tested (1×10^{-6} mL/sec) and ready to use to reduce your equipment prep time. In addition, all units are individually laser etched with a serial number and flow rate, so there's no more guesswork from lost tags or labels.

www.restek.com/SoilGas

Restek Corporation
www.restek.com
1-814-353-1300, ext. 3

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2023 Restek Corporation. All rights reserved. Printed in the U.S.A.

EVPR4130-UNV