

Mitigate Reruns and Keep Your PFAS Samples Flowing



New Resprep PFAS SPE cartridges are a two-in-one solution that simplifies sample preparation for complex matrices. They contain both polymeric WAX and carbon sorbents in a single, factory-packed cartridge, which simplifies and speeds up workflows for busy PFAS testing labs. Every lot of WAX material is tested for the 40 PFAS listed in EPA Method 1633 to ensure method detection limits can be met, and the ultra-clean CarboPrep Plus carbon bed ensures effective cleanup. In addition, an optional integrated filter aid reduces time-consuming reruns because it significantly mitigates clogging, so samples are processed more completely and consistently compared to using manually packed glass wool.

www.restek.com/ResprepPFAS

Restek Corporation
www.restek.com
1-814-353-1300, ext. 3