

Improve Productivity for Pesticides Analysis with More Effective SPE



New Resprep FL + CarboPrep Plus SPE cartridges enhance Florisil cleanup, boosting sample throughput for pesticides >10x. The dual-bed cartridge—containing both Florisil and carbon layers—drastically reduces instrument downtime by eliminating contaminants that degrade DDT and cause inlet performance requirements to fail. Cleaner extracts mean less instrument downtime, more accurate data, and increased lab productivity. Seamlessly integrating into existing workflows, Resprep FL + CarboPrep Plus SPE cartridges use the same solvents, volumes, and steps as Florisil cleanup procedures, making it simple to switch to an SPE that significantly improves lab productivity.

www.restek.com

Restek Corporation
www.restek.com
1-814-353-1300, ext. 3

For information on Restek patents and trademarks, visit www.restek.com/patents-trademarks. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2025 Restek Corporation. All rights reserved. Printed in the U.S.A.

EVPR4469-UNV