


Minimize Maintenance Time with Finger-Tight GC Fittings



Spending too much time on column installation and maintenance? Finger-tight fittings are faster and easier to install than traditional fittings, and they give you leak-tight connections with much less effort. In addition, finger-tight Hot Swap column nuts dissipate heat quickly, so you can change the column while the injector is still hot without burning your fingers or waiting for the instrument to cool down. Pair your Hot Swap fittings with finger-tight capillary gauges to ensure consistent insertion distances with every installation, so you can reduce variation and improve chromatographic performance.

Restek finger-tight fittings are now available for Agilent, PerkinElmer, Shimadzu, and Thermo GCs, so you can streamline operations and ordering for all your instruments. Switch today to increase sample throughput and improve lab productivity.

Get details at www.restek.com/FingerTight

Restek Corporation
www.restek.com
1-814-353-1300, ext. 3

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2023 Restek Corporation. All rights reserved. Printed in the U.S.A.

GNPR4115-UNV