

Restek LC Accessories

Simplify Your Routine LC Maintenance with These Handy Kits

For Agilent and Waters Systems

- Significant savings over instrument manufacturer prices.
- High-quality components in every kit.
- Wide range of options for HPLC systems and pumps.

www.restek.com/LC-Maintenance




RESTEK

Pure Chromatography

www.restek.com



25271



25798



25915

LC Maintenance Kits

Simplify routine maintenance with these handy kits for LC.

Description for Agilent LCs	Model #	Similar to Agilent part #	qty.	cat.#
Autosampler PM Kit Includes: rotor seal; piston seals (2); needle assembly; needle seat; finger caps (3)	1100, 1200	G1313-68709	kit	25271
Pump PM Kit Includes: PTFE frits (2); outlet cap; gold disk seal; active inlet cartridge; piston seals (4); glass solvent filters (2)	1050, 1100, 1200	G1311-68710	kit	25270
Pump PM Kit Includes: sapphire plunger (2); piston seals (2); outlet ball check valve; active inlet cartridge; PTFE frits (5)	1100, 1200	5065-4499	kit	25915
Description for Waters LCs	Model #	Similar to Waters part #	qty.	cat.#
ACQUITY BSM Pump Kit Includes: in-line filter assembly, stainless steel frit; tube assembly, transducer to check valve; tube assembly, SSV to in-line filter; primary check valve (2); wash seal, float flanged, (2); head plunger seal kit (2); sapphire plungers, (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µl; accumulator check valve, double ball & seat (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)	ACQUITY UPLC BSM	201000173	kit	25797
ACQUITY I2V BSM Pump Kit Includes: tube assembly, SSV to I2v; transducer to check valve, tube assembly, I2v; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); air filter, pump; air filter, pump handle; mixer assembly, 50 µl; filter frit cartridge, stainless steel; accumulator check valve, double ball & seat (2); check valve cartridge (2); solvent bottle filter, stainless steel (7); pump O-ring, (PTFE)	ACQUITY UPLC I2V BSM	201000197	kit	25798
ACQUITY H-Class QSM Pump Kit Includes: solvent bottle filter, stainless steel (5); 20 micron frit holder assembly; tube assembly, transducer to check valve; wash seal, float flanged (2); head plunger seal kit (2); sapphire plungers (2); check valve, double ball & seat (2); I2 check valve cartridge; air filter, door; pump O-ring, (PTFE); mixer assembly, 100 µl	ACQUITY H-CLASS QSM	201000233	kit	25799
2690/2695 Pump & Autosampler PM Kit Includes: sapphire plungers (2); seal wash plunger seals (4); head plunger seals (4); wash tube seals (4); battery, sparge diffuser; filter insert; face seals (4); solvent reservoir 10 µm filters (4); 250 µL syringe, check valve cartridges (2); wash tube seal; seal wash tube; PTFE washer; filter retainer; lower wash seal frit; needle wash frit; TFE washer; needle assembly; gold injector seals (2); stainless steel ferrule; compression screw	Alliance 2690, 2695	WAT270944	kit	25143
600 Pump PM Kit Includes: PerformancePLUS cartridges (4); sparge diffusers (4); Super Seals (2); solvent reservoir 10 µm filters (4); sapphire plungers (2); reference valve button; valve disk spacer; valve disk; TFE ball plug; TFE seat; ruby ball; inlet tube body assembly manifold insert; insert seal; belleville washers (2); flat washer	600 Pump	WAT052675	kit	25144
717 Autosampler PM Kit Includes: seal pack assembly; tube assembly (0.020" ID); needle; needle compression screw; 0.062 stainless steel ferrule; precolumn filter assembly; filter insert; 250 µL WISP syringe	717 Autosampler	WAT052669	kit	25145
616 Pump PM Kit Includes: sapphire plungers (2); sparge diffusers (4); solvent reservoir 10 µm filters (4); plunger seals (2); plunger wash seals (2); O-rings (2); backup rings (2); check valve cartridges (4)	616 Pump	WAT052672	kit	25146
515 Pump PM Kit Includes: PerformancePLUS check valves (4); sparge diffuser; solvent reservoir 10 µm filter; sapphire plungers (2); plunger seals (2); pivot inserts (2); pivot guides (2); washer (2); plunger springs (2); retaining rings (2)	515 Pump	WAT052587	kit	26519
1525 Pump PM Kit Includes: sapphire plungers (4); check valve cartridges (8); plunger seals (4); solvent reservoir 10 µm filters (2); reference valve button; valve disk spacer; valve disk	1525 Pump	201000114	kit	26430

To order, visit www.restek.com/LC-Maintenance

RESTEK
Pure Chromatography

Questions about this or any other Restek product?

Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2018 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS1527A-UNV