


Inlet Supplies for Shimadzu GCs with Split/Splitless Inlet



Shimadzu Plug Thermolite Plus Septa

Diameter	qty.	Similar to Part #	cat.#
Shimadzu Plug	50-pk.	Agilent 8010-0215	23872
Shimadzu Plug	100-pk.	Agilent 8010-0216	23873

Shimadzu Plug Premium Non-Stick BTO Septa

Diameter	qty.	cat.#
Shimadzu Plug	50-pk.	27098
Shimadzu Plug	100-pk.	27099

Viton O-Rings

for Shimadzu 17A, 2010, 2014, and 2030 GCs

Instrument	Max Temp	qty.	Similar to Part #	cat.#
for Shimadzu 17A, 2010, 2014, and 2030 GCs	250 °C	10-pk.	Shimadzu 036-11203-84	24899

Topaz 3.5 mm ID Single Taper Inlet Liner w/ Wool

for Shimadzu 17A, 2010, 2014, and 2030 GCs equipped with split/splitless inlets

ID x OD x Length	Packing	qty	cat.#
Single Taper, Premium Deactivation, Borosilicate Glass			
3.5 mm x 5.0 mm x 95 mm	Quartz Wool	5-pk.	23336

Topaz 3.5 mm ID Single Taper Inlet Liner

for Shimadzu 17A, 2010, 2014, and 2030 GCs equipped with split/splitless inlets

ID x OD x Length	qty	Similar to Part #	cat.#
Single Taper, Premium Deactivation, Borosilicate Glass			
3.5 mm x 5.0 mm x 95 mm	5-pk.	Shimadzu 221-48335-01	23321

Topaz 3.5 mm ID Precision Inlet Liner w/ Wool

for Shimadzu 17A, 2010, 2014, and 2030 GCs equipped with split/splitless inlets

ID x OD x Length	Packing	qty	cat.#
Precision, Premium Deactivation, Borosilicate Glass			
3.5 mm x 5.0 mm x 95 mm	Quartz Wool	5-pk.	23320

Topaz 3.5 mm ID Cyclo Single Taper Inlet Liner

for Shimadzu 17A, 2010, 2014, and 2030 GCs equipped with split/splitless inlets

ID x OD x Length	qty	cat.#
Cyclo Single Taper, Premium Deactivation, Borosilicate Glass		
3.5 mm x 5.0 mm x 95 mm	5-pk.	23456

View our complete range of Shimadzu inlet parts at www.restek.com

Our tool-free inlet parts are designed to make column installation and inlet maintenance easier and quicker, reducing the downtime of your Shimadzu GC.



Tool-free Column Installation Kit

Description	Includes	Instrument	qty.	cat.#
Tool-free Column Installation Kit	Inlet adaptor, retaining nut, back washer, and Hot Swap column nut	for Shimadzu 17A, 2010, 2014, and 2030 GCs	kit	27220
Tool-free Column Installation Kit (Siltek Treated)	Siltek-treated inlet adaptor, retaining nut, back washer, and Hot Swap column nut	for Shimadzu 17A, 2010, 2014, and 2030 GCs	kit	27221

Inlet Adaptor

Description	Instrument	qty.	Similar to Part #	cat.#
Inlet Adaptor	for Shimadzu 17A, 2010, 2014, and 2030 GCs	ea.	Shimadzu 221-42998-00	27222
Inlet Adaptor (Siltek Treated)	for Shimadzu 17A, 2010, 2014, and 2030 GCs	ea.	Shimadzu 221-42998-00	27223

Back Washer

Description	Instrument	qty.	Similar to Part #	cat.#
Back Washer	for Shimadzu 17A, 2010, 2014, and 2030 GCs	ea.	Shimadzu 201-30051-84	27226

Retaining Nut

Description	Instrument	qty.	Similar to Part #	cat.#
Retaining Nut	for Shimadzu 17A, 2010, 2014, and 2030 GCs	ea.	Shimadzu 201-30008-00	27225

Graphite Capillary Ferrules for 1/16-Inch Compression-Type Fittings

Ferrule ID	Fits Column ID	Material	qty.	cat.#
0.4 mm	0.025/0.05/0.075/0.10/0.15/0.18/0.20/0.25 mm (fused silica); 0.18 mm (MXT)	Graphite	10-pk.	20200
0.5 mm	0.32 mm (fused silica); 0.25/0.32 mm (MXT)	Graphite	10-pk.	20201
0.8 mm	0.45/0.53 mm (fused silica); 0.53 mm (MXT)	Graphite	10-pk.	20202

Vespal/Graphite Capillary Ferrules for 1/16-Inch Compression-Type Fittings

Ferrule ID	Fits Column ID	Material	qty.	cat.#
0.4 mm	0.025/0.05/0.075/0.10/0.15/0.18/0.20/0.25 mm (fused silica); 0.18 mm (MXT)	VG1, 85% Vespal/15% Graphite	10-pk.	27062
0.5 mm	0.32 mm (fused silica); 0.25/0.32 mm (MXT)	VG1, 85% Vespal/15% Graphite	10-pk.	27064
0.8 mm	0.45/0.53 mm (fused silica); 0.53 mm (MXT)	VG1, 85% Vespal/15% Graphite	10-pk.	27067

Hot Swap Column Nut for 1/16" Fittings

Description	qty.	cat.#
Hot Swap Column Nut for 1/16" Fittings	ea.	27224

Capillary Installation Gauge for 1/16" Fittings

Description	qty.	cat.#
Capillary Installation Gauge for 1/16" Fittings	ea.	27227

For other Hot Swap configurations for Shimadzu GCs, visit www.restek.com/FingerTight

RESTEK
Pure Chromatography

Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS2129D-UNV