


Inlet Supplies for PerkinElmer GCs with CAP Inlet



Septa (11 mm) for PerkinElmer GCs

Septum Diameter	50-pk. cat.#	100-pk. cat.#
Thermolite Plus Septa (usable to 350 °C inlet temp.)		
11 mm (7/16"), with CenterGuide	23864	23865
Premium Non-Stick BTO Septa (usable to 400 °C inlet temp.)		
11 mm (7/16"), with CenterGuide, predrilled	27090	27091

Silicone O-Rings

Description	Max. Temp.	Similar to PE Part #	qty.	cat.#
Silicone O-Rings for PerkinElmer Auto SYS XL or Clarus w/CAP Injector	250 °C	N6101374	10-pk.	20262

Topaz Inlet Liner

ID x OD x Length	qty.	cat.#
Precision, Premium Deactivation, Borosilicate Glass with Quartz Wool 4.0 mm x 6.2 mm x 92.1 mm	5-pk.	23450

Capillary Ferrules for 1/16-Inch Compression-Type Fittings

Ferrule ID	Fits Column ID	qty.	VG2 (60/40)
0.4 mm	0.10/0.15/0.18/0.25/0.28 mm	10-pk.	20211
0.5 mm	0.32 mm	10-pk.	20212
0.8 mm	0.45/0.53 mm	10-pk.	20213

Visit www.restek.com for a complete product listing.

RESTEK
Pure Chromatography

Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2018 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS2130B-UNV