


Inlet Supplies for Thermo TRACE Ultra GCs



Septa (17 mm)

Septum Diameter

48-pk.
cat.#

Premium Non-Stick BTO Septa (usable to 330 °C inlet temp.)

17 mm with CenterGuide

27096

Graphite Sealing Ring

Description

Max.
Temp.

Similar to
TS Part #

qty.

cat.#

Graphite Sealing Ring
for Thermo Scientific TRACE, 8000,
8000 TOP & Focus SSL

450 °C
450 °C

290-334-06
290-334-06

ea.
2-pk.

21898
21899

Topaz Inlet Liner

ID x OD x Length

qty.

cat.#

Single Taper, Premium Deactivation, Borosilicate Glass with Quartz Wool
3.0 mm x 8.0 mm x 105 mm

5-pk.

23442

ID x OD x Length

qty.

cat.#

Single Taper, Premium Deactivation, Borosilicate Glass
3.0 mm x 8.0 mm x 105 mm

5-pk.

23441

ID x OD x Length

qty.

cat.#

Precision, Premium Deactivation, Borosilicate Glass with Quartz Wool
5.0 mm x 8.0 mm x 105 mm

5-pk.

23327

M4 Graphite Ferrules

Ferrule
ID

Fits Column ID

Similar to
TS part #

2-pk.
cat.#

10-pk.
cat.#

0.4 mm 0.025/0.05/0.075/0.10/0.15/0.18/
0.20/0.25 mm (fused silica);
0.18 mm (MXT)

29013488 (2-pk.)
29053488 (10-pk.)

20280

20281

0.5 mm 0.32 mm (fused silica);
0.25/0.32 mm (MXT)

29013487 (2-pk.)
29053487 (10-pk.)

20282

20283

0.8 mm 0.45/0.53 mm (fused silica);
0.53 mm (MXT)

29013486 (2-pk.)
29053486 (10-pk.)

20284



20285

Visit www.restek.com for
a complete product listing.

RESTEK

Pure Chromatography

www.restek.com



Reducing Nuts

Description	qty.	cat.#
Adaptor for Capillary Column on Split/Splitless Injector for Thermo TRACE and Focus SSL; for use w/Standard 1/16" type Ferrules	ea.	24917
Siltek-Treated Adaptor for Capillary Column on Split/Splitless Injector for Thermo TRACE and Focus SSL; for use w/Standard 1/16" type Ferrules	ea.	20543

Standard 1/16" Ferrules

Ferrule ID	Fits Column ID	qty.	Graphite
0.4 mm	0.025/0.05/0.075/0.10/0.15/0.18/0.20/0.25 mm (fused silica); 0.18 mm (MXT)	10-pk.	20200
0.5 mm	0.32 mm (fused silica); 0.25/0.32 mm (MXT)	10-pk.	20201
0.8 mm	0.45/0.53 mm (fused silica); 0.53 mm (MXT)	10-pk.	20202

Column Nuts

Description	qty.	cat.#
Hot Swap Column Nut for 1/16" Fittings for Agilent GCs (Except Intuvo); PerkinElmer Clarus 590/690 and GC2400 GCs; Shimadzu 17A, 2010, 2010 PLUS, 2014, and 2030 GCs; Thermo TRACE 1300/1310 and 1600/1610 GCs	ea.	27224

Installation Gauge

Description	qty.	cat.#
Capillary Installation Gauge for 1/16" Fittings for Agilent GCs (Except Intuvo); PerkinElmer Clarus 590/690 and GC2400 GCs; Shimadzu 17A, 2010, 2010 PLUS, 2014, and 2030 GCs; Thermo TRACE 1300/1310 and 1600/1610 GCs	ea.	27227

Visit www.restek.com for a complete product listing.

RESTEK
Pure Chromatography

Questions? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe. To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS2132D-UNV