


Inlet Supplies for PerkinElmer 590/690 and GC2400 GCs



Septa (11 mm)

Septum Diameter	50-pk. cat.#	100-pk. cat.#
Thermolite Plus Septa (usable to 350 °C inlet temp.)		
11 mm (7/16") with CenterGuide	23864	23865
Premium Non-stick BTO Septa (usable to 400 °C inlet temp.)		
11 mm (7/16") with CenterGuide, predrilled	27090	27091

Viton O-Rings

Description	Max Temp	Similar to Agilent Part #	10-pk. cat.#	50-pk. cat.#
Viton O-Rings for Agilent GCs	300 °C	5188-5365	22241	22242

Topaz Inlet Liner

ID x OD x Length	qty.	cat.#
Single Taper, Premium Deactivation, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23800
Single Taper, Premium Deactivation, Borosilicate Glass		
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23796
Straight, Premium Deactivation, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23801
Precision, Premium Deactivation, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23799



Dual Vespel Ring Inlet Seals Washerless, Leak-Tight Seals

	2-pk. cat.#	10-pk. cat.#	50-pk. cat.#
0.8 mm ID Dual Vespel Ring Inlet Seal			
Gold-Plated	21240	21241	23418

Replacement Inlet Seals with Washers

Single-column 0.8 mm ID (Opening)	Similar to Agilent part #	2-pk. cat.#	10-pk. cat.#	50-pk. cat.#
Gold-Plated	5188-5367	21317	21318	23415

View our complete range of PerkinElmer inlet parts at www.restek.com



Reducing Nut

Description	qty.	cat.#
Direct Replacement Reducing Nut for Agilent 5890/6850/6890/7890/8890 GCs	ea.	22078

Compact Ferrules

Ferrule ID	Fits Column ID	qty.	Graphite
0.4 mm	0.025/0.05/0.075/0.10/0.15/0.18/0.20/0.25 mm (fused silica); 0.18 mm (MXT)	10-pk.	20250
0.5 mm	0.32 mm (fused silica); 0.25/0.32 mm (MXT)	10-pk.	21007
0.8 mm	0.45/0.53 mm (fused silica); 0.53 mm (MXT)	10-pk.	20252

Standard 1/16" Ferrules

Ferrule ID	Fits Column ID	qty.	Graphite
0.4 mm	0.025/0.05/0.075/0.10/0.15/0.18/0.20/0.25 mm (fused silica); 0.18 mm (MXT)	10-pk.	20200
0.5 mm	0.32 mm (fused silica); 0.25/0.32 mm (MXT)	10-pk.	20201
0.8 mm	0.45/0.53 mm (fused silica); 0.53 mm (MXT)	10-pk.	20202

Column Nuts

Description	qty.	cat.#
Capillary Column Nut, Stainless Steel, for Use w/Agilent Compact Ferrules	2-pk.	21884
Capillary Column Nut, Brass, for Use w/Agilent Compact Ferrules	2-pk.	21878
Capillary Column Nut, Stainless Steel, for Use w/Standard 1/16" type Ferrules	2-pk.	20883
Capillary Column Nut, Brass, for Use w/Standard 1/16" type Ferrules	2-pk.	21879

Visit www.restek.com for a complete product listing.

RESTEK

For information on Restek patents and trademarks, visit www.restek.com/patents-trademarks To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.
© 2025 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS4192A-UNV