

Streamline Nitrosamines Testing with Restek's New Multicomponent Standard



Comprised of the most frequently targeted nitrosamine compounds for pharmaceutical impurities testing, Restek's new, multicomponent Nitrosamines 7 Standard contains seven key nitrosamines in a single ampul, significantly reducing the complexity and challenges in sourcing the compounds you need. Explore the benefits of this seven-in-one nitrosamine certified reference standard (CRM) and new internal standard to simplify sourcing, reduce standard preparation, and to help ensure the accuracy of your data.

- Seven critical nitrosamines in a single ampul.
- Formulated for optimal stability with a concentration of 100 µg/mL for each analyte, allowing for flexibility in constructing calibration curves.
- Two independently manufactured lots available with verified lot-to-lot agreement.
- Internal standard available at a concentration of 1000 µg/mL.

www.restek.com/p/36025

Restek Corporation
www.restek.com
1-814-353-1300, ext. 3

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries. To unsubscribe from future Restek communications or to update your preferences, visit www.restek.com/subscribe To update your status with an authorized Restek distributor or instrument channel partner, please contact them directly.

© 2024 Restek Corporation. All rights reserved. Printed in the U.S.A.

PHPR4225-UNV